Instructions For Use

SAS-MX Serum Protein SB

Cat. No. 100200

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Fiche technique

Réf. 100200

SAS-MX Serum-protein SB

Anleitung

Kat. Nr. 100200

SAS-MX Sieroproteine SB

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Cod. 100200

SAS-MX Proteínas Séricas SB

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Contents

English	1
Français	5
Deutsch	10
Italiano	15
Español	20

INTENDED PURPOSE

The SAS-MX Serum Protein SB kit is intended for the separation and quantitation of serum proteins by agarose gel electrophoresis.

Serum contains over 100 individual proteins, each with a specific set of functions which are subject to specific variation in concentration under different pathological conditions¹.

Since the introduction of moving boundary electrophoresis by Tiselius², and the subsequent use of zone electrophoresis, serum proteins have been fractionated on the basis of their charge at a particular pH. The SAS-MX Serum Protein kit separates serum proteins into 6 main classes (albumin, alpha 1-globulin, alpha 2-globulin, beta 1-globulin, beta 2-globulin, and gamma globulin) according to charge in an agarose gel. The proteins are then stained to allow visualisation and quantitative interpretation. Each of the classical electrophoretic zones, with the exception of albumin, normally contains 2 or more components. The relative proportions of these fractions have proven to be useful aids in the diagnosis and prognosis of certain disease states^{3,5}.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION**1. SAS-MX Serum Protein SB Gel**

Contains agarose in a Tris / Barbitol buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

2. Tris / Barbitol Buffer Concentrate

Contains barbital and sodium barbital with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well. Buffer salts may crystallize slightly on standing. Wash any crystals from the bottle with diluted buffer to ensure complete dissolution.

3. Acid Blue Stain Concentrate

Contains concentrated Acid Blue stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

4. Destain Solution Concentrate

Contains concentrated Destain Solution. Dilute the contents of the bottle to 2 litres with purified water. Store in a tightly stoppered bottle.

5. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates and Blotters A and C to complete 10 gels.

STORAGE AND SHELF-LIFE**1. SAS-MX Serum Protein SB Gel**

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Barbital Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

3. Acid Blue Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration of the stain solution.

4. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. The diluted destain solution is stable for 6 months at 15...30°C. Cloudiness may indicate deterioration of the destain solution.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Drying oven with forced air capable of 60...70°C

Methanol

Purified water

SAMPLE COLLECTION AND PREPARATION

Freshly collected serum is the specimen of choice. Whilst samples can be stored at 15...30°C for up to 4 days, 2...6°C for up to 2 weeks or 6 months at -20°C for standard 5 band patterns⁶, the beta 2-globulin band (Complement C3) will degrade rapidly on storage. Urine and CSF can also be used following a suitable concentration step (50 - 100X). The use of plasma will result in a fibrinogen band between the beta and gamma fractions.

Interfering Factors: 1) Haemolysis may cause false elevation in the alpha 2 and beta fractions.
2) Inaccurate results may be obtained on specimens left uncovered, due to evaporation.

Samples / controls should be diluted 1:4 (1+3) with buffer before use. For more accurate albumin quantitations by densitometry, the sample dilution can be increased to 1:8 (1+7). The exact amount should be determined in each individual laboratory. The use of increased sample dilution will increase the quantitation of albumin, but will reduce the sensitivity of the method for minimonoclonal bands.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Remove the overlay and blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 5.
3. Apply 3µl of sample to each slit and allow to absorb for 4 minutes.
4. Whilst the samples are absorbing, pour 25ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from step 2 and remove both blotter and template.

6. Position the gel in the chamber agarose side down, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 80 volts, 25 minutes.
8. Following electrophoresis, dry the gel at 60...70°C. **NOTE:** The drying of the gel should take no more than 5-10 minutes to prevent diffusion of bands. If this cannot be achieved, fix the gel for 5 minutes in methanol prior to drying.
9. Immerse the dry gel in stain solution for 10 minutes.
10. Destain the gel in 2 x 2 minute washes of destain solution or until the background is clear.
11. Wash the gel briefly in purified water and dry.

INTERPRETATION OF RESULTS

It is recommended that any evaluation of the gels is performed against normal values produced for this method in each individual laboratory.

For a complete review of serum protein evaluation, see Ritzmann, S.E, 1982⁵. Studies show that the values are the same for both males and non-pregnant females. Some differences are seen in pregnant females at term and women on oral contraceptives.

Age has some effect on normal levels. Cord blood has a decreased total protein, albumin, alpha2 and beta fractions; slightly increased alpha I and normal or increased gamma fraction (largely of maternal origin). The gamma globulins drop rapidly until about 3 months of age, while other fractions have reached adult levels by this time. Adult levels of the gamma globulins are not reached until 10-16 years of age. The albumin decreases and beta globulin increases over the age of 40.

1. Qualitative Evaluation:

The gels may be visually inspected for the presence or absence of particular bands of interest.

2. Quantitative Evaluation:

Scan the gels, side down, at 595nm.

In either case, an elevation or decrease in particular serum components or the detection of unusual serum components require further investigation. The completed SAS-MX Serum Protein SB gel is stable for an indefinite period of time.

QUALITY CONTROL

Kemtrol Serum Controls (Cat. No. 7024 and 7025) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert provided for appropriate assay values.

LIMITATIONS

Since all electrophoresis procedures are non-linear, it is important to follow these Instructions For Use closely to ensure optimal resolution and reproducible results. Failure to follow these Instructions For Use may affect the results obtained.

REFERENCE VALUES

Using 20 normal specimens from male and female donors with an age range of 20-59 years, the following normal ranges were obtained (these are presented as a guideline only):

Protein Fraction	Mean (%)	\pm 2SD Range
Albumin	56.8	52.1 - 61.5
Alpha-1	3.6	2.5 - 4.7
Alpha-2	9.8	6.7 - 12.9
Beta-1	8.8	5.9 - 11.7
Beta-2	5.6	3.9 - 7.3
Gamma	15.4	11.9 - 18.8

PERFORMANCE CHARACTERISTICS

Reproducibility

	Within Gel(n=10)		Between Gel(n=100)	
	Mean (%)	CV(%)	Mean(%)	CV(%)
Albumin	59.8	3.6	59.2	1.9
Alpha-1	2.9	5.2	3.5	0.3
Alpha-2	7.5	3.8	8.0	0.6
Beta-1	8.5	3.4	9.2	0.8
Beta-2	6.3	9.6	6.0	0.6
Gamma	15.0	8.8	14.0	1.1

Sensitivity

The method is sensitive to 0.04g/L per band, determined as the lowest concentration of protein which was evident as a discrete band on the completed gel.

Linearity

The linearity of the method is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.

BIBLIOGRAPHY

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Apparatus for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
4. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, pages 579-582, 1986.
5. Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
6. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, (3rd Edition), W.B. Saunders Co., Philadelphia, page 524, 1995.

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LITERATUR

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290. 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
4. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, Seite 579-582, 1986.
5. Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
6. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, (3rd Edition), W.B. Saunders Co., Philadelphia, page 524, 1995.

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BIBLIOGRAFIA

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Apparatus for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Ritzmann, S.E and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins, Qualitative and Quantitative Assays' in Laboratory Medicine, Harper and Row, Inc., Hagerstown, 1979.
4. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, pages 579-582, 1986.
5. Ritzmann, S.E. (Ed.), Protein Abnormalities Volume I : Physiology of Immunoglobulins - Diagnostic and Clinical Aspects', Allen R. Liss, Inc., New York, 1982.
6. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, (3rd Edition), W.B. Saunders Co., Philadelphia, page 524, 1995.

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