SAS-MX Immunofix

Instructions For Use REF. I00300

SAS-MX Immunofix
Fiche technique

SAS-MX Immunofixation
Anleitung

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INTENDED PURPOSE

The SAS-MX IFE kit is intended for the separation and identification of monoclonal gammopathies by agarose gel electrophoresis with the Helena BioSciences SAS-MX electrophoresis chamber.

Immunofixation electrophoresis (IFE) is a two stage procedure using high resolution agarose electrophoresis in the first stage, and immunoprecipitation in the second phase.

The greatest demand for IFE is in the clinical laboratory, where it is used primarily for the detection of monoclonal gammopathies. A monoclonal gammopathy is a primary disease state in which a single clone of plasma cells produce elevated levels of an immunoglobulin of a single class and type. Such immunoglobulins are referred to as monoclonal proteins, M-proteins or paraproteins. Their presence may be of a benign nature or of uncertain significance. In some cases, they are indicative of a malignancy, such as multiple myeloma or Waldenström's macroglobulinaemia. Differentiation must be made between polyclonal and monoclonal gammopathies, as polyclonal gammopathies are a secondary disease state due to clinical disorders such as chronic liver disease, collagen disorders, rheumatoid arthritis and chronic infection.

Alfonso first described immunofixation in the literature in 1964¹. Alper and Johnson published a more practical procedure in 1969, and published a number of studies utilising this technique^{2,4}. Immunofixation has been used as a procedure for the investigation of immunoglobulins since 1976^{5,6}.

The SAS-MX IFE kit separates serum proteins according to charge in an agarose gel. The proteins are then incubated with monospecific antisera, washed and stained to allow visualization of the immunoprecipitate for qualitative interpretation.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-MX IFE Gel

Contains agarose in a Tris / Barbitol / Aspartate buffer with sodium azide as preservative. The gel is ready for use as packaged.

2. Tris / Barbitol Buffer Concentrate

Contains concentrated Tris / Barbitol buffer with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well.

3. Acid Blue Stain Concentrate

Contains concentrated Acid Blue stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

4. SAS-MX IFE Antisera Kit

Contains SP protein fixative (containing acetic acid and sulphosalicylic acid), and monospecific antisera to human immunoglobulins - IgG, IgA, IgM, Kappa Light Chain (free & bound) and Lambda Light Chain (free & bound). All antisera contain sodium azide as a preservative. The antisera are ready for use as packaged.

5. Destain Solution Concentrate

Contains concentrated Destain Solution. Dilute the contents of the bottle in 2 litres of purified water. Store in a tightly stoppered bottle.

6. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates, Antisera Templates and Blotters A, B, C, D and X to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-MX IFE Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Barbital Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

3. Acid Blue Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration of the stain solution.

4. SAS-MX IFE Antisera Kit

The antisera kit should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

5. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted destain solution is stable for 6 months at 15...30°C. Cloudiness may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 9400 IFE Control Kit

Cat. No. 4062 Incubation Chamber

Cat. No. 9035 Immuno SuperPress / Cat. No. 5014 Development Weight

Drying oven with forced air capable of 60...70°C

Saline Solution (0.85% NaCl)

Purified water

SAMPLE COLLECTION AND PREPARATION

Freshly collected serum is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 72 hours or 2 weeks at -20°C. Urine and CSF can also be used following a suitable concentration step:

Total Protein (mg/dL)	Conc. Factor
<50	100 x
50 - 100	50 x
100 - 300	25 x
300 - 600	10 x
600 - 1000	5 x

These concentrated CSF and urine samples should be used neat on the gel.

Serum samples should be diluted 1 + 1 in the SP Lane, and diluted 1 + 9 with saline solution for all 5 immunoglobulin lanes.

When typing minimono-clonal bands, the sample for the immunoglobulin lanes should be undiluted serum.

When typing high concentration monoclonal bands, the sample for the immunoglobulin lanes should be diluted more than 1 + 9 to prevent prozoning.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Remove the overlay and blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 5.
3. Apply 3µl of sample to each slit and allow to absorb for 5 minutes.
4. Whilst the samples are absorbing, pour 40ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from Step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side up, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 120 volts, 25 minutes.
8. Following electrophoresis, place the gel into the incubation chamber containing a wet blotter and position the antiserum application template onto the gel surface. **NOTE:** Ensure good contact of template and gel.
9. Apply 1µl of the appropriate IFE Control to the wells in the gel. Close the incubation chamber and allow to absorb completely (approximately 3 minutes).
10. Apply 2 drops of Protein Fixative to the SP lane and 2 drops of the appropriate antiserum into the immunoglobulin lanes. Ensure even distribution of antiserum in the lane (including the control wells) by rocking the gel.
11. Incubate the gel at 15...30°C for 10 minutes.
12. Following incubation, remove the antiserum template by briefly washing the gel in saline solution with gentle agitation.
13. Place the gel on a blotter D agarose side up. Place a blotter B (wetted in saline) onto the surface of the gel followed by a blotter X. Press the gel in an IFE SuperPress for 5 minutes or use a

development weight for 10 minutes.

14. Remove the blotters and place the gel in saline solution for 4 minutes with gentle agitation.
15. Place the gel on a blotter D agarose side up. Place a blotter B (wetted in saline) onto the surface of the gel followed by a blotter D. Press the gel in an IFE Superpress for 1 minute or use a development weight for 3 minutes.
16. Remove the blotters and dry the gel at 60...70°C.
17. Immerse the dry gel in stain solution for 4 minutes.
18. Destain the gel in 2 x 2 minute washes of destain solution.
19. Wash the gel briefly in purified water and dry.

INTERPRETATION OF RESULTS

The majority of monoclonal proteins migrate in the cathodic, gamma region of the protein pattern, but due to their abnormal nature, they may migrate anywhere within the globulin region on protein electrophoresis. The monoclonal protein band on the immunofixation pattern will occupy the same position and shape as the abnormal band on the serum protein pattern. The abnormal protein is identified by the antiserum type it reacts with.

When low concentrations of abnormal protein are present, the abnormal band may appear as a band within the normal polyclonal immunoglobulin. A band can also be seen within a polyclonal background when there is a large polyclonal immunoglobulin presence also.

The publication 'Immunofixation for the Identification of Monoclonal Gammopathies' is available from Helena BioSciences on request.

LIMITATIONS

1. Antigen Excess

Antigen excess will occur if there is not a slight antibody excess or antigen / antibody equivalence at the site of precipitation. Antigen excess in IFE is usually due to an excess of the immunoglobulin in the patient sample. Antigen excess is characterised by prozoning (unstained areas in the centre of the immunofixed protein band, with staining around the edges). A higher dilution of the sample should be used in this event to optimise the immunoglobulin concentration.

2. Non-Specific Precipitation in All Immunoglobulin Lanes

Occasionally a completed IFE plate exhibits a precipitate band in the same position in every pattern across the plate. This may result from:

a) IgM monoclonal immunoglobulins.

IgM monoclonal proteins can adhere to the gel matrix. A band will appear in all 5 antiserum lanes of the gel. However, where the band reacts with a specific antiserum for the heavy chain and light chain, there will be an increase in size and staining intensity of the band, allowing the immunoglobulin type to be identified. Additional dilution of the sample will improve the discrimination between the IgM-antibody reaction and the non-specific staining of precipitated IgM protein in other lanes, simplifying the diagnosis.

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Samples with high titres of Rheumatoid Factor or other immune complexes may show a precipitate band at the sample application point. Reducing the sample with DTT or β -2-mercaptoethanol can eliminate this non-specific reaction (Mix 190 μ L of diluted serum to 10 μ L of 1% (w/v) DTT in 0.85% saline solution or mix 100 μ L of serum with 10 μ L of a 1:10 dilution of β -2-mercaptoethanol in water. Perform the IFE as usual. **NOTE:** Always work in a fume hood when using β -2-mercaptoethanol).

c) Fibrinogen.

Fibrinogen, if present in the sample, will show as a discrete band in all lanes of the immunofixation pattern. Fibrinogen is present in plasma, and sometimes in the serum of patients on anticoagulant therapy.

3. Reaction With Kappa or Lambda Light Chain Antisera but No Reaction with IgG, IgA or IgM Heavy Chain Antisera.

Samples showing this pattern may either have a free light chain monoclonal gammopathy or they may have an IgD or IgE monoclonal protein. In this situation, the IFE should be repeated, substituting IgD and IgE antisera for two of the other heavy chain antisera. Failure to obtain a reaction with IgD or IgE antisera would be indicative of free light chain disease.

4. Band In Gamma Region Showing No Reactivity With IFE Antisera.

C Reactive Protein (CRP) may be detected in patients with acute inflammatory response^{7,8}. CRP appears as a narrow band at the cathodic end of the serum protein pattern. Elevated Alpha-1-Antitrypsin and Haptoglobin are supportive evidence for CRP. Patients with a CRP band will probably have an elevated level when assayed for CRP. A narrow band on the point of sample application can sometimes be seen which can be caused by chylomicrons in the serum or precipitated protein in samples which have been stored frozen.

5. Non-Reactivity With Kappa and Lambda Antisera.

Occasionally a sample will have a reaction with a heavy chain antiserum but no light chain reaction is obvious. In this situation, the following need to be ruled out - a) Heavy chain disease, b) Very high concentrations of light chains, leading to antigen excess, c) Low concentrations of light chains, d) A typical light chain molecule that does not react with the antiserum, e) Light Chains with 'hidden' light chain determinants (as sometimes seen with IgA and IgD). To obtain definitive results, testing may include a) A higher or lower dilution of the sample to optimise the antibody / antigen equivalence, b) Antisera from more than one manufacturer to aid in the identification of atypical immunoglobulins, and c) Treat the sample with β -2-mercaptoethanol to 'reveal' the light chains.

PERFORMANCE CHARACTERISTICS

A series of samples were tested and compared to another commercially available test kit - both kits showed equivalent results.

OPTIONAL EXTRA MATERIALS

Cat. No. 9249 Antiserum to Human IgD

Cat. No. 9250 Antiserum to Human IgE

Cat. No. 9412 Antiserum to Human Free Kappa Light Chain

Cat. No. 9413 Antiserum to Human Free Lambda Light Chain

BIBLIOGRAPHY

1. Afonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
2. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
3. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.

4. Johnson, A.M., 'Genetic Typing of α -1-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of PiM.', J. Lab. Clin. Med., 1976; 87 : 152-163.
5. Cawley, L.P., Minard, B.J., Tourtellotte, W.W., Ma, B.I. and Chelle, C. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
6. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
7. Jeppsson, J.O., Laurell, C.B., and Franzén, B. 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
8. Killingsworth, L.M., Cooney, S.K., and Tylia, M.M. 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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BIBLIOGRAPHIE

1. Afonso, E., 'Quantitation Immuno-electrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
2. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
3. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
4. Johnson, A.M., 'Genetic Typing of a1-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of PiM.', J. Lab. Clin. Med., 1976; 87 : 152-163.
5. Cawley, L.P., Minard, B.J., Tourtellotte, W.W., Ma, B.I. and Chelle, C. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
6. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
7. Jeppsson, J.O., Laurell, C.B., and Franzén, B. 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
8. Killingsworth, L.M., Cooney, S.K., and Tyllia, M.M. 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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LITERATUR

1. Alfonso, E., 'Quantitation Immuno-electrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
2. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vo. Sang., 1969; 17 : 445-452.
3. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
4. Johnson, A.M., 'Genetic Typing of Alpha(1)-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of P.M.', J. Lab. Clin. Med., 1976; 87 : 152-163.
5. Cawley, L.P. et al. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
6. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
7. Jeppsson, J.E. et al., 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
8. Killingsworth, L.M. et al., 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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BIBLIOGRAFIA:

1. Afonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
2. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
3. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
4. Johnson, A.M., 'Genetic Typing of α 1-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of PiM.', J. Lab. Clin. Med., 1976; 87 : 152-163.

5. Cawley, L.P., Minard, B.J., Tourtellotte, W.W., Ma, B.I. and Chelle, C. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
6. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
7. Jeppsson, J.O., Laurell, C.B., and Franzén, B. 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
8. Killingsworth, L.M., Cooney, S.K., and Tyllia, M.M. 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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BIBLIOGRAFIA

1. Afonso, E., 'Quantitation Immuno-electrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
2. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
3. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
4. Johnson, A.M., 'Genetic Typing of a I-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of PiM.', J. Lab. Clin. Med., 1976; 87 : 152-163.
5. Cawley, L.P., Minard, B.J., Tourtellotte, W.W., Ma, B.I. and Chelle, C. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
6. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
7. Jeppsson, J.O., Laurell, C.B., and Franzén, B. 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
8. Killingsworth, L.M., Cooney, S.K., and Tyllia, M.M. 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1277P 2007/09 (5)