helena | BioSciences
Europe
www.helena-biosciences.com

Helena BioSciences Europe
Colima Avenue
Sunderland Enterprise Park
Sunderland
SR5 3XB
tel: +44 (0) 191 549 6064
fax: +44 (0) 191 549 6271
email: info@helena-biosciences.com

HL-2-1301P 2003.08 (8)

Other Helena BioSciences Europe offices:

Helena BioSciences Europe
6 Rue Charles Cros-ZAE
95320 Saint Leu La Foret
France
tel: +33 13 995 9292
fax: +33 13 995 6891
email: helena@helena.fr

Helena BioSciences Europe
Via Enrico Fermi, 24
20090 Assago (Milano)
Italy
tel: +39 02 488 1951
or +39 02 488 2141
fax: +39 02 488 2677

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SAS-MX URINE PROTEIN

INTENDED PURPOSE

The SAS-MX Urine Protein Kit is intended for the screening of unconcentrated urine samples by agarose gel electrophoresis.

Urinary proteins are derived primarily from plasma proteins that filter through the kidney. The appearance of abnormal plasma proteins in the urine is of great value in evaluating renal function. The appropriate study of proteinuria should include quantitative and qualitative assessment of the type and amount of proteins excreted^{1,4}. Electrophoresis is a technique which allows the differentiation of several types of proteinuria - Physiological, Glomerular (selective and non-selective), Tubular and proteinuria associated with Dysglobulinaemias^{1,5}.

The SAS-MX Urine Protein kit separates urine proteins according to charge in an agarose gel. The proteins are then stained to allow visualisation and quantitative or qualitative interpretation. The high sensitivity stain used in the kit allows most urine samples to be tested without prior concentration and the use of the Urine Protein Internal Standard allows an approximation of urine total protein, or individual band protein concentrations to be determined by densitometry

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information

COMPOSITION

1. SAS-MX Urine Protein Gel

Contains agarose in a Tris / Barbitol buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

2. Tris / Barbitol Buffer Concentrate

Contains barbitol and sodium barbitol with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well. Buffer salts may crystallize slightly on standing. Wash any crystals from the bottle with diluted buffer to ensure complete dissolution.

3. Urine Protein Stain

Contains Coomassie Blue stain powder. The powdered stain is stable until the expiry date indicated on the label. Dissolve the contents of the vial in 1 litre of 50% methanol in water and stir overnight. Acidify by adding 20ml of glacial acetic acid. Mix well and filter before use.


4. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates and Blotters A and C to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-MX Urine Protein Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Barbital Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

3. Urine Protein Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber
 Cat. No. 1525 EPS600 Power Supply
 Cat. No. 3032 Urine Protein Internal Standard
 Drying oven with forced air capable of 50...60°C
 Fixative Solution: Mix 500ml of purified water and 500ml methanol. Add 100ml of glacial acetic acid. Store in a tightly stoppered bottle.
 Destain Solution: Mix 50ml of glacial acetic acid and 950ml of purified water. Store in a tightly stoppered bottle.
 Purified water

SAMPLE COLLECTION AND PREPARATION

Freshly collected urine is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 72 hours or 2 weeks at -20°C. Samples should initially be used without concentration or dilution, unless the total protein is known to exceed 1000mg/L. If scanning the completed gel by densitometry, an indication of the Urine Total Protein, or the concentration of individual protein zones can be performed by the incorporation of the Urine Protein Internal Standard (Cat. No. 3032) into the sample preparation procedure. Full Instructions For Use are provided with each Urine Protein Internal Standard.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 5.
3. Apply 5µl of sample to each slit and allow to absorb completely. This may take up to 20 minutes.
4. Whilst the samples are absorbing, pour 25ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side down, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 80 volts, 25 minutes.
8. Following electrophoresis fix the gel for 5 minutes in fixative solution.

SAS-MX URINE PROTEIN

9. Dry the gel at 50...60°C. **NOTE:** Cloudiness of the gel, leading to problems with destaining may be observed if the gel is dried at temperatures above 60°C.
10. Immerse the dry gel in stain solution for 10 minutes.
11. Destain the gel in 2 x 1 minute washes of destain solution or until the background is clear.
12. Wash the gel briefly in purified water and dry.

QUALITY CONTROL

A normal serum sample (diluted 1:100 in purified water) can be used to identify the position of each band.

INTERPRETATION OF RESULTS

1. Qualitative Evaluation:

Visual interpretation of bands present on the gel:

Type of Proteinuria	Bands Observed On Gel	Proteins Present
Normal urine	small albumin band	albumin
Glomerular	albumin, alpha-1, beta, gamma	albumin, alpha-1 antitrypsin, transferrin, gamma globulins
Tubular	alpha-1, alpha-2, beta	retinol binding protein, beta2-microglobulin, alpha-2 microglobulin
Overflow	gamma or variable	immunoglobulins, free light chains

2. Quantitative Evaluation:


Scan the gels gel, side down, at 595nm.

In either case, the detection of unusual urine components may require further investigation. The completed SAS-MX Urine Protein gel is stable for an indefinite period of time.

The identification of the major urinary protein outlined can be achieved by immunofixation using the Helena BioSciences SAS-MX Urine IFE Kit (Cat. No. 100600).

LIMITATIONS

Samples with obvious crystalline precipitation present should be centrifuged prior to use.

SAS-MX URINE PROTEIN

PERFORMANCE CHARACTERISTICS

a) Reproducibility

Samples containing 400mg/L and 20mg/L of protein were assayed 5 times each on three gels to determine the within-gel and between-gel reproducibility of elevated and low protein concentrations by densitometry:

Band	Within-Run		Between-Run	
	Mean (%)	CV (%)	Mean	CV (%)
400mg/L	62.8	1.3	62.4	1.6
20mg/L	10.8	7.0	11.9	12.1

b) Sensitivity

25ng per band (equivalent to 5mg/L in the sample), determined as the lowest concentration of protein which was evident as a discrete band on the completed gel.

c) Linearity

Linearity is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.


d) Correlation of Total Protein Quantitation

A range of urine samples were tested on the SAS-MX Urine Protein Gel using the Urine Protein Internal Standard (Cat. No. 3032) to determine the total protein by densitometry. The values obtained (in mg/L) were compared to values obtained for urine total protein using a commercially available assay system based upon a pyrogallol red-molybdate complex assay¹.

BIBLIOGRAPHY

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab., 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Microprotein-PRTM (Procedure No. 611) Sigma Diagnostics, September 1995.

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1. Fauchier, P. and Catalan, F. Interpretive Guide to Clinical Electrophoresis, Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980, Mai/Juni : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Microprotein-PRTM (Procedure No. 611) Sigma Diagnostics, September 1995.

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
BIBLIOGRAFIA

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M. Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980; May/June: 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab. 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Microprotein-PRTM (Procedure No. 61 I) Sigma Diagnostics, September 1995.

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PROTEÍNA DE LA ORINA SAS-MX