helena | BioSciences
Europe
www.helena-biosciences.com

Helena BioSciences Europe
Colima Avenue
Sunderland Enterprise Park
Sunderland
SR5 3XB
tel: +44 (0) 191 549 6064
fax: +44 (0) 191 549 6271
email: info@helena-biosciences.com

Other Helena BioSciences Europe offices:

Helena BioSciences Europe
6 Rue Charles Cros-ZAE
95320 Saint Leu La Foret
France
tel: +33 13 995 9292
fax: +33 13 995 6891
email: helena@helena.fr

Helena BioSciences Europe
Via Enrico Fermi, 24
20090 Assago (Milano)
Italy
tel: +39 02 488 1951
or +39 02 488 2141
fax: +39 02 488 2677

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SAS-MX URINE IFE

INTENDED PURPOSE

The SAS-MX Urine IFE Kit is intended for the identification of urinary proteins by agarose gel electrophoresis followed by immunofixation of the proteins in the gel.

Urinary proteins are derived primarily from plasma proteins that filter through the kidney. The appearance of abnormal plasma proteins in the urine is of great value in evaluating renal function. The appropriate study of proteinuria should include quantitative and qualitative assessment of the type and amount of proteins excreted^{1,5}.

Alfonso first described immunofixation in the literature in 1964⁴. Alper and Johnson published a more practical procedure in 1969, and published a number of studies utilising this technique^{7,8}. Immunofixation has been used as a procedure for the investigation of immunoglobulins since 1976⁹⁻¹¹. The combination of electrophoretic separation of urine proteins, coupled to the identification of specific protein types by immunoprecipitation allows the differentiation of several types of proteinuria - Physiological, Glomerular (selective and non-selective), Tubular and proteinuria associated with Dysglobulinaemias⁵.


The SAS-MX Urine IFE kit separates urine proteins according to charge in an agarose gel. The proteins are then incubated with monospecific antisera, washed and stained to allow visualization of the immunoprecipitate for qualitative interpretation.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

- 1. SAS-MX Urine IFE Gel**
Contains agarose in a barbital buffer with sodium azide as preservative. The gel is ready for use as packaged.
- 2. Tris / Barbital Buffer Concentrate**
Contains concentrated Tris / Barbital buffer with sodium azide as preservative. Dilute the contents of the bottle in 1 litre of purified water and mix well.
- 3. Urine Protein Stain**
Contains Coomassie Blue stain powder. The powdered stain is stable until the expiry date indicated on the label. Dissolve the contents of the vial in 1 litre of 50% methanol and stir overnight. Acidify by adding 20ml of glacial acetic acid. Mix well and filter before use.
- 4. SAS-MX Urine IFE Antisera Kit**
Contains antisera to a) Human Serum/Urine Proteins, b) Urine Micro Proteins (beta-2 microglobulin and retinol binding protein), c) Urine Macro Proteins (Transferrin and alpha-1-antitrypsin), d) IgG,A,M heavy chains, e) Kappa light chains (free and bound) and f) Lambda light chains (free and bound). All antisera contain sodium azide as a preservative. The antisera is ready for use as packaged.
- 5. Other Kit Components**
Each kit contains Instructions For Use and sufficient Sample Application Templates, Antisera Templates and Blotters A, B, C, D and X to complete 10 gels.

SAS-MX URINE IFE

STORAGE AND SHELF-LIFE

1. SAS-MX Urine IFE Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Barbital Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

3. Urine Protein Stain

The powdered stain should be stored at 15...30°C and is stable until the expiry date indicated on the label. Dissolved stain is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration.

4. SAS-MX Urine IFE Antisera Kit

The antisera kit should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 9400 IFE Control Kit

Cat. No. 4062 Incubation Chamber

Cat. No. 9035 Immuno SuperPress / Cat. No. 5014 Development Weight

Drying oven with forced air capable of 60...70°C

Destain solution: Mix 1300ml of purified water and 300ml methanol. Add 400ml of glacial acetic acid.

Store in a tightly stoppered bottle.

Saline Solution (0.85% NaCl)

Purified water

The following items are not required for standard Urine IFE but may be required for further investigation:

Cat. No. 220500 Antiserum to Human Urine Free Kappa Light Chain (2ml)

Cat. No. 220600 Antiserum to Human Urine Free Lambda Light Chain (2ml)

Cat. No. 220900 Antiserum to Human Urine Pentavalent/Albumin (2ml)




Cat. No. 221000 Antiserum to Human Pentavalent (2ml)

SAMPLE COLLECTION AND PREPARATION

Freshly collected urine is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 72 hours or 2 weeks at -20°C. Samples for the Total Lane should have a concentration of 1g/L or less the other lanes should be 5g/L or less. Some samples may require special treatment such as concentration if the total protein is low or dilution if the total protein is high, or if the light chains are poor reacting. See LIMITATIONS for more information.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 5.
3. Apply 3µl of sample to each slit and allow to absorb for 5 minutes.
4. Whilst the samples are absorbing, pour 25ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained in Step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side up, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 120 volts, 40 minutes.
8. Following electrophoresis, place the gel in the incubation chamber containing a wet blotter and position the antiserum application template onto the gel surface. **NOTE:** Ensure good contact of template and gel.
9. Apply 2 drops of the appropriate antiserum / fixative into each lane. Ensure even distribution of antiserum in the lane by rocking the gel.
10. Incubate the gel at 15...30°C for 15 minutes.
11. Following incubation, remove the antiserum template by washing the gel in saline solution for 5 minutes with gentle agitation.
12. Place the gel on a blotter D agarose side up. Place a blotter B (wetted in saline) onto the surface of the gel followed by a blotter X. Press the gel in an IFE SuperPress for 5 minutes or use a development weight for 10 minutes.
13. Remove the blotters and place the gel in saline solution for 10 minutes with gentle agitation.
14. Place the gel on a blotter D agarose side up. Place a blotter B (wetted in saline) onto the surface of the gel followed by a blotter D. Press the gel in an IFE Superpress for 1 minute or use a development weight for 3 minutes.
15. Remove the blotters and dry the gel at 60...70°C.
16. Immerse the dry gel in stain solution for 5 minutes.
17. Destain the gel in 2 x 15 second washes of destain solution.
18. Wash the gel briefly in purified water and dry.

SAS-MX URINE IFE**INTERPRETATION OF RESULTS**

The possible identity of proteins present in the sample can be determined by comparing bands present in the total protein pattern with those present in the other lanes.

Qualitative Evaluation / visual interpretation of bands present on the gel:

Type of Proteinuria	Bands Observed On Gel	Proteins Present
Normal urine	small albumin band	albumin
Glomerular	albumin, alpha-1, beta, gamma	albumin, alpha-1 antitrypsin, transferrin, gamma globulins
Tubular	alpha-1, alpha-2, beta	Retinol Binding Protein, beta-2 microglobulin, alpha-2 microglobulin.
Overflow	gamma or variable	immunoglobulins, free light chains

LIMITATIONS**1. Antigen Excess**

Antigen excess will occur if there is not a slight antibody excess or antigen/antibody equivalence at the site of precipitation. Antigen excess in IFE is usually due to an excess of the immunoglobulin in the patient sample. Antigen excess is characterised by prozoning (unstained areas in the centre of the immunofixed protein band, with staining around the edges). A higher dilution of the sample should be used in this event to optimise the immunoglobulin concentration.

2. Band In Cathodic End of Gamma Region Showing No Reactivity With IFE Antisera

C Reactive Protein (CRP) may be detected in patients with acute inflammatory response^{7,8}. CRP appears as a narrow band at the cathodic end of the serum protein pattern. Elevated Alpha-1-Antitrypsin and Haptoglobin are supportive evidence for CRP. Patients with a CRP band will probably have an elevated level when assayed for CRP.

3. Non-Reactivity With Kappa and Lambda Antisera


Occasionally a sample will have a reaction with a heavy chain antiserum but no light chain reaction is obvious. In this situation, the following need to be ruled out - a) Heavy chain disease, b) Very high concentrations of light chains, leading to antigen excess, c) Low concentrations of light chains, d) Atypical light chain molecule that does not react with the antiserum, e) Light Chains with 'hidden' light chain determinants (as sometimes seen with IgA and IgD). To obtain definitive results, testing may include a) A higher or lower dilution of the sample to optimise the antibody/antigen equivalence, b) Antisera from more than one manufacturer to aid in the identification of atypical immunoglobulins, and c) Treat the sample with β -2-mercaptoethanol to 'reveal' the light chains.

PERFORMANCE CHARACTERISTICS

A series of samples were tested and compared to another commercially available test kit - both kits showed equivalent results.

BIBLIOGRAPHY

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Alfonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
7. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
8. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
9. Johnson, A.M., 'Genetic Typing of Alpha(1)-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of P.M.', J. Lab. Clin. Med., 1976; 87 : 152-163.
10. Cawley, L.P. et al. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
11. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.

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BIBLIOGRAFIA

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M. Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980; May/June: 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab. 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Alfonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
7. Alper, C.A. and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox Sang., 1969; 17 : 445-452.
8. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
9. Johnson, A.M., 'Genetic Typing of α 1-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of P.M.', J. Lab. Clin. Med., 1976; 87 : 152-163.
10. Cawley, L.P. et al. 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
11. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.

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IFE PARA ORINA SAS-MX**INTERPRETACIÓN DE RESULTADOS**