helena | BioSciences
Europe
www.helena-biosciences.com

Helena BioSciences Europe
Colima Avenue
Sunderland Enterprise Park
Sunderland
SR5 3XB
tel: +44 (0) 191 549 6064
fax: +44 (0) 191 549 6271
email: info@helena-biosciences.com

Other Helena BioSciences Europe offices:

Helena BioSciences Europe
6 Rue Charles Cros-ZAE
95320 Saint Leu La Foret
France
tel: +33 13 995 9292
fax: +33 13 995 6891
email: helena@helena.fr

Helena BioSciences Europe
Via Enrico Fermi, 24
20090 Assago (Milano)
Italy
tel: +39 02 488 1951
or +39 02 488 2141
fax: +39 02 488 2677

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SAS-MX LIPOPROTEIN**INTENDED PURPOSE**

The SAS-MX Lipoprotein Kit is intended for the separation and quantitation of lipoproteins in serum or plasma by agarose gel electrophoresis.

Since Fredrickson and Lees proposed a system for phenotyping hyperlipoproteinaemia in 1965¹, the concept of coronary artery disease detection and prevention utilizing lipoprotein electrophoresis has become a relatively common test.

Epidemiological studies have related dietary intake of fats, especially cholesterol and blood levels of the lipids with the incidences of atherosclerosis, major manifestations of which are cardiovascular disease and stroke. Ischemic heart disease has also been related to hypercholesterolaemia^{2,3}. The need for accurate determination of lipoprotein phenotypes resulted from the recognition that hyperlipoproteinaemia is symptomatic of a group of disorders dissimilar in clinical features, prognosis and responsiveness to treatment. Since treatments of the disorders vary with the different phenotypes, it is absolutely necessary that the correct phenotype be established before therapy is begun⁴. In the classification system proposed by Fredrickson and Lees, only types II, III and IV have a proven relationship to atherosclerosis. Plasma lipids do not circulate freely in the plasma, but are transported bound to protein and can thus be classified as lipoproteins. The various fractions are made of different combinations of protein, cholesterol, glycerides, cholesterol esters, phosphatides and free fatty acids⁵. Several techniques have been employed to separate the plasma lipoproteins, including ultracentrifugation, thin layer chromatography, immunological techniques, and electrophoresis. Electrophoresis and ultracentrifugation are two of the most widely used methods and each has given rise to its own terminology. Table I shows the correlation of these classifications and the relative lipid and protein composition of each fraction.


Classification according to:		Composition - % in each fraction			
Electrophoretic Mobility	Ultra-centrifuge	Protein	Glyceride	Cholesterol	Phospholipids
Chylomicrons		2%	98%		
Beta	LDL*	21%	12%	45%	22%
pre-Beta	VLDL*	10%	55%	13%	22%
Alpha	HDL*	50%	6%	18%	26%

*Nonstandard abbreviations: LDL (low density lipoprotein), VLDL (very low density lipoprotein), HDL (high density lipoprotein).

Various exceptions to the above classifications inevitably exist. One of these is the "sinking pre-beta", which is pre-beta migrating material which "sinks" in the ultracentrifuge along with the LDL (beta migrating) fraction⁶. This is the Lp(a) lipoprotein reported by Dahien⁷. It is considered a normal variant found in 10% of the population.

Another exception is the "floating beta", which is beta migrating materials "floating" in the ultracentrifuge with the VLDL.

This abnormal lipoprotein appears in Type III hyperlipoproteinaemias. Various types of support media have been used for the electrophoretic separation of lipoproteins. When Fredrickson originally devised the classification system, he used paper electrophoresis⁸. More recently agarose-gel, starch block and polyacrylamide gel have been used^{5,7}.

The SAS-MX Lipoprotein Kit separates serum / plasma lipoproteins according to charge in agarose gel. The lipoproteins are then fixed and stained for visualisation.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-MX Lipoprotein Gel

Contains agarose in a Tris / Barbitol buffer with sodium azide and thiomersal as preservative. The gel is ready for use as packaged.

2. Tris / Barbitol Buffer Concentrate

Contains concentrated Tris / Barbitol buffer with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well.

3. SAS-MX Lipoprotein Stain

Contains Fat Red 7B stain. Dissolve the contents of the vial in 1 litre of Methanol, stir for 24 hours and filter before use. Preparation of working stain: Immediately prior to use, add 5ml of purified water to 25ml of the stock stain. Add the water drop-wise with stirring.

4. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates and Blotters A and C to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-MX Lipoprotein Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Tris / Barbitol Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C.

3. SAS-MX Lipoprotein Stain

The powdered stain should be stored at 15...30°C and is stable until the expiry date indicated on the label. Dissolved stain is stable for 6 months at 15...30°C. Store in a tightly stoppered bottle.

ITEMS REQUIRED BUT NOT PROVIDED


Cat. No. 4063 Chamber

Cat. No. 1525 EPS600 Power Supply

Drying oven with forced air capable of 60...70°C

Destain solution: mix 75ml of methanol and 25ml of purified water immediately before use.

Purified water

SAS-MX LIPOPROTEIN

SAMPLE COLLECTION AND PREPARATION

Fresh serum or EDTA anticoagulated plasma is the specimen of choice. Samples can be stored at 2...6°C for up to 5 days. DO NOT FREEZE.

Patient Preparation: For the most accurate phenotyping of lipoprotein patterns, the following precautions should be observed before sampling:

- 1) The patient should fast for a 12-14 hour period prior to sampling to prevent interference from meal-induced chylomicrons.
- 2) Discontinue all drugs for 3-4 weeks if possible.
- 3) The patient should be maintaining a stable weight and be on a normal diet for at least 1 week.
- 4) Wait 4-8 weeks after a myocardial infarction or similar traumatic episode.

Interfering Factors:

- 1) Heparin therapy can lead to alterations in the migration of the lipoproteins, particularly beta lipoprotein.
- 2) Samples should not be collected into heparin anticoagulant for similar reasons.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard the blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in step 5.
3. Apply 2µl of sample to each slit and allow to absorb for 7 minutes.
4. Whilst the sample is absorbing, pour 25ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from Step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side up, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 80 volts, 45 minutes.
8. Following electrophoresis, dry the gel at 60...70°C.
9. Place the dry gel in a staining dish and carefully pour the 30ml of freshly prepared working stain onto the gel. Stain for 2 minutes.
10. Destain the gel in 2 x 15-30 seconds washes of destain solution.
11. Wash the gel briefly in purified water and dry.

INTERPRETATION OF RESULTS

It is recommended that any evaluation of the gels is performed against normal values produced for this method in each individual laboratory.

1. **Qualitative Evaluation:** Visually inspect the gels for the presence or absence of particular bands of interest.
2. **Quantitative Evaluation:** Scan the gels gel side down at 525nm.

The alpha-lipoprotein (HDL) is the fastest moving fraction and migrates furthest towards the anode. The beta-lipoprotein (LDL) band is usually the most prominent fraction, migrating closest to the application point. Pre-beta lipoprotein (VLDL) migrates between the alpha and beta lipoproteins. The mobility of the pre-beta lipoprotein varies with the degree of resolution obtained, the type of pre-beta present, and the amount of beta-lipoproteins present. Sometimes, the pre-beta will appear as a smear just in front of the beta-lipoproteins, other times it may split in to 2 separate fractions or may be lacking altogether. The integrity of the pre-beta fraction decreases with sample age. Chylomicrons, when present, remain at the application point.

Calculating the amount of each lipid fraction as mg/dL or mmol/L is not recommended (see LIMITATIONS).

A normal fasting serum can be defined as a clear serum with negligible chylomicrons and normal cholesterol and triglyceride levels. On electrophoresis, the beta-lipoprotein appears as the major fraction with the pre-beta lipoprotein faint or absent and the alpha-lipoprotein band definite but less intense than the beta.

A patient must have an elevated cholesterol or triglycerides to have hyperlipoproteinaemia. The elevation must be determined to be primary or secondary to metabolic disorders such as hypothyroidism, obstructive jaundice, nephrotic syndrome, dysproteinemias, or poorly controlled insulinopaenic diabetes mellitus.

Primary lipidaemia arises from genetically determined factors or environmental factors of unknown mechanism such as diet, alcohol intake and drugs, especially oestrogen or steroid hormones¹². Also considered primary are those lipoproteinaemias associated with ketosis-resistant diabetes, pancreatitis and obesity. Diabetes mellitus and pancreatitis can be confusing, for it is often difficult to tell whether the hyperlipoproteinaemia or the disease is the causative factor.

For a complete review of Lipoprotein phenotyping, with descriptions of the criteria, see Fredrickson, D.S. and Lees, R.S.^{1,8}.

Marked increases in the alpha lipoproteins are seen in obstructive liver disease and cirrhosis. Marked decreases are seen in parenchymal liver disease. Tangier's disease is a rare genetic disorder characterised by the total absence of normal alpha lipoproteins. Heterozygotes exhibit decreased levels of alpha lipoproteins. It should be noted that hyperoestrogenaemia (pregnancy and oral contraceptive use) may cause moderate elevations in the alpha lipoproteins¹².

Abetalipoproteinaemia is a primary inherited defect characterised by severe deficiency of all lipoproteins of density less than 1.063 (all but the alpha lipoproteins). It is accompanied by numerous clinical symptoms and life expectancy is limited. A few cases of familial hypobetalipoproteinaemia have been reported. There is some evidence that the mutation is different from that producing abetalipoproteinaemia.

Lipoprotein-X is an abnormal lipoprotein often seen in patients with obstructive liver disease. It consists of unesterified (free) cholesterol, phospholipids, and VLDL protein. It migrates slower than LDL. Because of its particular lipid content, it stains poorly or not at all with the usual lipid stains and so is not usually detected by standard lipoprotein electrophoresis. Lipoprotein-X is clearly visible when using cholesterol-specific enzymatic staining methods.

SAS-MX LIPOPROTEIN

QUALITY CONTROL

The Lipotrol Control (Cat. No. 5069) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert provided with the control for acceptable assay values.

LIMITATIONS

Fat Red 7B, as well as the Sudan fat stains, has a much greater affinity for triglycerides and cholesterol esters than it has for free cholesterol and phospholipids. Bands seen after staining with these dyes do not reflect a true quantitation of the total plasma lipids¹⁰.

Since the lipid composition of each lipoprotein fraction is variable, it is essential to determine the total cholesterol and triglyceride levels before attempting to classify a pattern^{8,9}. When it comes to diagnosing or ruling out a Type III hyperlipoproteinaemia, a more definitive quantitation of the lipoproteins such as ultracentrifugation¹¹ or PAGE electrophoresis¹¹ is essential.

REFERENCE VALUES

It is recommended that any evaluation of the gels is performed against normal values which have been produced for this test in each individual laboratory.

A normal range study was performed using samples from 48 apparently healthy male and female volunteers:

Fraction	Range
Alpha Lipoproteins	14 - 46%
Pre-Beta Lipoproteins	6 - 40%
Beta Lipoproteins	28 - 61.7%
Chylomicrons	0 - 2%

PERFORMANCE CHARACTERISTICS

Within-Run Precision: 8 replicates of the same sample on a single gel.

Fraction	Mean (%)	CV (%)
Alpha Lipoprotein	26.3	5.2
Pre-Beta Lipoprotein	29.4	4.0
Beta Lipoprotein	44.4	4.3

Between-Run Precision: A single sample run on 10 different gels.

Fraction	Mean (%)	CV (%)
Alpha Lipoprotein	25.2	9.0
Pre-Beta Lipoprotein	34.7	4.2
Beta Lipoprotein	44.0	2.7

BIBLIOGRAPHY

1. Fredrickson, D.S and Lees, R.S. 'A System For Phenotyping Hyperlipoproteinemias', Circulation, 1965; 31(3) : 321-327.
2. Henry, R.J. Ed., 'Clinical Diagnosis and Management of Laboratory Methods', 17th Ed., W.B. Saunders & Co., New York, 194-201, 1984.
3. Lewis, L.A. and Oppet, J.J. Ed., 'CRC Handbook of Electrophoresis Vol II Lipoproteins in Disease', CRC Press Inc., Florida, 63-239, 1980.
4. Levy, R.I. and Fredrickson, D.S. 'Diagnosis and Management of Hyperlipoproteinemia', Am. J. Cardiol., 1968; 22(4) : 576-583.
5. Houstmuller, A.J., 'Agarose-gel Electrophoresis of Lipoproteins: A Clinical Screening Test', Koninklijke Van Gorcum and Comp., The Netherlands, p5, 1969.
6. Stonde, N.J. and Levy, R.I. 'The Hyperlipidemias and Coronary Artery Disease', Disease-A-Month, 1972.
7. Dahlen, G. 'The Pre-Beta Lipoprotein Phenomenon in Relation to Serum Cholesterol and Triglyceride Levels: The Lp(a) Lipoprotein and Coronary Heart Disease', Umea University Medical Dissertations, Sweden, No. 20, 1974.
8. Fredrickson, D.S., Levy, R.I. and Lees, R.S., 'Fat Transport in Lipoproteins - An Integrated Approach To Mechanisms and Disorders' N. Eng. J. Med., 1967; 276 : 34-42, 94-103, 148-156, 215-226, 273-281.
9. Fredrickson, D.S. 'When To Worry About Hyperlipidemia', Consultant, December 1974.
10. Davidsohn, I. And Henry, J.B., Todd-Sanford: 'Clinical Diagnosis by Laboratory Methods', 15th ed., p 639, 1974.
11. Masket, B.H., Levy, R.I and Fredrickson, D.S., 'The Use Of Polyacrylamide Gel Electrophoresis in Differentiating Type III Hyperlipoproteinemia', J. Lab. Clin. Med., 1973; 81(5) : 794-802.
12. World Health Organisation Memorandum: Classification of Hyperlipidemias and Hyperlipoproteinemias', Circulation, 1972; 45 : 501-508.

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BIBLIOGRAPHIE

1. Fredrickson, D.S et Lees, R.S. 'A System For Phenotyping Hyperlipoproteinemias', Circulation, 1965 ; 31(3) : 321-327.
2. Henry, R.J. Ed., 'Clinical Diagnosis and Management of Laboratory Methods', 17^e éd., W.B. Saunders & Co., New York, 194-201, 1984.
3. Lewis, L.A. et Oppet, J.J. Ed., 'CRC Handbook of Electrophoresis Vol II Lipoproteins in Disease', CRC Press Inc., Floride, 63-239, 1980.
4. Levy, R.I. et Fredrickson, D.S. 'Diagnosis and Management of Hyperlipoproteinemia', Am. J. Cardiol., 1968 ; 22(4) : 576-583.
5. Houstmuller, A.J., 'Agarose-gel Electrophoresis of Lipoproteins: A Clinical Screening Test', Koninklijke Van Gorcum and Comp., Hollande, p 5, 1969.
6. Stonde, N.J. et Levy, R.I. 'The Hyperlipidemias and Coronary Artery Disease', Disease-A-Month, 1972.
7. Dahlen, G. 'The Pre-Beta Lipoprotein Phenomenon in Relation to Serum Cholesterol and Triglyceride Levels: The Lp(a) Lipoprotein and Coronary Heart Disease', Umea University Medical Dissertations, Suède, n°. 20, 1974.
8. Fredrickson, D.S., Levy, R.I. et Lees, R.S., 'Fat Transport in Lipoproteins - An Integrated Approach To Mechanisms and Disorders' N. Eng. J. Med., 1967 ; 276 : 34-42, 94-103, 148-156, 215-226, 273-281.
9. Fredrickson, D.S. 'When To Worry About Hyperlipidemia', Consultant, décembre 1974.
10. Davidsohn, I. et Henry, J.B., Todd-Sanford: 'Clinical Diagnosis by Laboratory Methods', 15^e éd., p 639, 1974.
11. Masket, B.H., Levy, R.I et Fredrickson, D.S., 'The Use Of Polyacrylamide Gel Electrophoresis in Differentiating Type III Hyperlipoproteinemia', J. Lab. Clin. Med., 1973 ; 81(5) : 794-802.
12. Mémoire de l'Organisation mondiale de la santé : Classification of Hyperlipidemias and Hyperlipoproteinemias', Circulation, 1972 ; 45 : 501-508.

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LITERATUR

1. Fredrickson, D.S and Lees, R.S. 'A System For Phenotyping Hyperlipoproteinemias', Circulation, 1965; 31(3) : 321-327.
2. Henry, R.J. Ed., 'Clinical Diagnosis and Management of Laboratory Methods', 17th Ed., W.B. Saunders & Co., New York, 194-201, 1984.
3. Lewis, L.A. and Oppet, J.J. Ed., 'CRC Handbook of Electrophoresis Vol II Lipoproteins in Disease', CRC Press Inc., Florida, 63-239, 1980.
4. Levy, R.I. and Fredrickson, D.S. 'Diagnosis and Management of Hyperlipoproteinemia', Am. J. Cardiol., 1968; 22(4) : 576-583.
5. Houstmuller, A.J., 'Agarose-gel Electrophoresis of Lipoproteins: A Clinical Screening Test', Koninklijke Van Gorcum and Comp., The Netherlands, p5, 1969.
6. Stonde, N.J. and Levy, R.I. 'The Hyperlipidemias and Coronary Artery Disease', Disease-A-Month, 1972.
7. Dahlen, G. 'The Pre-Beta Lipoprotein Phenomenon in Relation to Serum Cholesterol and Triglyceride Levels: The Lp(a) Lipoprotein and Coronary Heart Disease', Umea University Medical Dissertations, Sweden, No. 20, 1974.
8. Fredrickson, D.S., Levy, R.I. and Lees, R.S., 'Fat Transport in Lipoproteins - An Integrated Approach To Mechanisms and Disorders' N. Eng. J. Med., 1967; 276 : 34-42, 94-103, 148-156, 215-226, 273-281.
9. Fredrickson, D.S. 'When To Worry About Hyperlipidemia', Consultant, December 1974.
10. Davidsohn, I. And Henry, J.B., Todd-Sanford: 'Clinical Diagnosis by Laboratory Methods', 15th ed., p 639, 1974.
11. Masket, B.H., Levy, R.I and Fredrickson, D.S., 'The Use Of Polyacrylamide Gel Electrophoresis in Differentiating Type III Hyperlipoproteinemia', J. Lab. Clin. Med., 1973 ; 81 (5) : 794-802.
12. World Health Organisation Memorandum: Classification of Hyperlipidemias and Hyperlipoproteinemias', Circulation, 1972; 45 : 501-508.

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1. Fredrickson, D.S and Lees, R.S. 'A System For Phenotyping Hyperlipoproteinemias', *Circulation*, 1965; 31(3) : 321-327.
2. Henry, R.J. Ed., 'Clinical Diagnosis and Management of Laboratory Methods', 17th Ed., W.B. Saunders & Co., New York, 194-201, 1984.
3. Lewis, L.A. and Oppet, J.J. Ed., 'CRC Handbook of Electrophoresis Vol II Lipoproteins in Disease', CRC Press Inc., Florida, 63-239, 1980.
4. Levy, R.I. and Fredrickson, D.S. 'Diagnosis and Management of Hyperlipoproteinemia', *Am. J. Cardiol.*, 1968; 22(4) : 576-583.
5. Houstmuller, A.J., 'Agarose-gel Electrophoresis of Lipoproteins: A Clinical Screening Test', Koninklijke Van Gorcum and Comp., The Netherlands, p5, 1969.
6. Stonde, N.J. and Levy, R.I. 'The Hyperlipidemias and Coronary Artery Disease', *Disease-A-Month*, 1972.
7. Dahlen, G. 'The Pre-Beta Lipoprotein Phenomenon in Relation to Serum Cholesterol and Triglyceride Levels: The Lp(a) Lipoprotein and Coronary Heart Disease', *Umea University Medical Dissertations, Sweden, No. 20, 1974.*
8. Fredrickson, D.S., Levy, R.I. and Lees, R.S., 'Fat Transport in Lipoproteins - An Integrated Approach To Mechanisms and Disorders' *N. Eng. J. Med.*, 1967; 276 : 34-42, 94-103, 148-156, 215-226, 273-281.
9. Fredrickson, D.S. 'When To Worry About Hyperlipidemia', *Consultant*, December 1974.
10. Davidsohn, I. And Henry, J.B., Todd-Sanford: 'Clinical Diagnosis by Laboratory Methods', 15th ed., p 639, 1974.
11. Masket, B.H., Levy, R.I and Fredrickson, D.S., 'The Use Of Polyacrylamide Gel Electrophoresis in Differentiating Type III Hyperlipoproteinemia', *J. Lab. Clin. Med.*, 1973 ; 81 (5) : 794-802.
12. World Health Organisation Memorandum: Classification of Hyperlipidemias and Hyperlipoproteinemias', *Circulation*, 1972; 45 : 501-508.

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LIPOPROTEÍNAS SAS-MX