**helena** | BioSciences  
Europe  
www.helena-biosciences.com

Helena BioSciences Europe  
Colima Avenue  
Sunderland Enterprise Park  
Sunderland  
SR5 3XB  
tel: +44 (0) 191 549 6064  
fax: +44 (0) 191 549 6271  
email: info@helena-biosciences.com

Other Helena BioSciences Europe offices:

Helena BioSciences Europe  
6 Rue Charles Cros-ZAE  
95320 Saint Leu La Foret  
France  
tel: +33 13 995 9292  
fax: +33 13 995 6891  
email: helena@helena.fr

Helena BioSciences Europe  
Via Enrico Fermi, 24  
20090 Assago (Milano)  
Italy  
tel: +39 02 488 1951  
or +39 02 488 2141  
fax: +39 02 488 2677

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## **SAS-MX HIGH RESOLUTION**

### **INTENDED PURPOSE**

The SAS-MX High Resolution Kit is intended for the separation of serum or plasma proteins by agarose gel electrophoresis.

Serum contains over 100 individual proteins, each with a specific set of functions which are subject to specific variation in concentration under different pathological conditions<sup>1</sup>.

Since the introduction of moving boundary electrophoresis by Tiselius<sup>2</sup>, and the subsequent use of zone electrophoresis, serum proteins have been fractionated on the basis of their charge at a particular pH. The SAS-MX High Resolution kit separates serum proteins beyond the classical 5 band pattern into 10-15 discrete fractions, thereby increasing the diagnostic usefulness of the protein patterns<sup>3,5</sup>. Approximately 15 serum proteins have been studied extensively because they may be measured easily<sup>6,9</sup>.

The SAS-MX High Resolution Kit separates serum / plasma proteins according to charge in a buffered agarose gel. The protein bands are then fixed and stained for visualisation and qualitative interpretation.

### **WARNINGS AND PRECAUTIONS**

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

### **COMPOSITION**

#### **1. SAS-MX High Resolution Gel**

Contains agarose in a barbital buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

#### **2. Tris / Barbital Buffer Concentrate**

Contains concentrated Tris / Barbital buffer with sodium azide added as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well.

#### **3. High Resolution Stain Concentrate**

Contains concentrated high resolution stain in a solution containing 5 parts methanol : 5 parts water : 1 part glacial acetic acid. Dilute the contents of the bottle to 700ml with destain solution.

#### **4. Other Kit Components**

Each kit contains Instructions For Use and sufficient Sample Application Templates, Wicks and Blotters A and C to complete 10 gels.

### **STORAGE AND SHELF-LIFE**

#### **1. SAS-MX High Resolution Gel**

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

#### **2. Tris / Barbital Buffer**

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

**3. High Resolution Stain**

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain is stable for 6 months at 15...30°C. Poor staining performance may indicate deterioration of the stain solution.

**ITEMS REQUIRED BUT NOT PROVIDED**

Cat. No. 4063 SAS-MX Chamber  
 Cat. No. 1525 EPS600 Power Supply  
 Cat. No. 3039 Chamber Cooling Device  
 Cat. No. 5141 High Resolution Protein Marker  
 Drying oven with forced air capable of 60...70°C  
 Destain Solution: Mix 500ml of methanol and 500ml of purified water. Add 100ml of glacial acetic acid.  
 Store in a tightly stoppered bottle.  
 Purified water

**SAMPLE COLLECTION AND PREPARATION**

Fresh serum or plasma is the specimen of choice. Samples may be stored covered for a maximum of 48 hours at 2...6°C. The use of plasma will result in a fibrinogen band between the beta and gamma fractions. The use of aged samples will cause the C3 band to migrate in the transferrin region. Serum and Plasma samples should be diluted 1 + 1 with buffer (1 part sample + 1 part buffer). Urine and CSF samples can be used as samples following a suitable concentration step:

Total Protein (g/L)	Concentration Factor
<0.5	100 x
0.5 - 1.0	50 x
1.0 - 3.0	25 x
3.0 - 6.0	10 x
6.0 - 10.0	5 x

Following concentration, urine and CSF samples should be applied directly to the gel.

**STEP-BY-STEP PROCEDURE**

**NOTE:** Ensure that the Chamber Cooling Device has been stored at 2...6°C for at least 60 minutes prior to use.

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 5.
3. Apply 2µl of sample to each slit and allow to absorb for 5 minutes. (A 4µl sample can be applied if darker bands are preferred.)
4. Whilst the samples are absorbing, pour 150ml of buffer into each outer section of the SAS-MX Chamber. Place the cooling device into the central trough of the chamber and wet the surface with a few drops of buffer.
5. Following sample absorption, blot the template with the blotter A retained from Step 2 and remove both blotter and template.

2

**SAS-MX HIGH RESOLUTION**

6. Position the gel on top of the cooling device agarose side up, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber. Take care to avoid trapped air bubbles under the gel.
7. Prepare a wick for each side of the gel by placing 3 wicks together in 2 sets, to make 2 thick wicks. Immerse each wick in the buffer, gently squeeze out excess buffer and attach each wick to the edge of the gel, parallel with the edge of the cooling device. Ensure one edge of the wick is immersed in the buffer and rub one finger across the portion of the wick on the gel to ensure good contact.
8. Electrophorese the gel: 250 volts, 25 minutes.
9. Following electrophoresis, immerse the gel in stain solution for 15 minutes.
10. Destain the gel in destain solution for 30 seconds.
11. Dry the gel at 60...70°C.
12. Destain the gel in 2 x 30 seconds washes of destain solution.
13. Dry the gel.

**ALTERNATIVE STAINING PROCEDURES**

Two alternative staining protocols are available depending on the desired staining requirement of the user.


1. Substitute Amido Black Stain (Cat. No. 3048) for the High Resolution Stain. **FOR SERUM OR PLASMA SAMPLES ONLY.** Do not dilute the sample prior to electrophoresis and follow the staining protocol outlined above. The amido black stain is less sensitive than the coomassie stain, but provides a more uniform staining quality, allowing a better qualitative evaluation of differences in band intensity.
2. Double Staining Technique. **FOR CSF AND URINE SAMPLES ONLY.** After electrophoresis perform amido black staining through to the completed gel then perform the coomassie staining technique through to the completed gel. The double stain technique may allow clearer visualisation of urine protein bands and oligoclonal bands in CSF.

**INTERPRETATION OF RESULTS**

It is recommended that any evaluation of the gel is performed against normal patterns produced for this method in each individual laboratory.

Visually inspect the gel for the presence or absence of particular protein bands:

Figure 1 illustrates the migration patterns of the main plasma proteins which can be identified using the SAS-MX High Resolution Procedure.

3

English

**SAS-MX HIGH RESOLUTION**

High Resolution protein electrophoresis patterns are primarily interpreted by comparing the relative intensities of the bands obtained on unknown samples with those obtained from normal individuals. One of the most common abnormal serum protein patterns is that observed in the non-specific inflammatory response, which is characterized by an increase in  $\alpha$ 1-antitrypsin and haptoglobin with decreased prealbumin, albumin and transferrin. While it is not useful in establishing a general diagnosis, it is useful in monitoring a patient's response to therapy.

Other examples of clinically important variations are:

1. Elevation of the transferrin band, suggesting low levels of iron.
2. Presence of monoclonal proteins, suggesting abnormalities in the immune system.
3. Low haptoglobin, suggesting elevated red cell turnover, or in-vitro haemolysis.
4. CRP presence, indicating an acute inflammatory response.
5. Low prealbumin, albumin and transferrin with diffuse hypergammaglobulinaemia, suggesting chronic inflammation, infection or antigenic stimulation.
6. Low C3 on fresh samples, suggesting complement consumption.

High Resolution protein electrophoresis is an excellent analytical tool to gain a broad overview of urine proteins<sup>3,10</sup>. Normal urine normally contains a trace of albumin, and sometimes a faint transferrin band. Glomerular-type proteinuria usually consists of strong bands of albumin, both  $\alpha$ 1-acid glycoprotein and  $\alpha$ 1-antitrypsin in a broad  $\alpha$ 1 zone and transferrin in the  $\beta$ 1 region. The serum pattern shows decreases in these proteins, with increases in the proteins retained by the glomerulus. The urine pattern in Tubular proteinuria usually consists of a faint albumin band, a double band in the  $\alpha$ 2 region due to  $\alpha$ 2-microglobulin, a strong band in the mid-beta region due to  $\beta$ 2-microglobulin and sometimes diffuse background staining in the gamma region due to free light chains. Chronic renal disease, or renal failure can lead to a mixed type pattern, caused by damage to both the tubules and glomerulus.

An elevation or decrease in particular sample components or the detection of unusual sample components require further investigation. The completed SAS-MX High Resolution gel is stable for an indefinite period of time.

**BIBLIOGRAPHY**

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Laurell, C.B. 'Composition and Variation of the Gel Electrophoretic Fractions of Plasma, Cerebrospinal Fluid and Urine' Scand. J. Clin. Lab. Invest., 1972; 29 (Suppl. 24) : 71.
4. Killingsworth, L.M. et al. 'Protein Analysis, Deciphering Cerebrospinal Fluid Patterns' Diag. Med., 1980; March/April : 1-7.
5. Killingsworth, L.M. et al. 'Protein Analysis, Finding Clues to Disease in Urine' Diag. Med., 1980; May/June : 69-75.
6. Killingsworth, L.M. 'Clinical Applications of Protein Determinations in Biological Fluids Other Than Blood' Clin. Chem., 1982; 28(5) : 1093-1102.
7. Ritzmann, S.E and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins. Qualitative and Quantitative Assays' in Laboratory Medicine, Harper & Row, Inc., Hagerstown, 1979.
8. Killingsworth, L.M. et al. 'Protein Analysis' Diag. Med., 1980; Jan/Feb : 3-15.

9. Killingsworth, L.M. 'Plasma Protein Patterns in Health and Disease' CRC Critical Reviews in Clinical Laboratory Sciences, August 1979.
10. Peterson, P.A., Ervin, P.E. and Beggard, I. 'Differentiation of Glomerular, Tubular and Normal Proteinuria: determination of Urinary Excretion of  $\alpha$ 2-microglobulin, Albumin and Total Protein' J. Clin. Invest., 1969; 48 : 1189-1198.

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#### BIBLIOGRAPHIE

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974 ; 291 : 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937 ; 33 : 524.
3. Laurell, C.B. 'Composition and Variation of the Gel Electrophoretic Fractions of Plasma, Cerebrospinal Fluid and Urine' Scand. J. Clin. Lab. Invest., 1972 ; 29 (suppl. 24) : 71.
4. Killingsworth, L.M. et al. 'Protein Analysis, Deciphering Cerebrospinal Fluid Patterns' Diag. Med., 1980; mars/avril : 1-7.
5. Killingsworth, L.M. et al. 'Protein Analysis, Finding Clues to Disease in Urine' Diag. Med., 1980 ; mai/juin : 69-75.
6. Killingsworth, L.M. 'Clinical Applications of Protein Determinations in Biological Fluids Other Than Blood' Clin. Chem., 1982 ; 28(5) : 1093-1102.
7. Ritzmann, S.E et Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins. Qualitative and Quantitative Assays' in Laboratory Medicine, Harper & Row, Inc., Hagerstown, 1979.
8. Killingsworth, L.M. et al. 'Protein Analysis' Diag. Med., 1980; jan/fév : 3-15.
9. Killingsworth, L.M. 'Plasma Protein Patterns in Health and Disease' CRC Critical Reviews in Clinical Laboratory Sciences, août 1979.
10. Peterson, P.A., Ervin, P.E. et Beggard, I. 'Differentiation of Glomerular, Tubular and Normal Proteinuria: determination of Urinary Excretion of ?2-microglobulin, Albumin and Total Protein' J. Clin. Invest., 1969; 48 : 1189-1198.

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**LITERATUR**

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Laurell, C.B. 'Composition and Variation of the Gel Electrophoretic Fractions of Plasma, Cerebrospinal Fluid and Urine' Scand. J. Clin. Lab. Invest., 1972; 29 (Suppl. 24) : 71.
4. Killingsworth, L.M. et al. 'Protein Analysis, Deciphering Cerebrospinal Fluid Patterns' Diag. Med., 1980; March/April : 1-7.
5. Killingsworth, L.M. et al. 'Protein Analysis, Finding Clues to Disease in Urine' Diag. Med., 1980; May/June : 69-75.
6. Killingsworth, L.M. 'Clinical Applications of Protein Determinations in Biological Fluids Other Than Blood' Clin. Chem., 1982; 28(5) : 1093-1102.
7. Ritzmann, S.E and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins. Qualitative and Quantitative Assays' in Laboratory Medicine, Harper & Row, Inc., Hagerstown, 1979.
8. Killingsworth, L.M. et al. 'Protein Analysis' Diag. Med., 1980; Jan/Feb : 3-15.
9. Killingsworth, L.M. 'Plasma Protein Patterns in Health and Disease' CRC Critical Reviews in Clinical Laboratory Sciences, August 1979.
10. Peterson, P.A., Ervin, P.E. and Beggard, I. 'Differentiation of Glomerular, Tubular and Normal Proteinuria: determination of Urinary Excretion of  $\alpha_2$ -microglobulin, Albumin and Total Protein' J. Clin. Invest., 1969; 48 : 1189-1198.

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#### BIBLIOGRAFIA

1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Laurell, C.B. 'Composition and Variation of the Gel Electrophoretic Fractions of Plasma, Cerebrospinal Fluid and Urine' Scand. J. Clin. Lab. Invest., 1972; 29 (Suppl. 24) : 71.
4. Killingsworth, L.M. et al. 'Protein Analysis, Deciphering Cerebrospinal Fluid Patterns' Diag. Med., 1980; March/April : 1-7.
5. Killingsworth, L.M. et al. 'Protein Analysis, Finding Clues to Disease in Urine' Diag. Med., 1980; May/June : 69-75.
6. Killingsworth, L.M. 'Clinical Applications of Protein Determinations in Biological Fluids Other Than Blood' Clin. Chem., 1982; 28(5) : 1093-1102.
7. Ritzmann, S.E. and Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins. Qualitative and Quantitative Assays' in Laboratory Medicine, Harper & Row, Inc., Hagerstown, 1979.
8. Killingsworth, L.M. et al. 'Protein Analysis' Diag. Med., 1980; Jan/Feb : 3-15.
9. Killingsworth, L.M. 'Plasma Protein Patterns in Health and Disease' CRC Critical Reviews in Clinical Laboratory Sciences, August 1979.
10. Peterson, P.A., Ervin, P.E. and Beggard, I. 'Differentiation of Glomerular, Tubular and Normal Proteinuria: determination of Urinary Excretion of  $\alpha$ 2-microglobulin, Albumin and Total Protein' J Clin. Invest., 1969; 48 : 1189-1198.

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1. Alper, C.A. 'Plasma Protein Measurements as a Diagnostic Aid', N. Eng. J. Med., 1974; 291 : 287-290.
2. Tiselius, A. 'A New Approach for Electrophoretic Analysis of Colloidal Mixtures', Trans. Faraday Soc., 1937; 33 : 524.
3. Laurell, C.B. 'Composition and Variation of the Gel Electrophoretic Fractions of Plasma, Cerebrospinal Fluid and Urine' Scand. J. Clin. Lab. Invest., 1972; 29 (Suppl. 24) : 71.
4. Killingsworth, L.M. et alter. 'Protein Analysis, Deciphering Cerebrospinal Fluid Patterns' Diag. Med., 1980; Marzo/Abril : 1-7.
5. Killingsworth, L.M. et alter. 'Protein Analysis, Finding Clues to Disease in Urine' Diag. Med., 1980; Mayo/Junio: 69-75.
6. Killingsworth, L.M. 'Clinical Applications of Protein Determinations in Biological Fluids Other Than Blood' Clin. Chem., 1982; 28(5) : 1093-1102.
7. Ritzmann, S.E y Daniels, J.C. 'Diagnostic Proteinology: Separation and Characterization of Proteins. Qualitative and Quantitative Assays' en Laboratory Medicine, Harper & Row, Inc., Hagerstown, 1979.
8. Killingsworth, L.M. et alter. 'Protein Analysis' Diag. Med., 1980; Ene/Feb: 3-15.
9. Killingsworth, L.M. 'Plasma Protein Patterns in Health and Disease' CRC Critical Reviews in Clinical Laboratory Sciences, agosto 1979.
10. Peterson, P.A., Ervin, P.E. y Beggard, I. 'Differentiation of Glomerular, Tubular and Normal Proteinuria: determination of Urinary Excretion of  $\alpha_2$ -microglobulin, Albumin and Total Protein' J. Clin. Invest., 1969; 48 : 1189-1198.

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### ALTA RESOLUCIÓN SAS-MX