SAS-MX CK VIS-10 Isoenzyme

Instructions For Use

REF 101800

REF 101801 (gel kit only) REF 101802 (reagent kit only)

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
INTENDED PURPOSE

The SAS-MX CK Vis-10 Isoenzyme Kit is intended for the separation and quantitation of Creatine Kinase Isoenzymes in serum or plasma by agarose gel electrophoresis.

Creatine phosphokinase (CK) (EC 2.7.3.2.) is an energy transfer enzyme, which catalyses the reversible reaction


CK exists primarily in skeletal muscle, cardiac muscle and the brain, with small amounts in several other tissues. A number of diverse clinical episodes such as surgical procedures, intramuscular injections and myocardial infarct induce increased CK activity in the serum^{2,3}. The source of elevated CK activity may be narrowed by isoenzyme assessment. There are two molecular CK subunits, designated M and B, the combinations of which produce three isoenzymes. CK-MM (isolated primarily from skeletal muscle), CK-MB (myocardium) and CK-BB (primarily from the brain).

CK isoenzyme analysis is one of the most important procedures used in the early detection of myocardial damage⁴. After an acute myocardial infarction (MI), CK-MB appears in the serum in approximately 4 to 6 hours, reaches peak activity at 18-24 hours, and may disappear completely within 72 hours. Within the first 48 hours after MI, CK-MB is present in 100% of the patients with MI, as well as in some cases of severe coronary insufficiency^{1,3,5}.

The most definitive laboratory testing in the diagnosis of MI is accomplished by performing studies of CK isoenzyme in conjunction with lactate dehydrogenase (LD) isoenzymes^{3,5,6}. The specificity and sensitivity achieved with these two tests has eliminated the necessity for additional enzyme studies in accurately diagnosing MI⁷. The most important consideration in the interpretation of CK and LD isoenzyme patterns is the detection of the characteristic change of pattern of multiple examinations (the relatively fast appearance and disappearance of CK-MB and the flip of LD1 over LD2)^{1,3,9}. Persistent elevation in CK-MB is not indicative of myocardial infarct. CK-MB may be useful in diagnosing a small infarct in which total CK never exceeds the upper limit of normal¹⁰.

CK produced by myocardium is only 25-40% CK-MB, the remainder being CK-MM. Therefore, an elevation in CK due to myocardial infarction produces not only a rise in CK-MB but in CK-MM as well. The isoenzymes of CK have been assessed by various methods^{11,20}. Electrophoresis offers the distinct advantage of complete separation of the isoenzymes without risk of carryover³.

The isoenzymes of CK are separated according to their electrophoretic mobility on an agarose gel, followed by incubation with a reagent which visualizes the CK isoenzymes according to the following reaction sequence:

The patterns may be qualitatively examined or the patterns may be quantitated by densitometry.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-MX CK Isoenzyme Gel (x10)

Contains agarose in a AMP buffer with sodium azide as preservative. The gel is ready for use as packaged.

2. SAS-MX CK Isoenzyme Buffer (1x 100ml)

Contains MOPSO buffer pH 7.6 - 8.0. Dilute the contents of the bottle to 1 litre with purified water.

3. CK Isoenzyme Reagent (10x 1ml)

Contains those components required for the identification of creatine kinase isoenzymes according to the reaction sequence shown above. See STEP-BY-STEP PROCEDURE for reconstitution instructions.

4. CK Isoenzyme Diluent (1x 10.5ml)

Contains MES, sucrose and sodium azide as preservative. The diluent is ready for use as packaged.

5. CK Isoenzyme Chromogen (1x 1ml)

Contains tetranitro blue tetrazolium (TNBT) and dimethylformamide (DMF). The chromogen is ready for use as packaged.

6. CK Isoenzyme Activator (1x 1ml)

Contains 2-Mercaptoethanol in water. The activator is ready for use as packaged.

7. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates (Template 2 = 10 x 8 sample, Template 4 = 10 x 10 sample), and Blotters A and D to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-MX CK Isoenzyme Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. SAS-MX CK Isoenzyme Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C.

3. CK Isoenzyme Reagent

The reagent vials should be stored at 2...6°C and is stable until the expiry date indicated on the label.

4. CK Isoenzyme Diluent

The diluent should be stored at 2...6°C and is stable until the expiry date indicated on the label.

5. CK Isoenzyme Chromogen

The chromogen should be stored at 2...6°C and is stable until the expiry date indicated on the label.

6. CK Isoenzyme Activator

The activator should be stored at 2...6°C and is stable until the expiry date indicated on the label. Avoid contamination.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 4062 Incubation Chamber

Drying Oven with forced air capable of 60...70°C

Incubator capable of 45°C

Glass Plate

Destain Solution: Mix 100ml of glacial acetic acid and 900ml of purified water. Store in a tightly stoppered bottle.

Purified water

SAMPLE COLLECTION AND PREPARATION

Fresh serum is the specimen of choice. Samples should be tested as soon as possible after collection. If required, samples can be stored refrigerated at 2...6°C for up to 48 hours²³ or frozen at -20°C for up to 2 weeks²³. Allow to thaw at room temperature.

Sample Collection: Proper timing of sample collection is critical for the accurate interpretation of CK Isoenzyme patterns in the assessment of myocardial infarction (MI). A minimum of 4 samples should be collected - one on admission, followed by samples every 8-12 hours for a minimum of 36 hours.

Interfering factors:

- 1) Whilst red cells contain no CK activity, the products of haemolysis can interfere with the reaction sequence for the enzyme, so haemolysed samples are not recommended for use²¹.
- 2) Repeated freezing and thawing of samples destroys CK activity.
- 3) Certain drugs can affect CK isoenzyme activity²².
- 4) CK is inactivated by heat²¹.
- 5) In some samples, albumin may produce a visible band between MB and BB band regions.

STEP-BY-STEP PROCEDURE

1. To every 100µl of patient sample or control add 1µl of CK Isoenzyme Activator and incubate at 15...30°C for 10 minutes before use.
2. Remove the gel from the packaging and place on a paper towel. Blot the gel surface at the application point, (between the '-' marks at the edges of the gel) with a blotter A, discard the blotter. **NOTE:** From the time the gel is taken out of the plastic shell, no more than 5-7 minutes should elapse before electrophoresis is started.
3. Align the sample application template (either 8 sample or 10 sample according to preference) with the '-' marks at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in step 6.
4. Apply 1µl of sample to each slit and allow to absorb for 4 minutes.
5. Whilst the samples are absorbing, pour 45ml of buffer into each inner section of the SAS-MX Chamber.
6. Following sample absorption, blot the template with the blotter A retained from Step 3 and remove both blotter and template.
7. Position the gel in the chamber agarose side down, aligning the positive (+) and the negative (-) sides with the corresponding positions on the chamber.
8. Electrophorese the gel: 100 volts, 20 minutes.
9. Whilst the samples are electrophoresing, place a blotter D moistened with water into the incubation chamber and place in the incubator set at 45°C to equilibrate.

10. Approximately 3-4 minutes before the end of electrophoresis, reconstitute 1 vial of CK Isoenzyme Reagent by adding 1 ml of CK Isoenzyme Diluent and 70 μ l of CK Isoenzyme Chromogen and mix well. Do not vortex.
11. Following electrophoresis, place the gel agarose side up on the glass plate. Pour the contents of the CK Isoenzyme Reagent vial along the anode (+) side of the gel.
12. Using a serological pipette, spread the reagent from the anode edge to the cathode (-) edge of the gel. Wait 15 seconds and spread the reagent from cathode to anode edge. Wait 15 seconds and roll the excess off the gel.
13. Remove the gel from the glass plate and place in the incubation chamber. Incubate the gel: 30 minutes, 45°C.
14. Immerse the gel in destain solution for 10 minutes.
15. Wash the gel in purified water and dry at 60...70°C.

INTERPRETATION OF RESULTS

It is recommended that any evaluation of the gels is performed against normal values which have been produced for this test in each individual laboratory. Completed gels should be stored in the dark at all times.

Qualitative Evaluation: The SAS-MX CK Isoenzyme gel can be visually inspected for the presence or absence of bands of interest.

Quantitative Evaluation: Scan the gels at 590nm within 1 hour of completion.

CK-BB is the fastest moving, most anodic band and CK-MM is the slowest moving most cathodic band. CK-MB migrates intermediate to CK-MM and CK-BB¹. A number of atypical CK bands have been reported. Atypical bands migrating between CK-MB have been attributed to CK-BB complexed to IgG and CK-MM complexed to lipoprotein as well as others without positive identification. Mitochondrial CK migrates cathodally to CK-MM and a band designated -macro- CK isolated from a cancer patient also migrated cathodally to CK-MM.

CK-MM

Often the only isoenzyme of CK found in normal serum¹.

Elevated in skeletal muscle injury, myocardial injury or brain injury^{1,3}.

CK-MB

May be present in serum from normal subjects in the amount of 0-4%²⁴. Small amounts of CK-MB activity have been interpreted as an alert to possible myocardial infarct and should be followed by serial CK and LDH Isoenzyme studies.

Positive indication of myocardial infarct when the following criteria are met: The proper clinical setting², CK-MB activity >5% of the total CK activity and a minimum of 10IU/L^{1,15,25}, CK-MB shows the characteristic change in pattern (relatively rapid appearance and disappearance)^{1,3,9}.

Positive identification of second myocardial infarct: after the first MI the CK-MB increases after starting to decline. The total CK may not show an increase after starting to decline.

Values following open heart surgery³ - CK and LD isoenzymes are less specific following open heart surgery than in most diagnostic situations. The CK-MB will be elevated due to myocardial damage from the operative procedure as well as trauma to the heart from manipulation and cannulation. The LD is flipped secondary to haemolysis from extracorporeal circulation. Infarct patients have higher levels of CK-MB activity, but the wide range of isoenzyme activity seen in non-MI patients overlaps that noted

in patients with MI, making complete discrimination impossible. Despite this difficulty, accuracy in diagnosing MI can be increased by doing serial determinations of CK-MB in the post-operative period and monitoring the trends. Peri-operative infarct patients will usually have a progressive rise in CK-MB levels, while non-MI patients exhibit a more rapid decrease in CK-MB post-operative^{2,26}. Elevation in diseases other than MI: Severe coronary insufficiency, Duchenne's muscular dystrophy, Rocky Mountain Spotted Fever, Rhabdomyolysis, Dermatomyositis, Myoglobinuria, Polymyositis, Reye's Syndrome.

CK-BB

Often seen in the serum of patients with prostatic carcinoma and occasionally in the serum of patients with other carcinomas and malignant tumours¹ (See LIMITATIONS). Rarely seen in the serum of patients with brain injury due to damage to the blood-brain barrier^{1,27}. Occasionally seen in the serum of patients with severe shock syndrome (probably due to lung or small bowel involvement). Occasionally seen in the serum of patients with chronic renal failure, women in labour¹.

A number of atypical CK bands have been reported. Atypical bands migrating between CK-MB and CK-MM have been attributed to CK-BB complexed to IgG^{28,29} and CK-MM complexed to lipoprotein³⁰, as well as others without positive identification³¹⁻³³. Mitochondrial CK migrates cathodal to CK-MM³⁴, and a band designated 'macro' CK, isolated from a cancer patient also migrated cathodal to CK-MM³⁵.

QUALITY CONTROL

The CK/LD Control (Cat.No.5134) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert provided for acceptable assay values.

LIMITATIONS

The SAS-MX CK Isoenzyme procedure is not designed to identify tumour markers.

REFERENCE VALUES

Normal range studies performed using apparently healthy male and female volunteers gave the following results:

CK-MM = 100%

CK-MB = 0%

CK-BB = 0%

PERFORMANCE CHARACTERISTICS**Linearity**

The method gives a linear response to CK activities up to 1300 IU/L per band.

Sensitivity

The method is sensitive to 8 IU/L per band.

Reproducibility

	Within Gel (n=8)		Between Gel (n=4)	
	Mean (%)	CV (%)	Mean (%)	CV (%)
CK-MM	76.5	0.5	76.4	0.8
CK-MB	11.1	1.8	11.2	4.5
CB-BB	12.5	2.7	12.5	2.5

BIBLIOGRAPHY

1. Brish, L.K., 'CK & LD Isoenzyme. A Self-Instructional Text', Am. Soc. Of Clin. Path. Press, Chicago, 1984 : 30-47.
2. Wolf, P.L., Griffiths, J.C and Koett, J.W., 'Interpretation of Electrophoretic Patterns of Proteins and Isoenzymes', Mason Publishers, New York, 1981 : 60.
3. Tietz, N.W., Ed. 'Textbook of Clinical Chemistry', W.B. Saunders Co., Philadelphia, 1986 : 678-700.
4. Galen, R.S., 'The Enzyme Diagnosis of Myocardial Infarction', Human Path., 1975; 6(2): 141-155.
5. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
6. Marmor, A. and Alpan, G., 'Specificity of Creatine Kinase MB Isoenzyme for Myocardial Injury', Clin. Chem., 1978; 24(12): 2206.
7. Galen, R.S., Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
8. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978; 11(6): 232-234.
9. Yasmineh, W.G., Pyle, R.B, Cohn, J.N., Nicoloff, D.M., Hanson, N.Q. and Steele, B.W., 'Serial Serum Creatine Phosphokinase MB Isoenzyme Activity After Myocardial Infarction', Circulation, 1977; 55(5): 733-737.
10. D'Souza, J.P., Sine, H.E., Horvitz, R.A., Kubasik, N.P., Brody, B.B. and Barold, S.S., 'The Significance of the MB Isoenzyme in Patients with Acute Cardiovascular Disease with a Normal or Borderline Total CPK Activity', Clin. Biochem., 1978; 11(5): 204-209.
11. Mercer, D.W., 'Separation of Tissue and Serum Creatine Kinase Isoenzymes by Ion-Exchange Column Chromatography', Clin. Chem., 1974; 20(1): 36-40.
12. Rao, P.S., Lukes, J.J., Ayres, S.M. and Mueller, H., 'New manual and Automated Method for Determining Activity of Creatine Kinase Isoenzyme MB, by Use of Dithiothreitol: Clinical Applications', Clin. Chem., 1975; 21(11): 1612-1618.
13. Jockers-Wretou, E. and Pfeleiderer, G., 'Quantitation of Creatine Kinase Isoenzymes in Human Tissues and Sera by an Immunological Method', Clin. Chim. Acta, 1975; 58: 223-232.
14. Roberts, R., Sobel, B.E. and Parker, C.W., 'Radioimmunoassay for Creatine Kinase Isoenzymes', Science, 1976; 194: 855-857.
15. Van Der Veen, K.J. and Willebrands, A.F., 'Isoenzymes of Creatine Phosphokinase in Tissue Extracts and in Normal and Pathological Sera', Clin. Chim. Acta, 1966; 13: 312-316.
16. Wong, R. and Swallen, T.O., 'Cellulose Acetate Electrophoresis of Creatine Phosphokinase Isoenzymes in the Diagnosis of Myocardial Infarction', Am. J. Clin. Pathol., 1975; 64: 209-216.
17. Burger, A., Richterich, R. und Aebi, H., 'Die Heterogenitat der Kreatin-Kinase', Biochemische Zeitschrift, 1964; 339: 305-314.
18. Deul, D.H. and van Breeman, J.F.L., 'Electrophoresis of Creatine Phosphokinase from Various Organs', Clin. Chim. Acta, 1964; 10: 276-283.
19. Rosalki, S.B., 'Creatine Phosphokinase Isoenzymes', Nature, 1965; 207: 414.
20. Trainer, T.D. and Gruenig, D., 'A Rapid Method for the Analysis of Creatine Phosphokinase Isoenzymes', Clin. Chim. Acta, 1968; 21: 151-154.
21. Henry, R.J. and Cannon, D.C., Ed's, 'Clinical Chemistry Principles and Techniques', 2nd Edition, Harper and Row, New York, 1974: 901-903.

22. Young, D.S. , 'Effects of Drugs on Clinical Laboratory Tests', 3rd Edition, AACC Press, Washington D.C., 1990.
23. Cho, H.W. and Meltzer, H.Y., 'Factors Affecting Stability of Isoenzymes of Creatine Phosphokinase', *Am. J. Clin. Pathol.*, 1979; 71(1): 75-82.
24. Marmor, A., Alpan, G., Keidar, S., Grenadier, E. and Palant, A., 'The MB Isoenzyme of Creatine Kinase as an Indicator of Severity of Myocardial Ischaemia', *Lancet*, 1978; Oct 14: 812-814.
25. Galen, R.S., Personal Communication, June 1980.
26. Galen, R.S., 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
27. Bayer, P.M., Gabl, F., Granditsch, G., Widhalm, K., Zyman, H. and Deutsch, E., 'Creatine Kinase Isoenzymes in Cerebral Fluid in a Case of Brain Damage', *Clin. Chem.*, 1976; 22(8): 1405-1407.
28. Stein, W. and Bohner, 'Immunoglobulin-Bound Creatine Kinase BB ("Macro CK") in Three Patients with Different Diseases', *Clin. Chem.*, 1979; 25(8): 1513-1514.
29. Pudek, M.R., Jacobson, B.E., Urquhart, N. and Rabkin, S.W., 'Effect of Macro Creatine Kinase BB on Results of Different Creatine Kinase Isoenzyme Methods', *Clin. Chem.*, 1982; 28(6): 1400.
30. Velletri, K., Griffiths, W.C. and Diamond, I., 'Abnormal Electrophoretic Mobility of a Creatine Kinase MM Isoenzyme', *Clin. Chem.*, 1975; 21(12): 1837-1838.
31. Sax, S.M., Moore, J.J., Giegel, J.L. and Welsh, M., 'Atypical Increase in Serum Creatine Kinase Activity in Hospital Patients', *Clin. Chem.*, 1976; 22(1): 87-91.
32. Ljungdahl, L. and Gerhardt, W., 'Creatine Kinase Isoenzyme Variants in Human Serum', *Clin. Chem.*, 1978; 24(5): 832-834.
33. Lott, J.A. and Heinz, J.W., 'A Stable, Unusual Isoenzyme of Creatine Kinase Observed in a Case of Acute Anoxia', *Clin. Chem.*, 1978; 24(6): 1047.
34. James, G.P. and Harrison, R.L., 'Creatine Kinase Isoenzymes of Mitochondrial Origin in Human Serum', *Clin. Chem.*, 1979; 25(6): 943-947.
35. Yuu, Hoyuo, et al., *Clin. Chem.*, 1978; 24(11): 2054-2057.

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LITERATUR

1. Brish, L. K., 'CK & LD Isoenzyme. A Self-Instructional Text', Am. Soc. Of Clin. Path. Press, Chicago, 1984 : 30-47.
2. Wolf, P.L., Griffiths, J.C and Koett, J.W., 'Interpretation of Electrophoretic Patterns of Proteins and Isoenzymes', Mason Publishers, New York, 1981 : 60.
3. Tietz, N.W., Ed. 'Textbook of Clinical Chemistry', W.B. Saunders Co., Philadelphia, 1986 : 678-700.
4. Galen, R. S., 'The Enzyme Diagnosis of Myocardial Infarction', Human Path., 1975; 6(2): 141-155.
5. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
6. Marmor, A. and Alpan, G., 'Specificity of Creatine Kinase MB Isoenzyme for Myocardial Injury', Clin. Chem., 1978; 24(12) : 2206.
7. Galen, R.S., Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
8. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978; 11(6): 232-234.
9. Yasmineh, W.G., Pyle, R.B, Cohn, J.N., Nicoloff, D.M., Hanson, N.Q. and Steele, B.W., 'Serial Serum Creatine Phosphokinase MB Isoenzyme Activity After Myocardial Infarction', Circulation, 1977; 55(5): 733-737.
10. D'Souza, J.P., Sine, H.E., Horvitz, R.A., Kubasik, N.P., Brody, B.B. and Barold, S.S., 'The Significance of the MB Isoenzyme in Patients with Acute Cardiovascular Disease with a Normal or Borderline Total CPK Activity', Clin. Biochem., 1978; 11(5): 204-209.
11. Mercer, D.W., 'Separation of Tissue and Serum Creatine Kinase Isoenzymes by Ion-Exchange Column Chromatography', Clin. Chem., 1974; 20(1) : 36-40.
12. Rao, P.S., Lukes, J.J., Ayres, S.M. and Mueller, H., 'New manual and Automated Method for Determining Activity of Creatine Kinase Isoenzyme MB, by Use of Dithiothreitol: Clinical Applications', Clin. Chem., 1975; 21(11) : 1612-1618.
13. Jockers-Wretou, E. and Pfeleiderer, G., 'Quantitation of Creatine Kinase Isoenzymes in Human Tissues and Sera by an Immunological Method', Clin. Chim. Acta, 1975; 58: 223-232.
14. Roberts, R., Sobel, B.E. and Parker, C.W., 'Radioimmunoassay for Creatine Kinase Isoenzymes', Science, 1976; 194: 855-857.
15. Van Der Veen, K.J. and Willebrands, A.F., 'Isoenzymes of Creatine Phosphokinase in Tissue Extracts and in Normal and Pathological Sera', Clin. Chim. Acta, 1966; 13: 312-316.
16. Wong, R. and Swallen, T.O., 'Cellulose Acetate Electrophoresis of Creatine Phosphokinase Isoenzymes in the Diagnosis of Myocardial Infarction', Am. J. Clin. Pathol., 1975; 64: 209-216.
17. Burger, A., Richterich, R. und Aebi, H., 'Die Heterogenitat der Kreatin-Kinase', Biochemische Zeitschrift, 1964; 339: 305-314.
18. Deul, D.H. and van Breeman, J.F.L., 'Electrophoresis of Creatine Phosphokinase from Various Organs', Clin. Chim. Acta, 1964; 10: 276-283.
19. Rosalki, S.B., 'Creatine Phosphokinase Isoenzymes', Nature, 1965; 207: 414.
20. Trainer, T.D. and Gruenig, D., 'A Rapid Method for the Analysis of Creatine Phosphokinase Isoenzymes', Clin. Chim. Acta, 1968; 21: 151-154.
21. Henry, R.J. and Cannon, D.C., Ed's, 'Clinical Chemistry Principles and Techniques', 2nd Edition, Harper and Row, New York, 1974: 901-903.

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BIBLIOGRAFIA

1. Brish, L.K., 'CK & LD Isoenzyme. A Self-Instructional Text', Am. Soc. Of Clin. Path. Press, Chicago, 1984 : 30-47.
2. Wolf, P.L., Griffiths, J.C and Koett, J.W., 'Interpretation of Electrophoretic Patterns of Proteins and Isoenzymes', Mason Publishers, New York, 1981 : 60.
3. Tietz, N.W., Ed. 'Textbook of Clinical Chemistry', W.B. Saunders Co., Philadelphia, 1986. 678-700.
4. Galen, R.S. 'The Enzyme Diagnosis of Myocardial Infarction', Human Path., 1975; 6(2): 141-155.
5. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
6. Marmor, A. and Alpan, G., 'Specificity of Creatine Kinase MB Isoenzyme for Myocardial Injury', Clin. Chem., 1978; 24(12): 2206.
7. Galen, R.S. Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
8. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978; 11(6): 232-234.
9. Yasmineh, W.G., Pyle, R.B, Cohn, J.N., Nicoloff, D.M., Hanson, N.Q. and Steele, B.W., 'Serial Serum Creatine Phosphokinase MB Isoenzyme Activity After Myocardial Infarction', Circulation, 1977; 55(5): 733-737.
10. D'Souza, J.P., Sine, H.E., Horvitz, R.A., Kubasik, N.P., Brody, B.B. and Barold, S.S., 'The Significance of the MB Isoenzyme in Patients with Acute Cardiovascular Disease with a Normal or Borderline Total CPK Activity', Clin. Biochem., 1978; 11(5): 204-209.
11. Mercer, D.W., 'Separation of Tissue and Serum Creatine Kinase Isoenzymes by Ion-Exchange Column Chromatography', Clin. Chem., 1974; 20(1): 36-40.
12. Rao, P.S., Lukes, J.J., Ayres, S.M. and Mueller, H., 'New manual and Automated Method for Determining Activity of Creatine Kinase Isoenzyme MB, by Use of Dithiothreitol: Clinical Applications', Clin. Chem., 1975; 21(11): 1612-1618.
13. Jockers-Wretou, E. and Pfeleiderer, G., 'Quantitation of Creatine Kinase Isoenzymes in Human Tissues and Sera by an Immunological Method', Clin. Chim. Acta, 1975; 58: 223-232.
14. Roberts, R., Sobel, B.E. and Parker, C.W., 'Radioimmunoassay for Creatine Kinase Isoenzymes', Science, 1976; 194: 855-857.
15. Van Der Veen, K.J. and Willebrands, A.F., 'Isoenzymes of Creatine Phosphokinase in Tissue Extracts and in Normal and Pathological Sera', Clin. Chim. Acta, 1966; 13: 312-316.
16. Wong, R. and Swallen, T.O., 'Cellulose Acetate Electrophoresis of Creatine Phosphokinase Isoenzymes in the Diagnosis of Myocardial Infarction', Am. J. Clin. Pathol., 1975; 64: 209-216.
17. Burger, A., Richerich, R. und Aebi, H., 'Die Heterogenitat der Kreatin-Kinase', Biochemische Zeitschrift, 1964; 339: 305-314.
18. Deul, D.H. and van Breeman, J.F.L., 'Electrophoresis of Creatine Phosphokinase from Various Organs', Clin. Chim. Acta, 1964; 10: 276-283.
19. Rosalki, S.B., 'Creatine Phosphokinase Isoenzymes', Nature, 1965; 207: 414.
20. Trainer, T.D. and Gruenig, D., 'A Rapid Method for the Analysis of Creatine Phosphokinase Isoenzymes', Clin. Chim. Acta, 1968; 21: 151-154.
21. Henry, R.J. and Cannon, D.C., Ed's, 'Clinical Chemistry Principles and Techniques', 2nd Edition, Harper and Row, New York, 1974: 901-903.

22. Young, D.S. , 'Effects of Drugs on Clinical Laboratory Tests', 3rd Edition, AACC Press, Washington D.C., 1990.
23. Cho, H.W. and Meltzer, H.Y., 'Factors Affecting Stability of Isoenzymes of Creatine Phosphokinase', *Am. J. Clin. Pathol.*, 1979; 71(1): 75-82.
24. Marmor, A., Alpan, G., Keidar, S., Grenadier, E. and Palant, A., 'The MB Isoenzyme of Creatine Kinase as an Indicator of Severity of Myocardial Ischaemia', *Lancet*, 1978; Oct 14: 812-814.
25. Galen, R.S. Personal Communication, June 1980.
26. Galen, R.S. 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
27. Bayer, P.M., Gabl, F., Granditsch, G., Widhalm, K., Zyman, H. and Deutsch, E., 'Creatine Kinase Isoenzymes in Cerebral Fluid in a Case of Brain Damage', *Clin. Chem.*, 1976; 22(8): 1405-1407.
28. Stein, W. and Bohner, 'Immunoglobulin-Bound Creatine Kinase BB ("Macro CK") in Three Patients with Different Diseases', *Clin. Chem.*, 1979; 25(8): 1513-1514.
29. Pudek, M.R., Jacobson, B.E., Urquhart, N. and Rabkin, S.W., 'Effect of Macro Creatine Kinase BB on Results of Different Creatine Kinase Isoenzyme Methods', *Clin. Chem.*, 1982; 28(6): 1400.
30. Velletri, K., Griffiths, W.C. and Diamond, I., 'Abnormal Electrophoretic Mobility of a Creatine Kinase MM Isoenzyme', *Clin. Chem.*, 1975; 21(12): 1837-1838.
31. Sax, S.M., Moore, J.J., Giegel, J.L. and Welsh, M., 'Atypical Increase in Serum Creatine Kinase Activity in Hospital Patients', *Clin. Chem.*, 1976; 22(1): 87-91.
32. Ljungdahl, L. and Gerhardt, W., 'Creatine Kinase Isoenzyme Variants in Human Serum', *Clin. Chem.*, 1978; 24(5): 832-834.
33. Lott, J.A. and Heinz, J.W., 'A Stable, Unusual Isoenzyme of Creatine Kinase Observed in a Case of Acute Anoxia', *Clin. Chem.*, 1978; 24(6): 1047.
34. James, G.P. and Harrison, R.L., 'Creatine Kinase Isoenzymes of Mitochondrial Origin in Human Serum', *Clin. Chem.*, 1979; 25(6): 943-947.
35. Yuu, Hoyuo, et al., *Clin. Chem.*, 1978; 24(11): 2054-2057.

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1. Brish, L.K., 'CK & LD Isoenzyme. A Self-Instructional Text', Am. Soc. Of Clin. Path. Press, Chicago, 1984 : 30-47.
2. Wolf, P.L., Griffiths, J.C and Koett, J.W., 'Interpretation of Electrophoretic Patterns of Proteins and Isoenzymes', Mason Publishers, New York, 1981 : 60.
3. Tietz, N.W. (Ed.), Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, 1986. 678-700.
4. Galen, R.S., 'The Enzyme Diagnosis of Myocardial Infarction', Human Path., 1975; 6(2): 141-155.
5. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
6. Marmor, A. and Alpan, G., 'Specificity of Creatine Kinase MB Isoenzyme for Myocardial Injury', Clin. Chem., 1978; 24(12) : 2206.
7. Galen, R.S., Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
8. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978, 11 (6): 232-234.
9. Yasmineh, W.G., Pyle, R.B, Cohn, J.N., Nicoloff, D.M., Hanson, N.Q. and Steele, B.W., 'Serial Serum Creatine Phosphokinase MB Isoenzyme Activity After Myocardial Infarction', Circulation, 1977; 55(5): 733-737.
10. D'Souza, J.P., Sine, H.E., Horvitz, R.A., Kubasik, N.P., Brody, B.B. and Barold, S.S., 'The Significance of the MB Isoenzyme in Patients with Acute Cardiovascular Disease with a Normal or Borderline Total CPK Activity', Clin. Biochem., 1978, 11 (5): 204-209.
11. Mercer, D.W., 'Separation of Tissue and Serum Creatine Kinase Isoenzymes by Ion-Exchange Column Chromatography', Clin. Chem., 1974; 20(1) : 36-40.
12. Rao, P.S., Lukes, J.J., Ayres, S.M. and Mueller, H., 'New manual and Automated Method for Determining Activity of Creatine Kinase Isoenzyme MB, by Use of Dithiothreitol: Clinical Applications', Clin. Chem., 1975; 21(11) : 1612-1618.
13. Jockers-Wretou, E. and Pfeleiderer, G., 'Quantitation of Creatine Kinase Isoenzymes in Human Tissues and Sera by an Immunological Method', Clin. Chim. Acta, 1975; 58: 223-232.
14. Roberts, R., Sobel, B.E. and Parker, C.W., 'Radioimmunoassay for Creatine Kinase Isoenzymes', Science, 1976; 194: 855-857.
15. Van Der Veen, K.J. and Willebrands, A.F., 'Isoenzymes of Creatine Phosphokinase in Tissue Extracts and in Normal and Pathological Sera', Clin. Chim. Acta, 1966; 13: 312-316.
16. Wong, R. and Swallen, T.O., 'Cellulose Acetate Electrophoresis of Creatine Phosphokinase Isoenzymes in the Diagnosis of Myocardial Infarction', Am. J. Clin. Pathol., 1975; 64: 209-216.
17. Burger, A., Richterich, R. und Aebi, H., 'Die Heterogenitat der Kreatin-Kinase', Biochemische Zeitschrift, 1964; 339: 305-314.
18. Deul, D.H. and van Breeman, J.F.L., 'Electrophoresis of Creatine Phosphokinase from Various Organs', Clin. Chim. Acta, 1964; 10: 276-283.
19. Rosalki, S.B., 'Creatine Phosphokinase Isoenzymes', Nature, 1965; 207: 414.
20. Trainer, T.D. and Gruenig, D., 'A Rapid Method for the Analysis of Creatine Phosphokinase Isoenzymes', Clin. Chim. Acta, 1968; 21: 151-154.
21. Henry, R.J. and Cannon, D.C., Ed's, 'Clinical Chemistry Principles and Techniques', 2nd Edition, Harper and Row, New York, 1974: 901-903.

22. Young, D.S. , 'Effects of Drugs on Clinical Laboratory Tests', 3rd Edition, AACC Press, Washington D.C., 1990.
23. Cho, H.W. and Meltzer, H.Y., 'Factors Affecting Stability of Isoenzymes of Creatine Phosphokinase', *Am. J. Clin. Pathol.*, 1979; 71(1): 75-82.
24. Marmor, A., Alpan, G., Keidar, S., Grenadier, E. and Palant, A., 'The MB Isoenzyme of Creatine Kinase as an Indicator of Severity of Myocardial Ischaemia', *Lancet*, 1978; Oct 14: 812-814.
25. Galen, R.S., Personal Communication, June 1980.
26. Galen, R.S., 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
27. Bayer, P.M., Gabl, F., Granditsch, G., Widhalm, K., Zyman, H. and Deutsch, E., 'Creatine Kinase Isoenzymes in Cerebral Fluid in a Case of Brain Damage', *Clin. Chem.*, 1976; 22(8) : 1405-1407.
28. Stein, W. and Bohner, 'Immunoglobulin-Bound Creatine Kinase BB ("Macro CK") in Three Patients with Different Diseases', *Clin. Chem.*, 1979; 25(8) : 1513-1514.
29. Pudek, M.R., Jacobson, B.E., Urquhart, N. and Rabkin, S.W., 'Effect of Macro Creatine Kinase BB on Results of Different Creatine Kinase Isoenzyme Methods', *Clin. Chem.*, 1982; 28(6) : 1400.
30. Velletri, K., Griffiths, W.C. and Diamond, I., 'Abnormal Electrophoretic Mobility of a Creatine Kinase MM Isoenzyme', *Clin. Chem.*, 1975; 21(12) : 1837-1838.
31. Sax, S.M., Moore, J.J., Giegel, J.L. and Welsh, M., 'Atypical Increase in Serum Creatine Kinase Activity in Hospital Patients', *Clin. Chem.*, 1976; 22(1) : 87-91.
32. Ljungdahl, L. and Gerhardt, W., 'Creatine Kinase Isoenzyme Variants in Human Serum', *Clin. Chem.*, 1978; 24(5) : 832-834.
33. Lott, J.A. and Heinz, J.W., 'A Stable, Unusual Isoenzyme of Creatine Kinase Observed in a Case of Acute Anoxia', *Clin. Chem.*, 1978; 24(6) : 1047.
34. James, G.P. and Harrison, R.L., 'Creatine Kinase Isoenzymes of Mitochondrial Origin in Human Serum', *Clin. Chem.*, 1979; 25(6) : 943-947.
35. Yuo, Hoyuo, et al., *Clin. Chem.*, 1978; 24(11) : 2054-2057.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1304P 2008/01 (8)