SAS-MX LD Isoenzyme

Instructions For Use REF 101900

SAS-MX Isoenzyme LDH
Fiche technique

SAS-MX LD-Isoenzym
Anleitung

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INTENDED PURPOSE

The SAS-MX LD Isoenzyme Kit is intended for the separation and quantitation of Lactate Dehydrogenase Isoenzymes in serum or plasma by agarose gel electrophoresis.

Lactate dehydrogenase (LD) (EC 1.1.1.27) is an enzyme found in virtually all human tissues, with the liver, skeletal muscle, heart and kidney having the greatest concentrations. The wide distribution of LD in body tissues limits the usefulness of total LD determinations in diagnosis. Testing for the source of elevated LD activity may be indicated with isoenzyme assessment¹.

Five isoenzymes of LD can be demonstrated in human serum. Each isoenzyme is designated by a number which is related to its electrophoretic mobility. The fastest moving fraction (most anodic) is designated LD1 and is found primarily in heart muscle. The slowest moving (most cathodic) is LD5 found primarily in liver and skeletal muscle. The others - LD2, LD3, and LD4 are found in varying degrees along with LD1 and LD5 in all tissues^{1,2}.

The most important use of LD isoenzymes is in the diagnosis of myocardial damage. LD2 is found in highest concentration in normal human serum. The ratio LD1/LD2 is therefore less than one. Following myocardial infarction (MI), there is substantial elevation in LD1 so that the LD1/LD2 ratio following MI will approach or even exceed 1, a phenomenon referred to as 'flipped LD'. The LD level begins to rise approximately 12-24 hours following myocardial infarction, frequently reaching levels two to three times (or greater) the upper limit of normal. Peak activity is usually reached on day 3-4 and activity may remain elevated for as long as two weeks after infarction³.

The most definitive testing in the diagnosis of MI is accomplished by performing creatine kinase (CK) isoenzyme studies in conjunction with LD isoenzyme studies^{1,5}. The specificity and sensitivity achieved with these two tests has eliminated the necessity for additional enzyme studies in accurately diagnosing MI. Studies have shown that CK and LD isoenzyme analyses, in conjunction with the proper clinical setting and electrocardiogram results, are virtually 100% accurate in properly diagnosing myocardial infarction^{2,6}.

The isoenzymes of LD have been determined by various methods⁷⁻¹¹. Electrophoresis provides far more information than the other methods because it allows complete separation of all five isoenzymes with no risk of carryover. The support media used in electrophoresis includes cellulose acetate, agar, agarose and acrylamide gels¹.

The SAS-MX LD Isoenzyme Kit utilises a modification of the method of Preston⁸ and separates the lactate dehydrogenase isoenzymes according to mobility in an agarose gel. The isoenzyme bands are then visualised with a colorimetric reagent which quantitatively visualises the isoenzymes according to the following reaction sequence:

The patterns may be quantitated by densitometry at 595nm.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-MX LD Isoenzyme Gel (x10)

Contains agarose in a barbital buffer with sodium azide as preservative. The gel is ready for use as packaged.

2. LD Buffer Concentrate (2x 100ml)

Contains concentrated AMPD, barbital, bicine and aspartate buffer with sodium azide added as preservative. Dilute the contents of the bottle to 500ml with purified water and mix well.

3. LD Isoenzyme Reagent (10x 1ml)

Contains those components required for the identification of lactate dehydrogenase isoenzymes according to the reaction sequence shown above. See STEP-BY-STEP PROCEDURE for reconstitution instructions.

4. LD Isoenzyme Diluent (1x 15ml)

Contains AMPD, bicine, barbital and aspartate buffer, pH 8.1 - 8.3 with sodium azide as preservative.

5. Other Kit Components

Each kit contains Instructions For Use and sufficient Sample Application Templates and Blotters A, B, C and D to complete 10 gels.

STORAGE AND SHELF LIFE

1. SAS-MX LD Isoenzyme Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. LD Buffer

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

3. LD Isoenzyme Reagent

The reagent vials should be stored at 2...6°C and are stable until the expiry date indicated on the label. Reconstituted reagent is stable for 48 hours at 2...6°C kept in the dark.

4. LD Isoenzyme Diluent

The diluent should be stored at 2...6°C and is stable until the expiry date indicated on the label. Avoid contamination of the diluent. Cloudiness or poor performance may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 5014 Development weight

Cat. No. 4062 Incubation chamber

Drying oven with forced air capable of 60...70°C

Incubator capable of 45°C

Glass plate

Destain Solution: mix 100ml of glacial acetic acid with 900ml of purified water. Store in a tightly stoppered bottle.

SAMPLE COLLECTION AND PREPARATION

Fresh serum is the specimen of choice. Heparin or EDTA anticoagulated plasma can also be used, Oxalate inhibits LD activity and should not be used¹². Plasma should be platelet-free as platelets contain LD13. Samples should be tested as soon as possible after collection as no one storage condition is optimal for all isoenzymes^{12, 14, 15, 16}. If required, samples can be stored refrigerated at 15...30°C or 2...6°C for up to 48 hours¹². DO NOT FREEZE - freezing destroys LD5 activity¹².

Sample Collection: Proper timing of sample collection is critical for the accurate interpretation of LD Isoenzyme patterns in the assessment of MI. A minimum of 3 samples should be collected - one on admission, one 6-13 hours later and the third 24-37 hours following admission²⁻⁴.

Interfering factors:

- 1) Haemolysis may affect LD1 and LD2 quantitations as erythrocytes contain high levels of LD activity^{1,2,12}.
- 2) LD activity is reduced in uraemic sera due to the presence of inhibitors such as urea and oxalate. LD5 is particularly affected¹⁷.
- 3) Acetone and chloroform inhibit all isoenzymes except LD1¹⁴.
- 4) Certain drugs can affect LD isoenzyme activity¹⁸.

STEP-BY-STEP PROCEDURE

1. Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard the blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in step 5.
3. Apply 3µl of sample to each slit and allow to absorb for 3 minutes.
4. Whilst the sample is absorbing, pour 35ml of buffer into each inner section of the SAS-MX Chamber.
5. Following sample absorption, blot the template with the blotter A retained from Step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side down, aligning the positive (+) and the negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 100 volts, 15 minutes.
8. Whilst the samples are electrophoresing, place a blotter D moistened with water into the incubation chamber and place in the incubator set at 45°C to equilibrate.
9. Approximately 3-4 minutes before the end of electrophoresis, reconstitute one vial of LD Isoenzyme Reagent by adding 1 ml of LD Isoenzyme Diluent and mix well.
10. Following electrophoresis, place the gel agarose side up on the glass plate. Pour the contents of the LD Isoenzyme Reagent vial along the anode (+) edge of the gel.
11. Using a serological pipette, spread the reagent from the anode edge to the cathode (-) edge of the gel. Wait 15 seconds and spread the reagent back to the anode edge. Wait 15 seconds and roll the excess off the gel.
12. Place the gel in the incubation chamber. Incubate the gel: 25 minutes, 45°C.
13. Wash the gel in destain solution for 2 minutes.
14. Place the gel agarose side up on the glass plate. Place a blotter B (wetted in destain solution), 3 folded paper towels on to the surface of the gel. Press the gel for 3-5 minutes.
15. Remove the towels and blotter B and place the gel in destain solution for 1 minute.
16. Dry the gel at 60...70°C.

INTERPRETATION OF RESULTS

It is recommended that any evaluation of the gels is performed against normal values which have been produced for this test in each individual laboratory. Completed gels should be stored in the dark at all times.

Qualitative Evaluation: The SAS-MX LD Isoenzyme gel can be visually inspected for the presence or absence of bands of interest.

Quantitative Evaluation: Scan the gels at 595nm, slit 5 within 2 hours of completion. Completed gels should be stored in the dark at all times.

Following electrophoresis, 5 zones of LD activity are normally detected. the fastest zone (LD1) moves with a migration similar to the alpha globulins and the slowest zone (LD5) moves with the gamma globulins. The other 3 zones have intermediate mobilities. The LD activity in serum reflects the breakdown of numerous cell types and all 5 zones can be seen. LD2 predominates, followed by LD1 and LD3. LD4 and LD5 occur in minor amounts.

1. LD2 is the LD Isoenzyme present in the largest amounts in normal serum^{1-4,12}.
2. LD1 is elevated and may be greater than LD2 in:
 - a) Myocardial infarction^{1-4,12}.
 - b) Duchenne's muscular dystrophy presents a pattern like MI but clinical symptoms help in easily differentiating the two diseases^{19,20}.
 - c) Haemolysis (including haemolytic anaemias) should be strongly considered whenever total serum LD reaches levels greater than 5x normal and the isoenzymes show a greater LD1 and LD2. Total LD is much higher in haemolytic anaemia than in MI unless the MI is accompanied by severe shock. Pernicious anaemia in relapse gives an LD pattern like haemolysis and some of the highest total serum LD levels are found with this condition^{2,14}.
 - d) Renal infarct^{2,12}.
3. LD3 is elevated in pulmonary infarctions^{7,12,21}.
4. LD4 elevation has not been associated with any particular pathology.
5. LD5 is elevated in hepatic and muscular damage and skin diseases.
6. When total LD is markedly elevated but all of the isoenzymes have normal percentages, this is called an isomorphic pattern. Widely divergent groups of clinical diagnoses have shown this type of pattern and include cardiorespiratory diseases, malignancy, fracture, diseases of the central nervous system, infection and inflammation, hepatic cirrhosis, alcoholism, trauma without fracture, infectious mononucleosis, hypothyroidism, uraemia, necrosis, pseudomononucleosis, viraemia and intestinal obstruction^{1,2,22}.
7. LD and CK isoenzymes are less specific following open heart surgery than they are in most diagnostic situations. The CK-MB will be elevated due to myocardial damage from the operative procedure and trauma caused by manipulation and cannulation of the heart. The LD2/LD1 may be elevated secondary to haemolysis from extracorporeal circulation²²⁻²⁴.

QUALITY CONTROL

The CK/LD Control (Cat.No.5134) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert provided for acceptable assay values.

LIMITATIONS

The SAS-MX LD Isoenzyme method is not designed to identify tumour markers.

Refer to SAMPLE COLLECTION AND HANDLING for other interfering factors.

Further testing required:

1. Total LD activity may be determined. Conflicting reports exist about the true value of total serum enzyme levels and the severity of a disease^{1,4,25}.
2. In diagnosing MI, CK Isoenzyme studies should also be performed^{1,4}.
3. Haptoglobin studies may be performed to rule out haemolysis as a cause of elevated LD1 and LD2.

BIBLIOGRAPHY

1. Brish, L.K. 'CK & LD Isoenzymes A Self-Instructional Text', Am. Soc. of Clin. Path. Press, Chicago, 1984: 85-120.
2. Hadden, D.M. and Prentiss, T., 'Cardiac Profiling by Electrophoresis', Lab. Mgt., 1977; May: 19-24.
3. Henry, J.D. 'Clinical Diagnosis and Management by Laboratory Methods', 16th Edition, 1979; 1: 366-369, W.B. Saunders Co., Philadelphia.
4. Galen, R.S., Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
5. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978; 11(6): 232-234.
6. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
7. Nerenberg, S.T., 'Electrophoretic Screening Procedures', 1973: 96-125, Lea & Febiger, Philadelphia.
8. Preston, J.A., Briere, R.O. and Batsakis, J.G., 'Rapid Electrophoretic Separation of Lactate Dehydrogenase Isoenzymes on Cellulose Acetate', Am. J. Clin. Pathol., 1965; 43(3): 256-260.
9. Hsu, M-Y., Kohler, M.M., Barolia, L. and Bondar, R.J.L., 'Separation of Five Isoenzymes of Serum Lactate Dehydrogenase by Discontinuous Gradient Elution from a Miniature Ion-Exchange Column', Clin. Chem., 1979; 25(8): 1453-1458.
10. Usategui-Gomez, M., Wicks, R.W. and Warshaw, M., 'Immunochemical Determination of the Heart Isoenzyme of Lactate Dehydrogenase(LDH1) in Human Serum', Clin. Chem., 1979; 25(5): 729-734.
11. Rotenberg, Z., Weinberger, I., Sagie, A., Fuchs, J., Sperling, O. and Agmon, J., 'Lactate Dehydrogenase Isoenzymes in Serum During Recent Acute Myocardial Infarction', Clin. Chem., 1987; 33(8): 1419-1420.
12. Tietz, N.W. et al., 'Textbook of Clinical Chemistry', 1986: 691-700, W.B. Saunders Co., Philadelphia.
13. Rothwell, D.J., Jendrzyczak, B., Becker, M. and Dumas, B.T., 'Lactate Dehydrogenase Activities in Serum and Plasma', Clin. Chem., 1976; 22(7): 1024-1026.
14. Latner, A.L. and Skillen, A.W., 'Isoenzymes in Biology and Medicine', 1968: 146-157, Academic Press, London.
15. Kreuzer, H.H. and Fennis, W.H.S., 'Lactic Dehydrogenase Isoenzymes in Blood Serum after Storage at Different Temperatures', Clin. Chim. Acta, 1964; 9: 64-68.

16. Galen, R.S., Personal Communication, Dec. 1981.
17. Clark, P.I., Kostuk, W.J. and Henderson, A.R., 'Time-Dependency of Human Lactate Dehydrogenase Isoenzyme 5 Inhibition by Urea', *Clin. Chem.*, 1976; 22(12): 2059.
18. Young, D.S., 'Effects of Drugs on Clinical Laboratory Tests', 3rd Ed., 1990, AACC Press, Washington D.C.
19. Roses, A.D., Roses, M.J., Nicholson, G.A. and Roe, C.R., 'Lactate Dehydrogenase Isoenzyme 5 in Detecting carriers of Duchenne Muscular Dystrophy', *Neurology*, 1977; 27(May): 414-421.
20. Yasmineh, W.G., Ibrahim, G.A., Abbasnezhad, M. and Awad, E.A., 'Isoenzyme Distribution of Creatine Kinase and lactate Dehydrogenase in Serum and Skeletal Muscle in Duchenne Muscular Dystrophy, Collagen Disease, and Other Muscular Disorders', *Clin. Chem.*, 1978; 24(11): 1985-1989.
21. Papadopoulos, N.M., 'Clinical Applications of Lactate Dehydrogenase Isoenzymes', *Annals of Clin. Lab. Sci.*, 1977; 7(6): 506-510.
22. Jacobs, D.S., Robinson, R.A., Clark, G.M. and Tucker, J.M., 'Clinical Significance of the Isomorphic Pattern of the Isoenzymes of Serum Lactate Dehydrogenase', *Annals of Clin. and Lab. Sci.*, 1977; 7(5): 411-421.
23. Galen, R.S., 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
24. Mohiuddin, S.M., Raffetto, J., Sketch, M.H., Lynch, J.D., Schultz, R.D. and Runco, V., 'LDH Isoenzymes and Myocardial Infarction in patients Undergoing Coronary Bypass Surgery: An Excellent Correlation', *Am. Heart J.*, 1976; 92(5): 584-588.
25. Chapelle, J-P., Albert, A., Smeets, J-P., Marechal, J-P., Heusghem, C. and Kulbertus, H.E., 'Does Lactate Dehydrogenase Isoenzyme-5 Contribute to the Predictive Power of Total Lactate Dehydrogenase in Myocardial Infarction?' *Clin. Chem.*, 1983; 29(5): 774-777.

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LITERATUR

1. Brish, L.K. 'CK & LD Isoenzymes A Self-Instructional Text', Am. Soc. of Clin. Path. Press, Chicago, 1984 : 85-120.
2. Hadden, D.M. and Prentiss, T., 'Cardiac Profiling by Electrophoresis', Lab. Mgt., 1977; May: 19-24.
3. Henry, J.D. 'Clinical Diagnosis and Management by Laboratory Methods', 16th Edition, 1979; 1: 366-369, W.B. Saunders Co., Philadelphia.
4. Galen, R.S., Reiffel, J.A. and Gambino, S.R., 'Diagnosis of Acute Myocardial Infarction', J. Am. Med. Assoc., 1975; 232(2): 145-147.
5. Frolich, J., Brosseuk, A., Grant, A. and McLennan, M., 'Study of the Value of CPK and LDH Isoenzyme Determinations in the Differential Diagnosis of Ischemic Chest Pain', Clin. Biochem., 1978; 11(6): 232-234.
6. Wagner, G.S., Roe, C.R., Limbird, L.E., Rosati, R.A. and Wallace, A.G., 'The Importance of Identification of the Myocardial-Specific Isoenzyme of Creatine Phosphokinase (MB Form) in the Diagnosis of Acute Myocardial Infarction', Circulation, 1973; 47: 263-269.
7. Nerenberg, S.T., 'Electrophoretic Screening Procedures', 1973: 96-125, Lea & Febiger, Philadelphia.

8. Preston, J.A., Briere, R.O. and Batsakis, J.G., 'Rapid Electrophoretic Separation of Lactate Dehydrogenase Isoenzymes on Cellulose Acetate', *Am. J. Clin. Pathol.*, 1965; 43(3): 256-260.
9. Hsu, M-Y., Kohler, M.M., Barolia, L. and Bondar, R.J.L., 'Separation of Five Isoenzymes of Serum Lactate Dehydrogenase by Discontinuous Gradient Elution from a Miniature Ion-Exchange Column', *Clin. Chem.*, 1979; 25(8) : 1453-1458.
10. Usategui-Gomez, M., Wicks, R.W. and Warshaw, M., 'Immunochemical Determination of the Heart Isoenzyme of Lactate Dehydrogenase(LDH1) in Human Serum', *Clin. Chem.*, 1979; 25(5) : 729-734.
11. Rotenberg, Z., Weinberger, I., Sagie, A., Fuchs, J., Sperling, O. and Agmon, J., 'Lactate Dehydrogenase Isoenzymes in Serum During Recent Acute Myocardial Infarction', *Clin. Chem.*, 1987; 33(8) : 1419-1420.
12. Tietz, N.W. et al., 'Textbook of Clinical Chemistry', 1986: 691-700, W.B. Saunders Co., Philadelphia.
13. Rothwell, D.J., Jendrzyczak, B., Becker, M. and Dumas, B.T., 'Lactate Dehydrogenase Activities in Serum and Plasma', *Clin. Chem.*, 1976; 22(7) : 1024-1026.
14. Latner, A.L. and Skillen, A.W., 'Isoenzymes in Biology and Medicine', 1968: 146-157, Academic Press, London.
15. Kreutzer, H.H. and Fennis, W.H.S., 'Lactic Dehydrogenase Isoenzymes in Blood Serum after Storage at Different Temperatures', *Clin. Chim. Acta*, 1964; 9: 64-68.
16. Galen, R.S., Personal Communication, Dec. 1981.
17. Clark, P.I., Kostuk, W.J. and Henderson, A.R., 'Time-Dependency of Human Lactate Dehydrogenase Isoenzyme 5 Inhibition by Urea', *Clin. Chem.*, 1976; 22(12) : 2059.
18. Young, D.S., 'Effects of Drugs on Clinical Laboratory Tests', 3rd Ed., 1990, AACC Press, Washington D.C.
19. Roses, A.D., Roses, M.J., Nicholson, G.A. and Roe, C.R., 'Lactate Dehydrogenase Isoenzyme 5 in Detecting carriers of Duchenne Muscular Dystrophy', *Neurology*, 1977; 27(May): 414-421.
20. Yasmineh, W.G., Ibrahim, G.A., Abbasnezhad, M. and Awad, E.A., 'Isoenzyme Distribution of Creatine Kinase and lactate Dehydrogenase in Serum and Skeletal Muscle in Duchenne Muscular Dystrophy, Collagen Disease, and Other Muscular Disorders', *Clin. Chem.*, 1978; 24(11) : 1985-1989.
21. Papadopoulos, N.M., 'Clinical Applications of Lactate Dehydrogenase Isoenzymes', *Annals of Clin. Lab. Sci.*, 1977; 7(6): 506-510.
22. Jacobs, D.S., Robinson, R.A., Clark, G.M. and Tucker, J.M., 'Clinical Significance of the Isomorphic Pattern of the Isoenzymes of Serum Lactate Dehydrogenase', *Annals of Clin. and Lab. Sci.*, 1977; 7(5): 411-421.
23. Galen, R.S., 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
24. Mohiuddin, S.M., Raffetto, J., Sketch, M.H., Lynch, J.D., Schultz, R.D. and Runco, V., 'LDH Isoenzymes and Myocardial Infarction in patients Undergoing Coronary Bypass Surgery: An Excellent Correlation', *Am. Heart J.*, 1976; 92(5): 584-588.
25. Chapelle, J-P., Albert, A., Smeets, J-P., Marechal, J-P., Heusghem, C. and Kulbertus, H.E., 'Does Lactate Dehydrogenase Isoenzyme-5 Contribute to the Predictive Power of Total Lactate Dehydrogenase in Myocardial Infarction?' *Clin. Chem.*, 1983; 29(5) : 774-777.

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7. Nerenberg, S.T., 'Electrophoretic Screening Procedures', 1973: 96-125, Lea & Febiger, Philadelphia.
8. Preston, J.A., Briere, R.O. and Batsakis, J.G., 'Rapid Electrophoretic Separation of Lactate Dehydrogenase Isoenzymes on Cellulose Acetate', *Am. J. Clin. Pathol.*, 1965; 43(3): 256-260.
9. Hsu, M-Y., Kohler, M.M., Barolia, L. and Bondar, R.J.L., 'Separation of Five Isoenzymes of Serum Lactate Dehydrogenase by Discontinuous Gradient Elution from a Miniature Ion-Exchange Column', *Clin. Chem.*, 1979; 25(8): 1453-1458.
10. Usategui-Gomez, M., Wicks, R.W. and Warshaw, M., 'Immunochemical Determination of the Heart Isoenzyme of Lactate Dehydrogenase(LDH1) in Human Serum', *Clin. Chem.*, 1979; 25(5): 729-734.
11. Rotenberg, Z., Weinberger, I., Sagie, A., Fuchs, J., Sperling, O. and Agmon, J., 'Lactate Dehydrogenase Isoenzymes in Serum During Recent Acute Myocardial Infarction', *Clin. Chem.*, 1987; 33(8): 1419-1420.
12. Tietz, N.W. et al., 'Textbook of Clinical Chemistry', 1986: 691-700, W.B. Saunders Co., Philadelphia.
13. Rothwell, D.J., Jendrzyczak, B., Becker, M. and Doumas, B.T., 'Lactate Dehydrogenase Activities in Serum and Plasma', *Clin. Chem.*, 1976; 22(7): 1024-1026.
14. Latner, A.L. and Skillen, A.W., 'Isoenzymes in Biology and Medicine', 1968: 146-157, Academic Press, London.
15. Kreutzer, H.H. and Fennis, W.H.S., 'Lactic Dehydrogenase Isoenzymes in Blood Serum after Storage at Different Temperatures', *Clin. Chim. Acta*, 1964; 9: 64-68.
16. Galen, R.S. Personal Communication, Dec. 1981.
17. Clark, P.I., Kostuk, W.J. and Henderson, A.R., 'Time-Dependency of Human Lactate Dehydrogenase Isoenzyme 5 Inhibition by Urea', *Clin. Chem.*, 1976; 22(12): 2059.
18. Young, D.S., 'Effects of Drugs on Clinical Laboratory Tests', 3rd Ed., 1990, AACC Press, Washington D.C.
19. Roses, A.D., Roses, M.J., Nicholson, G.A. and Roe, C.R., 'Lactate Dehydrogenase Isoenzyme 5 in Detecting carriers of Duchenne Muscular Dystrophy', *Neurology*, 1977; 27(May): 414-421.
20. Yasmineh, W.G., Ibrahim, G.A., Abbasnezhad, M. and Awad, E.A., 'Isoenzyme Distribution of Creatine Kinase and lactate Dehydrogenase in Serum and Skeletal Muscle in Duchenne Muscular Dystrophy, Collagen Disease, and Other Muscular Disorders', *Clin. Chem.*, 1978; 24(11): 1985-1989.
21. Papadopoulos, N.M., 'Clinical Applications of Lactate Dehydrogenase Isoenzymes', *Annals of Clin. Lab. Sci.*, 1977; 7(6): 506-510.
22. Jacobs, D.S., Robinson, R.A., Clark, G.M. and Tucker, J.M., 'Clinical Significance of the Isomorphic Pattern of the Isoenzymes of Serum Lactate Dehydrogenase', *Annals of Clin. and Lab. Sci.*, 1977; 7(5): 411-421.
23. Galen, R.S. 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
24. Mohiuddin, S.M., Raffetto, J., Sketch, M.H., Lynch, J.D., Schultz, R.D. and Runco, V., 'LDH Isoenzymes and Myocardial Infarction in patients Undergoing Coronary Bypass Surgery: An Excellent Correlation', *Am. Heart J.*, 1976; 92(5): 584-588.
25. Chapelle, J-P., Albert, A., Smeets, J-P., Marechal, J-P., Heusghem, C. and Kulbertus, H.E., 'Does Lactate Dehydrogenase Isoenzyme-5 Contribute to the Predictive Power of Total Lactate Dehydrogenase in Myocardial Infarction?' *Clin. Chem.*, 1983; 29(5): 774-777.

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7. Nerenberg, S.T., 'Electrophoretic Screening Procedures', 1973: 96-125, Lea & Febiger, Philadelphia.
8. Preston, J.A., Briere, R.O. and Batsakis, J.G., 'Rapid Electrophoretic Separation of Lactate Dehydrogenase Isoenzymes on Cellulose Acetate', *Am. J. Clin. Pathol.*, 1965; 43(3): 256-260.
9. Hsu, M-Y., Kohler, M.M., Barolia, L. and Bondar, R.J.L., 'Separation of Five Isoenzymes of Serum Lactate Dehydrogenase by Discontinuous Gradient Elution from a Miniature Ion-Exchange Column', *Clin. Chem.*, 1979; 25(8) : 1453-1458.
10. Usategui-Gomez, M., Wicks, R.W. and Warshaw, M., 'Immunochemical Determination of the Heart Isoenzyme of Lactate Dehydrogenase(LDH1) in Human Serum', *Clin. Chem.*, 1979; 25(5) : 729-734.
11. Rotenberg, Z., Weinberger, I., Sagie, A., Fuchs, J., Sperling, O. and Agmon, J., 'Lactate Dehydrogenase Isoenzymes in Serum During Recent Acute Myocardial Infarction', *Clin. Chem.*, 1987; 33(8) : 1419-1420.
12. Tietz, N.W. et al., 'Textbook of Clinical Chemistry', 1986: 691-700, W.B. Saunders Co., Philadelphia.
13. Rothwell, D.J., Jendrzyczak, B., Becker, M. and Doumas, B.T., 'Lactate Dehydrogenase Activities in Serum and Plasma', *Clin. Chem.*, 1976; 22(7) : 1024-1026.
14. Latner, A.L. and Skillen, A.W., 'Isoenzymes in Biology and Medicine', 1968: 146-157, Academic Press, London.
15. Kreutzer, H.H. and Fennis, W.H.S., 'Lactic Dehydrogenase Isoenzymes in Blood Serum after Storage at Different Temperatures', *Clin. Chim. Acta*, 1964; 9: 64-68.
16. Galen, R.S., Personal Communication, Dec. 1981.
17. Clark, P.I., Kostuk, W.J. and Henderson, A.R., 'Time-Dependency of Human Lactate Dehydrogenase Isoenzyme 5 Inhibition by Urea', *Clin. Chem.*, 1976; 22(12) : 2059.
18. Young, D.S., 'Effects of Drugs on Clinical Laboratory Tests', 3rd Ed., 1990, AACC Press, Washington D.C.
19. Roses, A.D., Roses, M.J., Nicholson, G.A. and Roe, C.R., 'Lactate Dehydrogenase Isoenzyme 5 in Detecting carriers of Duchenne Muscular Dystrophy', *Neurology*, 1977; 27(May): 414-421.
20. Yasmineh, W.G., Ibrahim, G.A., Abbasnezhad, M. and Awad, E.A., 'Isoenzyme Distribution of Creatine Kinase and lactate Dehydrogenase in Serum and Skeletal Muscle in Duchenne Muscular Dystrophy, Collagen Disease, and Other Muscular Disorders', *Clin. Chem.*, 1978; 24(11) : 1985-1989.
21. Papadopoulos, N.M., 'Clinical Applications of Lactate Dehydrogenase Isoenzymes', *Annals of Clin. Lab. Sci.*, 1977; 7(6): 506-510.
22. Jacobs, D.S., Robinson, R.A., Clark, G.M. and Tucker, J.M., 'Clinical Significance of the Isomorphic Pattern of the Isoenzymes of Serum Lactate Dehydrogenase', *Annals of Clin. and Lab. Sci.*, 1977; 7(5): 411-421.
23. Galen, R.S., 'Diagnostic Tests in Cardiac and Non-Cardiac Disorders' *Diag. Med.*, 1978 Feb., 74-87.
24. Mohiuddin, S.M., Raffetto, J., Sketch, M.H., Lynch, J.D., Schultz, R.D. and Runco, V., 'LDH Isoenzymes and Myocardial Infarction in patients Undergoing Coronary Bypass Surgery: An Excellent Correlation', *Am. Heart J.*, 1976; 92(5): 584-588.
25. Chapelle, J-P., Albert, A., Smeets, J-P., Marechal, J-P., Heusghem, C. and Kulbertus, H.E., 'Does Lactate Dehydrogenase Isoenzyme-5 Contribute to the Predictive Power of Total Lactate Dehydrogenase in Myocardial Infarction?' *Clin. Chem.*, 1983; 29(5) : 774-777.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1305P 2008/01 (7)