## SAS-MX HDL

### Instructions For Use REF. 101000

SAS-MX HDL  
Fiche technique

SAS-MX HDL  
Anleitung

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SAS-MX HDL  
Instrucciones de uso

### Contents


English .....	1
Français .....	6
Deutsch .....	11
Italiano .....	16
Español .....	21

**INTENDED PURPOSE**

The SAS-MX HDL Kit is intended for the separation and quantitation of high density lipoprotein (HDL) cholesterol from the other lipoprotein classes in serum or plasma by agarose gel electrophoresis and enzymatic staining of the cholesterol component of each lipoprotein according to the following reaction sequence:

The relationship of HDL Cholesterol to coronary heart disease (CHD) was reported by Barr et al. in 1951<sup>1</sup>, and Miller and Miller in 1975<sup>2</sup>. The work of Castelli et al.<sup>3,4</sup> focused attention on HDL Cholesterol assessment as the definitive laboratory test in determining the risk of CHD. The clinical recommendations from the National Cholesterol Education Program (NCEP) panels direct clinical laboratories to perform measurements of Total, HDL and LDL Cholesterols and Triglycerides<sup>7</sup>.

The cholesterol content of the lipoprotein fractions has been determined by ultracentrifugation<sup>8</sup>, selective precipitation<sup>9</sup>, and electrophoresis on various media<sup>10</sup>.

**WARNINGS AND PRECAUTIONS**

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Refer to product safety data sheet for risk and safety phrases and disposal information.

**COMPOSITION****1. SAS-MX HDL Gel (x10)**

Contains agarose in a Tris / Barbital buffer with sodium azide as preservative. The gel is ready for use as packaged.

**2. SAS-MX HDL Buffer (1x 100ml)**

Contains Tris, barbital and sodium barbital buffer with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well.

**3. SAS-MX HDL Reagent (10x 1ml)**

Contains those compounds required for the identification of cholesterol according to the reaction sequence shown above. See STEP-BY-STEP-PROCEDURE for reconstitution instructions.

**4. SAS-MX HDL Diluent (1x 27ml)**

Contains Triton X-100 in HEPES buffer, pH 6.9-7.1. The diluent is ready for use as packaged.

**5. Other Kit components**

Each kit contains Instructions For Use and sufficient Sample Application Templates, Blotters A, C and D to complete 10 gels.

## STORAGE AND SHELF LIFE

### 1. SAS-MX HDL Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

### 2. SAS-MX HDL Buffer

Buffer concentrate should be stored at 15...30°C until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C.

### 3. SAS-MX HDL Reagent

The reagent vials should be stored at 2...6°C and are stable until the expiry date indicated on the label. Reconstituted reagent is stable for 72 hours at 2...6°C.

### 4. SAS-MX HDL Diluent

The diluent should be stored at 2...6°C and is stable until the expiry date indicated on the label.

## ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS 600 Power Supply

Cat. No. 4062 Incubation Chamber

Drying oven with forced air capable of 60...70°C.

Incubator capable of 30°C.

Glass Plate

Destain Solution: Mix 50ml glacial acetic acid and 950ml of purified water. Store in a tightly stoppered bottle.

Purified water

## SAMPLE COLLECTION AND PREPARATION

Fresh serum or EDTA anticoagulated plasma is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 4 days. DO NOT FREEZE.

### Patient Preparation:

- 1) If HDL Cholesterol is the only parameter to be determined, the patient need not be fasting. If other lipid values are to be determined then the patient should fast for a 12-14 hour period prior to sampling.
- 2) Discontinue all drugs for 3-4 weeks if possible.
- 3) The patient should be maintaining a stable weight and be on normal diet for at least 1 week.
- 4) Wait 4-8 weeks after a myocardial infarction or similar traumatic episode.

### Interfering Factors:

- 1) Heparin therapy can lead to alterations in the migration of the lipoproteins, particularly beta lipoprotein.

**Sample preparation:** Samples should be diluted 1 + 1 with buffer (1 part sample + 1 part buffer).

**STEP-BY-STEP PROCEDURE**

1. Remove the gel from the packaging and place on a blotter D. Blot the gel surface with a blotter C, discard blotter.
2. Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in step 5.
3. Apply 5 $\mu$ l of sample to each slit and allow to absorb for 5 minutes.
4. Whilst the samples are absorbing, pour 35ml of buffer into each inner section of the SAS-MX Chamber. Place the plastic Heat Shield into the inner section of the chamber.
5. Following sample absorption, blot the template with the blotter A retained from step 2 and remove both blotter and template.
6. Position the gel in the chamber agarose side down, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber.
7. Electrophorese the gel: 80 volts, 35 minutes
8. Whilst the samples are electrophoresing, place the glass plate moistened with water into the incubation chamber and place in the incubator set at 30°C to equilibrate.
9. Approximately 2-4 minutes before the end of electrophoresis, reconstitute 1 vial of HDL reagent by adding 1ml of HDL diluent and mix well. Do not shake.
10. Following electrophoresis, place the gel agarose side up onto a blotter D and place into a drying oven with forced air at 60°C for 2 minutes.
11. Remove the gel from the oven and place on the pre-incubated glass plate. Pour the contents of the HDL reagent vial along the anode (+) edge of the gel.
12. Using a serological pipette, spread the reagent from the anode edge to the cathode (-) edge of the gel. Wait 30 seconds and spread the reagent back to the anode. Wait 30 seconds and roll the excess reagent off gel.
13. Remove the gel from the glass plate and place in the incubation chamber. Incubate the gel: 30 minutes, 30°C.
14. Immerse the gel in destain solution for 5 minutes with gentle agitation.
15. Dry the gel at 60...70°C.
16. Wash the gel in purified water for 5 minutes with gentle agitation.
17. Dry the gel.

**INTERPRETATION OF RESULTS**

It is recommended that any evaluation of the gels is performed against normal values which have been produced for this test in each individual laboratory.

Treatment decisions in the NCEP guidelines are based primarily on LDL Cholesterol levels<sup>7</sup>:

**Total Cholesterol**

Desirable blood cholesterol	<200 mg/dL
Borderline-High blood cholesterol	200-239 mg/dL
High blood cholesterol	>240 mg/dL

## HDL Cholesterol

Low HDL cholesterol	<35 mg/dL
Protective HDL cholesterol	>60 mg/dL

## Triglycerides

Desirable	<250 mg/dL
Borderline	250-500 mg/dL
Elevated	500-1000 mg/dL
Severely elevated/pancreatitis	> 1000 mg/dL

## LDL Cholesterol

Dietary Therapy	Initiation Level	LDL Goal
Without CHD and <2 risk factors	≥ 160 mg/dL	< 160 mg/dL
Without CHD and 2 risk factors	≥ 130 mg/dL	< 130 mg/dL
With CHD	> 100 mg/dL	≤ 100 mg/dL

## Drug Therapy

Without CHD and <2 risk factors	≥ 190 mg/dL	< 160 mg/dL
Without CHD and > 2 risk factors	≥ 160 mg/dL	< 130 mg/dL
With CHD	> 130 mg/dL	≤ 100 mg/dL

The risk factors considered in the classification scheme are Age (males ≥45 years old, females ≥55 years old), family history of premature CHD, smoking, hypertension and diabetes. HDL Cholesterol is considered high risk ≤35 mg/dL and is counted as a risk factor in the classification. HDL Cholesterol >60 mg/dL is considered protective and subtracts one from the total number of risk factors.

## QUALITY CONTROL

The Gel HDL Control (Cat. No. 3218) can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package insert for acceptable assay values.

## BIBLIOGRAPHY

1. Barr, D.P. et al 'Protein-lipid Relationships in Human Plasma', Am. J. Med., 1951 ; 11 : 480-493.
2. Miller, G.J and Miller, N.E. 'Plasma-High Density-Lipoprotein Concentration and Development of Ischemic Heart Disease', Lancet, 1976 ; 1 : 16-19.
3. Kannel, W.B. et al 'Serum Cholesterol, Lipoproteins and the Risk of Coronary Heart Disease' Ann. Inter. Med., 1971 ; 74 (1) : 1-12.
4. Gordon, T. et al 'High Density Lipoprotein as a Protective Factor Against Coronary Artery Disease. The Framingham Study', Am. J. Med., 1977 ; 62 : 707-714.
5. Galen, R.S. 'HDL Cholesterol, How Good a Risk Factor?', Diag. Med., 1979 ; Nov/Dec : 39-58.
6. Castelli, W.P. et al.'HDL Cholesterol and Other Lipids in Coronary Heart Disease - The Cooperative Lipoprotein Phenotyping Study' Circulation, 1977 ; 55 (5) : 767-772.
7. National Cholesterol Education Program, Second report of the expert panel on detection, evaluation and treatment of high blood cholesterol in adults (Adult Treatment Panel II). NCEP 1993.
8. Delalla, O.F. and Gofman, J.W. 'Ultracentrifugal Analysis of Serum Lipoprotein', in Methods of Biochemical Analysis, Vol. 1, D. Glick (Ed.), New York, Interscience, 459-478, 1954.

9. Burstein, M. and Scholnick, H.R., 'Precipitation of Chylomicrons and Very Low density Lipoproteins From Human Serum With Sodium Lauryl Sulfate', *Life Sci.*, 1972; 11 : 177-184.
10. Cobb, S.A. and Sanders J.L., 'Enzymic Determination of Cholesterol in Serum Lipoproteins Separated by Electrophoresis', *Clin. Chem.*, 1978; 24(7) : 1116-1120.

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**BIBLIOGRAPHIE**

1. Barr, D.P. et al. 'Protein-lipid Relationships in Human Plasma', Am. J. Med., 1951 ; 11 : 480-493.
2. Miller, G. J. et Miller, N. E. 'Plasma-High Density-Lipoprotein Concentration and Development of Ischemic Heart Disease', Lancet, 1976 ; 1 : 16-19.
3. Kannel, W.B. et al 'Serum Cholesterol, Lipoproteins and the Risk of Coronary Heart Disease' Ann. Inter. Med., 1971 ; 74 (1) : 1-12.
4. Gordon, T. et al. 'High Density Lipoprotein as a Protective Factor Against Coronary Artery Disease. The Framingham Study', Am. J. Med., 1977 ; 62 : 707-714.
5. Galen, R.S. 'HDL Cholesterol, How Good a Risk Factor?', Diag. Med., 1979 ; nov/déc : 39-58.
6. Castelli, W. P. et al. 'HDL Cholesterol and Other Lipids in Coronary Heart Disease - The Cooperative Lipoprotein Phenotyping Study', Circulation, 1977 ; 55 (5) : 767-772.
7. National Cholesterol Education Program, Second report of the expert panel on detection, evaluation and treatment of high blood cholesterol in adults (Adult Treatment Panel II). NCEP 1993.
8. Delalla, O. F. et Gofman, J. W. 'Ultracentrifugal Analysis of Serum Lipoprotein', in Methods of Biochemical Analysis, Vol. 1, D. Glick (éd.), New York, Interscience, 459-478, 1954.

9. Burstein, M. et Scholnick, H. R., 'Precipitation of Chylomicrons and Very Low density Lipoproteins From Human Serum With Sodium Lauryl Sulfate', *Life Sci.*, 1972 ; 11 : 177-184.
10. Cobb, S. A. et Sanders J. L., 'Enzymic Determination of Cholesterol in Serum Lipoproteins Separated by Electrophoresis', *Clin. Chem.*, 1978; 24(7) : 1116-1120.

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## LITERATUR

1. Barr, D.P. et al 'Protein-lipid Relationships in Human Plasma', Am. J. Med., 1951 ; 11 : 480-493.
2. Miller, G.J and Miller, N.E. 'Plasma-High Density-Lipoprotein Concentration and Development of Ischemic Heart Disease', Lancet, 1976 ; 1 : 16-19.
3. Kannel, W.B. et al 'Serum Cholesterol, Lipoproteins and the Risk of Coronary Heart Disease' Ann. Inter. Med., 1971 ; 74 (1) : 1-12.
4. Gordon, T. et al 'High Density Lipoprotein as a Protective Factor Against Coronary Artery Disease. The Framingham Study', Am. J. Med., 1977 ; 62 : 707-714.
5. Galen, R.S. 'HDL Cholesterol, How Good a Risk Factor?', Diag. Med., 1979 ; Nov/Dec : 39-58.
6. Castelli, W.P. et al.'HDL Cholesterol and Other Lipids in Coronary Heart Disease - The Cooperative Lipoprotein Phenotyping Study' Circulation, 1977; 55 (5): 767-772.
7. National Cholesterol Education Program, Second report of the expert panel on detection, evaluation and treatment of high blood cholesterol in adults (Adult Treatment Panel II). NCEP 1993.
8. Delalla, O.F. and Gofman, J.W. 'Ultracentrifugal Analysis of Serum Lipoprotein', in Methods of Biochemical Analysis, Vol. 1, D. Glick (Ed.), New York, Interscience, 459-478, 1954.

9. Burstein, M. and Scholnick, H.R., 'Precipitation of Chylomicrons and Very Low density Lipoproteins From Human Serum With Sodium Lauryl Sulfate', *Life Sci.*, 1972; 11 : 177-184.
10. Cobb, S.A. and Sanders J.L., 'Enzymic Determination of Cholesterol in Serum Lipoproteins Separated by Electrophoresis', *Clin. Chem.*, 1978; 24(7) : 1116-1120.

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**BIBLIOGRAFIA**

1. Barr, D.P. et al 'Protein-lipid Relationships in Human Plasma', Am. J. Med., 1951 ; 11 : 480-493.
2. Miller, G.J and Miller, N.E. 'Plasma-High Density-Lipoprotein Concentration and Development of Ischemic Heart Disease', Lancet, 1976 ; 1 : 16-19.
3. Kannel, W.B. et al 'Serum Cholesterol, Lipoproteins and the Risk of Coronary Heart Disease' Ann. Inter. Med., 1971 ; 74 (1) : 1-12.
4. Gordon, T. et al 'High Density Lipoprotein as a Protective Factor Against Coronary Artery Disease. The Framingham Study', Am. J. Med., 1977 ; 62 : 707-714.
5. Galen, R.S. 'HDL Cholesterol, How Good a Risk Factor?', Diag. Med., 1979 ; Nov/Dec : 39-58.
6. Castelli, W.P. et al 'HDL Cholesterol and Other Lipids in Coronary Heart Disease - The Cooperative Lipoprotein Phenotyping Study' Circulation, 1977 ; 55 (5) : 767-772.
7. National Cholesterol Education Program, Second report of the expert panel on detection, evaluation and treatment of high blood cholesterol in adults (Adult Treatment Panel II). NCEP 1993.

8. Delalla, O.F. and Gofman, J.W. 'Ultracentrifugal Analysis of Serum Lipoprotein', in *Methods of Biochemical Analysis*, Vol. 1, D. Glick (Ed.), New York, Interscience, 459-478, 1954.
9. Burstein, M. and Scholnick, H.R., 'Precipitation of Chylomicrons and Very Low density Lipoproteins From Human Serum With Sodium Lauryl Sulfate', *Life Sci.*, 1972; 11 : 177-184.
10. Cobb, S.A. and Sanders J.L., 'Enzymic Determination of Cholesterol in Serum Lipoproteins Separated by Electrophoresis', *Clin. Chem.*, 1978; 24(7) : 1116-1120.

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1. Barr, D.P. et al. 'Protein-lipid Relationships in Human Plasma', Am. J. Med., 1951 ; 11 : 480-493.
2. Miller, G.J y Miller, N.E. 'Plasma-High Density-Lipoprotein Concentration and Development of Ischemic Heart Disease', Lancet, 1976 ; 1: 16-19.
3. Kannel, W.B. et al. 'Serum Cholesterol, Lipoproteins and the Risk of Coronary Heart Disease' Ann. Inter. Med., 1971 ; 74 (1): 1-12.
4. Gordon, T. et al. 'High Density Lipoprotein as a Protective Factor Against Coronary Artery Disease. The Framingham Study', Am. J. Med., 1977 ; 62 : 707-714.
5. Galen, R.S. 'HDL Cholesterol, How Good a Risk Factor?', Diag. Med., 1979 ; Nov/Dic: 39-58.
6. Castelli, W.P. et al. 'HDL Cholesterol and Other Lipids in Coronary Heart Disease - The Cooperative Lipoprotein Phenotyping Study' Circulation, 1977 ; 55 (5) : 767-772.
7. Programa Nacional de Educación sobre el Colesterol, Segundo informe del panel de expertos sobre la detección, evaluación y tratamiento de la hipercolesterolemia en adultos (Panel de Tratamiento de Adultos II). NCEP 1993.

8. 'Ultracentrifugal Analysis of Serum Lipoprotein', in Methods of Biochemical Analysis, Vol. 1, D. Glick (Ed.), Nueva York, Interscience, 459-478, 1954.
9. Burstein, M. y Scholnick, H.R., 'Precipitation of Chylomicrons and Very Low density Lipoproteins From Human Serum With Sodium Lauryl Sulfate', Life Sci., 1972; 11 : 177-184.
10. Cobb, S.A. y Sanders J.L., 'Enzymic Determination of Cholesterol in Serum Lipoproteins Separated by Electrophoresis', Clin. Chem., 1978; 24(7) : 1116-1120.

Helena Biosciences Europe  
Queensway South  
Team Valley Trading Estate  
Gateshead  
Tyne and Wear  
NE11 0SD

Tel: +44 (0) 191 482 8440  
Fax: +44 (0) 191 482 8442  
Email: [info@helena-biosciences.com](mailto:info@helena-biosciences.com)  
[www.helena-biosciences.com](http://www.helena-biosciences.com)

HL-2-1309P 2008/01 (7)