**helena** | BioSciences  
Europe  
www.helena-biosciences.com

Helena BioSciences Europe  
Colima Avenue  
Sunderland Enterprise Park  
Sunderland  
SR5 3XB  
tel: +44 (0) 191 549 6064  
fax: +44 (0) 191 549 6271  
email: info@helena-biosciences.com

Other Helena BioSciences Europe offices:

Helena BioSciences Europe  
6 Rue Charles Cros-ZAE  
95320 Saint Leu La Foret  
France  
tel: +33 13 995 9292  
fax: +33 13 995 6891  
email: helena@helena.fr

Helena BioSciences Europe  
Via Enrico Fermi, 24  
20090 Assago (Milano)  
Italy  
tel: +39 02 488 1951  
or +39 02 488 2141  
fax: +39 02 488 2677

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## **SAS-MX ALP ISOENZYME**

### **INTENDED PURPOSE**

The SAS-MX ALP Isoenzyme kit is intended for the separation and quantitation of Alkaline Phosphatase Isoenzymes by agarose gel electrophoresis.

Alkaline Phosphatase (ALP) (EC 3.1.3.1) is an enzyme which catalyses the hydrolysis of phosphate esters at alkaline pH. The greatest concentrations of ALP are found in the liver, bone, intestine and placenta. However, practically every body tissue contains at least a small amount of ALP. Because of this wide distribution, the usefulness of a total ALP quantitation is limited. Fortunately, each source of ALP produces one predominant isoenzyme, and the tissue source of elevated ALP in serum can be determined by identifying the isoenzyme type.

The isoenzymes of ALP differ in their physicochemical and electrophoretic properties and, by taking advantage of these differences, the individual isoenzymes can be identified. In addition to the liver (including macrohepatic, or 'fast liver'), bone, intestinal and placental isoenzymes, other ALP isoenzymes have been identified in serum. These include Regan, Nagao, PA and renal isoenzymes.

The SAS-MX ALP Isoenzyme procedure is a high resolution electrophoresis method, utilizing differences in the response of the isoenzymes to neuraminidase<sup>2,3</sup> and a detergent in the gel to separate liver and macrohepatic (fast liver). The result is a separation capable of separating liver, bone, macrohepatic and intestinal ALP isoenzymes for use in the diagnosis and treatment of liver, bone, parathyroid and intestinal disorders. This system may separate the intestinal isoenzyme into 3 distinct bands - the clinical significance of these has not been determined.

### **WARNINGS AND PRECAUTIONS**

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

### **COMPOSITION**

#### **1. SAS-MX ALP Isoenzyme Gel**

Contains agarose in a Tris / Barbitol buffer with sodium azide and thiomersal as preservative. The gel is ready for use as packaged.

#### **2. Tris / Barbitol Buffer Concentrate**

Contains barbitol and sodium barbitol with sodium azide as preservative. Dilute the contents of the bottle to 1 litre with purified water and mix well. Buffer salts may crystallize slightly on standing. Wash any crystals from the bottle with diluted buffer to ensure complete dissolution.

#### **3. SAS-MX ALP Separation Enhancer**

Contains neuraminidase from *Vibrio Cholerae* (EC 3.2.1.18). The separation enhancer is ready for use as packaged.

#### **4. SAS-MX ALP Isoenzyme Substrate**

Contains BCIP (5-Bromo-4-chloro-3-indolyl phosphate) in 2-amino-2-methyl-1,3-propanediol buffer. The isoenzyme substrate is ready for use as packaged.

#### **5. SAS-MX ALP Isoenzyme Chromagen**

Contains NBT (nitroblue tetrazolium) powder. See STEP-BY-STEP PROCEDURE for reconstitution instructions.

#### **6. Other Kit Components**

Each kit contains Instructions For Use, sufficient Sample Application Templates, Blotters A, B, C, D and Wicks to complete 10 gels.

**STORAGE AND SHELF-LIFE****1. SAS-MX Alkaline Phosphatase Gel**

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE THE PLATES. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

**2. Tris / Barbital Buffer**

The buffer concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C. Cloudiness or poor performance of the diluted buffer may indicate deterioration.

**3. SAS-MX ALP Isoenzyme Separation Enhancer**

The separation enhancer should be stored at 2...6°C and is stable until the expiry date indicated on the label. Cloudiness or poor performance may indicate deterioration.

**4. Alkaline Phosphatase Substrate**

The alkaline phosphatase substrate should be stored at 2...6°C and is stable until the expiry date indicated on the label. Discolouration or poor performance may indicate deterioration.

**5. Alkaline Phosphatase Chromagen**

The alkaline phosphatase chromagen should be stored at 2...6°C and is stable until the expiry date indicated on the label. Avoid contamination. Reconstituted chromagen should be used within 30 minutes of preparation.

**ITEMS REQUIRED BUT NOT PROVIDED**

Cat. No. 4063 SAS-MX Chamber

Cat. No. 1525 EPS600 Power Supply

Cat. No. 3039 Chamber Cooling Device

Cat. No. 5014 Development Weight

Cat. No. 4062 Incubation chamber

Drying Oven with forced air capable of 60...70°C

Incubator capable of 45°C

Glass plate

Vortex

Destain solution: Mix 100ml of glacial acetic acid and 900ml of purified water. Store in a tightly stoppered bottle.

Purified water

**SAMPLE COLLECTION AND PREPARATION**

Freshly collected serum is the specimen of choice. Serum should be removed from the red cells as soon as possible after blood collection. Samples containing oxalate, citrate or EDTA cannot be used as these substances inhibit alkaline phosphatase activity<sup>4</sup>. Total ALP level should be determined by an appropriate method.

Patients should be fasting. Patients with blood groups O or B who are secretors may have an elevated Intestinal ALP around 2 hours after a fatty meal<sup>5,6</sup>.

It is recommended that serum be tested as soon as possible after collection and that the serum should be stored at 2...6°C. If necessary, the serum should be stored frozen at -20°C for no more than 24 hours<sup>7,8</sup>.

**SAS-MX ALP ISOENZYME**


**Interfering Factors:** 1) Phosphate, oxalate, citrate and cyanide can inhibit ALP activity<sup>4,8</sup>. 2) Excess glycine may inhibit ALP by complexing Magnesium ions, 3) EDTA can inhibit some ALP isoenzymes<sup>4</sup>, 4) Several drugs cause an enzymatic imbalance which may affect ALP level<sup>9</sup>.

Samples with high total ALP activity can be diluted if the high activity affects the ability to discriminate between isoenzyme types. **NOTE:** Dilutions should not be performed in saline or water, as the migration of the isoenzymes will be affected. Dilutions should be performed in serum from a male donor which has been heat treated (56°C for 15-20 minutes) to destroy endogenous ALP activity. The suitability of the serum for use as a diluent should be checked before use.

**STEP-BY-STEP PROCEDURE**

**NOTE:** Ensure that the Chamber Cooling Device has been stored at 2...6°C for at least 60 minutes prior to use.

- To every 25µl of sample / control add 5µl of separation enhancer and incubate at 15...30°C for 5 minutes. **NOTE:** For optimal separation it is important to keep to the times stated.
- Remove the gel from the packaging and place on a paper towel. Blot the gel surface with a blotter C, discard blotter.
- Align the sample application template with the arrows at the edge of the gel. Place a blotter A on top of the template and rub a finger across the slits to ensure good contact. Remove the blotter and retain for use in Step 6.
- Apply 5µl of sample to each slit and allow to absorb for 10 minutes.
- Whilst the sample is absorbing, pour 100ml of buffer into each outer section of the SAS-MX Chamber. Place the chamber cooling device into the central trough of the chamber and wet the surface with a few drops of buffer.
- Following sample absorption, blot the template with the blotter A retained from step 3 and remove both blotter and template.
- Position the gel on top of the cooling device agarose side up, aligning the positive (+) and negative (-) sides with the corresponding positions on the chamber. Take care to avoid trapped air bubbles under the gel.
- Prepare a wick for each side of the gel by placing 3 wicks together in 2 sets, to make 2 thick wicks. Immerse each wick in the buffer, gently squeeze out excess buffer and attach each wick to the edge of the gel, parallel with the edge of the cooling device. Ensure one edge of the wick is immersed in the buffer and rub one finger across the portion of the wick on the gel to ensure good contact.
- Electrophorese the gel: 250 volts, 30 minutes.
- Whilst the samples are electrophoresing, place a tissue moistened with water into the incubation chamber and place in the incubator set at 45°C to equilibrate.
- Approximately 10 minutes before the end of electrophoresis, reconstitute 1 vial of ALP isoenzyme chromagen by adding 1ml of ALP substrate. Vortex mix for 30 seconds and allow to stand undisturbed until use.
- Following electrophoresis, place the gel agarose side up, onto the glass plate. Pour the contents of the ALP isoenzyme vial along one edge of the gel.
- Using a serological pipette, spread the reagent across the gel surface. Wait 15 seconds and repeat spreading, wait 15 seconds and roll the excess off the gel.
- Place the gel in the incubation chamber. Incubate the gel: 30 minutes, 45°C.
- Wash the gel in destain solution for 5 minutes.

## SAS-MX ALP ISOENZYME

16. Place the gel on a blotter D agarose side up, and place a blotter B (wetted in destain solution) onto the surface of the gel followed by 2 blotter D. Press the gel for 5 minutes using the development weight.
17. Remove the weight and blotters and rinse the gel in purified water.
18. Dry the gel at 60...70°C.

### INTERPRETATION OF RESULTS.

#### 1. Qualitative Evaluation:

The gels may be visually inspected for the presence or absence of particular bands of interest.

#### 2. Quantitative Evaluation:

Scan the SAS-MX ALP Isoenzyme gels (gel side down) within 2 hours of completion at 595nm.


In either case, an increase or decrease in particular serum components or the detection of unusual serum components require investigation. The completed SAS-MX ALP Isoenzyme gel should be stored in the dark to prevent increased background colouration.

Interpretation of ALP isoenzyme patterns should not be attempted without knowledge of the total ALP level in the patient's serum. Serum from normal individuals may contain small amounts of liver, bone and intestinal ALP<sup>8,10,11</sup>. ALP levels are age and sex dependant<sup>8,12,13</sup>.

Pregnant females may show a placental band. The macrohepatic (fast liver) band seen in neoplasms should be interpreted as an alert to a disease state regardless of total ALP level.

The performance of the rare Nagao, Regan and PA isoenzymes with this system are not known at this time. Abnormal bands have been reported in patients with normal total ALP levels.

Figure 1 (below) shows a diagrammatic representation of the migration characteristics of the main ALP Isoenzyme bands.

The liver band migrates the most anodic of all the bands. The liver band on patients with a high total ALP will migrate more anodal than that on a normal level patient. Liver is the enzyme most frequently elevated when total ALP levels are elevated<sup>6,10</sup>. The liver ALP increases early in liver disease before most other liver function tests show abnormalities. The wide group of conditions leading to increased liver ALP include acute hepatitis, cirrhosis, fatty liver, drug induced liver disease, obstruction of biliary flow by carcinoma at the head of the pancreas, bile duct stricture, primary biliary cirrhosis and metastatic carcinoma of the liver<sup>9</sup>.

Next to the liver band, running more cathodal is the bone band. In samples containing high activities of the bone isoenzyme, the bone band runs more slowly than that of a normal patient. The bone band often has a characteristic square shape, as opposed to the more normal oblong shape of the electrophoretic bands. Bone isoenzyme is increased as a results of increased osteoblast activity. This isoenzyme is normally elevated in growing children and adults over the age of 50. The highest total ALP values have been attributed to increased bone isoenzyme in Paget's disease or renal rickets<sup>14</sup>. An abnormally high bone isoenzyme level may also be indicative of bone cancer, osteomalacia or celiac sprue<sup>5</sup>. A decreased bone ALP in children may be attributed to cretinism or hypophosphatasia.

Running cathodal to the bone band is the macrohepatic (fast liver) band. Macrohepatic ALP has been isolated in cases of metastatic carcinoma to the liver and has been suggested as a diagnostic tool in identifying such cases. It has also been isolated in patients with viral hepatitis, alcoholic cirrhosis and other liver diseases. Data generated by Viot and associates<sup>15</sup> suggest that macrohepatic ALP is highly correlated with the presence of liver metastases and that the presence of this isoenzyme could be predictive of the appearance of liver metastases. Viot also reports that macrohepatic ALP is seen occasionally in patients free of any disease state<sup>15</sup>.

The intestinal ALP isoenzyme band(s) run cathodal to the macrohepatic. The appearance of up to three bands is more common in samples from non-fasting patients.

The identification of ALP isoenzymes such as Regan, Nagao and PA may require the use of additional testing protocols such as heat inactivation<sup>16</sup>, amino acid inhibition<sup>17,18</sup>, urea denaturation<sup>17,18</sup> or separation on polyacrylamide gels<sup>20,21</sup>.

### QUALITY CONTROL

Alkaline Phosphatase Isoenzyme Control (Cat. No. 3059) contains normal levels of Liver and Bone isoenzymes and can be used to verify all phases of the procedure and should be used on each plate run. Refer to the package inserts for acceptable assay values.

### REFERENCE VALUES

It is recommended that any evaluation of the gels is performed against normal values which have been produced for this test in each individual laboratory.

An indication of normal ALP Isoenzymes can be found in Rhone, *et al*<sup>8</sup>, Fritsche *et al*<sup>10</sup> and Lee *et al*<sup>11</sup>.


### PERFORMANCE CHARACTERISTICS


#### Reproducibility

	Within-Run (n=15)		Between-Run (n=75)	
	Mean	CV (%)	Mean	CV (%)
Liver	51.8%	3.4	51.8%	4.5
Bone	48.3%	3.6	48.2%	4.8

#### Sensitivity

The method is sensitive to 51 +/- 5.7IU per band, determined as the lowest concentration of ALP isoenzyme which is evident as a discrete band on the completed gel.


### Linearity

The linearity of the method is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.

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14. Nerenburg, S.T. 'Medical Technology', Lea and Febiger Inc., Philadelphia, 122-123, 1973.