SAS-I Urine Analysis

Instructions For Use REF 200400

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INTENDED PURPOSE

The SAS-I Urine Analysis Kit is intended for the screening of unconcentrated urine samples by agarose gel electrophoresis.

Urinary proteins are derived primarily from plasma proteins that filter through the kidney. The appearance of abnormal plasma proteins in the urine is of great value in evaluating renal function. The appropriate study of proteinuria should include quantitative and qualitative assessment of the type and amount of proteins excreted^{1,5}. The combination of electrophoretic separation of urine proteins, coupled to the identification of specific protein types by immunoprecipitation allows the differentiation of several types of proteinuria - Physiological, Glomerular (selective and non-selective), Tubular and proteinuria associated with Dysglobulinaemias^{1,5}. The SAS-I Urine Analysis kit separates urine proteins according to charge in an agarose gel. The proteins are then stained to allow visualisation and quantitative or qualitative interpretation. The high sensitivity stain used in the kit allows most urine samples to be tested without prior concentration, and the use of the Urine Protein Internal Standard allows an approximation of urine total protein, or individual band protein concentrations to be determined by densitometry.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION**1. SAS-I Urine Analysis Gel**

Contains agarose in a Tris / Barbitol buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

2. Acid Violet Stain Concentrate

Contains concentrated Acid Violet stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

3. Destain Solution Concentrate

Dilute the contents of Destain A to 1 litre with purified water. Then add the contents of Destain B and add a further 1 litre of purified water, slowly.

4. Wash Solution

Contains concentrated Wash Solution. Dilute the contents of the bottle to 2 litres with saline solution.

5. Other Kit Components

Each kit contains Instructions For Use and sufficient Blotters B, C, D, X and Blotter Combs to complete 10 gels.

STORAGE AND SHELF-LIFE**1. SAS-I Urine Analysis Gel**

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources

2. Acid Violet Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration of the stain solution.

3. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted destain solution is stable for 6 months at 15...30°C.

4. Wash Solution

The Wash Solution concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted wash solution is stable for 6 months at 15...30°C. Cloudiness may indicate deterioration of the wash solution.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 210200 Sample Applicator Blades (1 x 10)

Cat. No. 210300 Sample Applicator Blades (5 x 10)

Cat. No. 210100 Disposable sample cups (100)

Cat. No. 5014 Development Weight

Cat. No. 3100 REP Prep

Drying oven with forced air capable of 60...70°C

Fixative Solution: Mix 500ml methanol with 500ml purified water. Store in a tightly stoppered bottle.

Saline solution (0.85% NaCl)

Purified water

The following items are not required for standard urine analysis but may be required for further investigation:

Cat. No. 220100 Antiserum to Human Urine Total Protein (2ml)

Cat. No. 220200 Antiserum to Human Urine Micro Proteins (2ml)

Cat. No. 220300 Antiserum to Human Urine Macro Proteins (2ml)

Cat. No. 220400 Antiserum to Human GAM Proteins (2ml)

Cat. No. 220700 Antiserum to Human Free & Bound Kappa Light Chain (2ml)

Cat. No. 220800 Antiserum to Human Free & Bound Lambda Light Chain (2ml)

Cat. No. 220500 Antiserum to Human Free Kappa Light Chain (2ml)

Cat. No. 220600 Antiserum to Human Free Lambda Light Chain (2ml)

Cat. No. 220900 Antiserum to Human Urine Pentavalent/Albumin (2ml)

Cat. No. 221000 Antiserum to Human Urine Pentavalent (2ml)

Cat. No. 9249 Antiserum to Human IgD

Cat. No. 9250 Antiserum to Human IgE

SAMPLE COLLECTION AND PREPARATION

Freshly collected urine is the specimen of choice. Samples can be stored refrigerated at 2...6°C for up to 72 hours or 2 weeks at -20°C. Samples should initially be used without concentration or dilution, unless the total protein is known to exceed 1000mg/L. If scanning the completed gel by densitometry, an indication of the Urine Total Protein, or the concentration of individual protein zones can be performed by the incorporation of the Urine Protein Internal Standard (Cat. No. 3032) into the sample preparation procedure. Full Instructions For Use are provided with each Urine Protein Internal Standard.

STEP-BY-STEP PROCEDURE

1. Pipette 35µl of the sample into the appropriate well of the SAS-I sample tray or disposable sample cups.
 - i) **SAS-I & SAS-I Plus users:** Carefully place the sample tray onto the applicator drawer. Ensure that the tray is pushed firmly down into position.
 - ii) **SAS-3 users:** Carefully locate the sample tray using the sample base locating pins. Ensure that the tray is positioned securely.
2. Remove the gel from the packaging and:
 - i) **SAS-I users:** place the gel in the SAS-I, agarose side up, aligning the positive and negative sides with the corresponding electrode posts.
 - ii) **SAS-I Plus users:** dispense 400µL of REP Prep onto the heat sink. Place the gel onto the heat sink, agarose side up, aligning the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
 - iii) **SAS-3 users:** place the alignment guide onto the pins and dispense 400µL of REP Prep onto the centre of the chamber. Place the gel into the chamber agarose side up, using the guide, align the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
3. Blot the surface of the gel with a blotter C, discard the blotter.
4.
 - i) **SAS-I users:** attach the electrodes onto the top side of the electrode posts so that they are in contact with the gel blocks.
 - ii) **SAS-I Plus users:** (as above). Place the cover over the gel and electrodes and press firmly for 5 seconds to ensure contact.
 - iii) **SAS-3 users:** attach the electrodes onto the the electrode posts so that they are in contact with the gel blocks.
5. Place two applicator blade assemblies in position on the instrument, (**SAS-3 users:** slot A and I0).
6. Perform the Urine Analysis electrophoresis:

- i) **SAS-I users:**
 Urine Analysis: (with internal standard, cat. no. 3032) 100 volts, 16 mins, 10 applications
 Urine Analysis: (without internal standard, cat. no. 3032) 100 volts, 17 mins, 10 applications
 Immunofixation: 80 volts, 17 mins, 10 applications
- ii) **SAS-I Plus users:**
 Urine Analysis: (with internal standard, cat. no. 3032) 100 volts, 20 mins, 20°C, 10 applications
 Urine Analysis: (without internal standard, cat. no. 3032) 100 volts, 21 mins, 20°C, 10 applications
 Immunofixation: 80 volts, 20 mins, 20°C, 10 applications
- iii) **SAS-3 users:**

Step	Time (mm:ss)	Temperature (°C)	Voltage	Other
Load Sample	00:10	21		Speed 1
Apply Sample	00:10	21		Speed 1* ¹
Repeat steps 'load sample, apply sample' for a total of 10 applications.				
Electrophoresis	18:00	20	100	
Dry	08:00	54* ²		

*¹ Use Location 2

*² Do not include dry step if immunofixation is required.

NOTE: For Urine Analysis, remove buffer blocks prior to drying.

7. Following electrophoresis, (**SAS-I Plus users:** remove the cover), remove the electrodes from the surface of the gel. **NOTE:** Gels can be stained to show all proteins present in the sample (see staining protocol), or the gels can be incubated with specific antisera to show particular proteins of interest (see immunofixation protocol).

STAINING PROTOCOL

8. **SAS-I and SAS-I Plus users:** Remove both gel blocks using the Gel Block Remover.
 9. Attach the gel to the staining chamber holder.
 10. Select the Urine test program on the staining unit and, following the prompts, Fix, Stain, Destain, Wash, and Dry the gel.
- a) **SAS-2 (Auto-Stainer)**

Step	Solution	Time (mm:ss)	Port	Temperature (°C)
Fix	Fixative solution	05:00	4	
Dry	—	15:00		65
Stain	Acid Violet Stain	10:00	5	
Wash	Purified water	00:10	1	
Wash	Purified water	00:10	1	
Destain	Destain solution	02:00	2	
Destain	Destain solution	02:00	2	
Wash	Water	01:00	1	
Dry	—	15:00		65

- b) **SAS-4 (Auto-Stainer)**

Step	Time (mm:ss)	Temperature (°C)	Other
Stain	04:00		Recirculate ON
Destain	02:00		Recirculate ON
Destain	02:00		Recirculate ON
Dry	12:00	64	

- c) **Manual**

Follow the sequence listed for the SAS-2 Auto-Stainer, using a staining bath for the Fix, Stain, Destain and Wash steps, and a Drying Oven with forced air at 60...70°C for the Dry steps.

11. At the end of the staining cycle, remove the gel from the staining chamber. The gel is now ready for examination.

IMMUNOFIXATION PROTOCOL

8. Program the following parameters into the instrument.
- i) **SAS-I Plus users:** Incubation step 1: 10 mins, 37°C (incubate)
 Incubation step 2: 8 mins, 40°C (D blot)
- ii) **SAS-3 users:**

Step	Time (mm:ss)	Temperature (°C)
Apply antisera	10:00	21
Insert combs	02:00	21
Blotter D	05:00	40
Dry	08:00	54

9. **(SAS-3 users: remove the alignment guide).** Position the antiserum application template onto the gel surface. **NOTE:** The milled antisera channels should be aligned centrally over the printed box on the gel in which the samples are applied.
10. Apply 2 drops (or 50µl) of the appropriate antiserum / fixative into the hole of the immunoglobulin lanes. Ensure that the antisera has completely filled the channels.
11. **SAS-I users:** Incubate the gel: 15...30°C, 10 minutes.
SAS-I Plus users: Incubate the gel.
12. Following incubation, place a blotter comb into the holes of the antiserum template. Allow 2 minutes for the excess antisera to be absorbed, then remove the blotter combs and the template from the gel surface.
13. **i) SAS-I users:** Remove the gel blocks using the Gel Block Remover and wash the gel in wash solution for 5 minutes.
ii) SAS-I Plus and SAS-3 users: Place a blotter D (smooth side down) onto the surface of the gel, leave for 10 seconds and remove.
14. **i) SAS-I users:** Place the gel on a blotter D agarose side up and place a blotter B (wetted in wash solution) onto the surface of the gel followed by two blotter X. Press the gel using the Development Weight for 10 minutes.
ii) SAS-I Plus and SAS-3 users: Place a blotter D (smooth side down) onto the surface of the gel and replace the antiserum template to hold the blotter flat. Blot the gel.
15. **i) SAS-I users:** Remove the blotters and place the gel in wash solution for 4 minutes with gentle agitation.
ii) SAS-I Plus users: Remove the blotter D.
iii) SAS-3 users: Remove the blotter D and dry the gel.
16. **i) SAS-I users:** Remove the gel from the wash solution and place on a blotter D agarose side up. Place a blotter B (wetted in wash solution) onto the surface of the gel followed by a blotter D. Press the gel for 3 minutes.
ii) SAS-I Plus and SAS-3 users: Remove the gel blocks using the Gel Block Remover. Go to Step 18.
17. **i) SAS-I users:** Remove the blotters.
18. Attach the gel to the staining chamber holder.
NOTE: Immediately after use, clean the antisera template with a mild, biocidal detergent. If possible, scrub the bottom of the template with a toothbrush or small test tube brush. Do not let the antisera dry on the template. Build up of antisera on the surface of the template will result in the formation of bubbles during the antisera application step. Dry the antisera template thoroughly. Water left in the holes will hinder the application of antisera for the next use. Store the template upside down to increase the air circulation and thus the drying potential of the applicator.
19. Select the Urine IFE test program on the staining unit and, following the prompts, Wash, Stain, Destain and Dry the gel.

a) SAS-2 (Auto-Stainer)

Step	Solution	Time (mm:ss)	Port	Temperature (°C)
Dry	—	10:00		55
Wash	Wash solution	07:00	4	
Stain	Acid Violet Stain	03:00	5	
Destain	Destain solution	02:00	2	
Dry	—	05:00		65
Wash	Wash solution	03:00	4	
Wash	Wash solution	03:00	4	
Dry	—	05:00		65

b) SAS-4 (Auto-Stainer)

Step	Time (mm:ss)	Temperature (°C)	Other
Wash	00:03		Recirculate ON
Wash	10:00		Recirculate ON
Stain	04:00		Recirculate ON
Destain	02:00		Recirculate ON
Destain	02:00		Recirculate ON
Dry	12:00	63	

c) Manual

Follow the sequence listed for the SAS-2 Auto-Stainer, using a staining bath for the Wash, Stain and Destain steps, and a Drying Oven with forced air at 60...70°C for the Dry steps.

20. At the end of the staining cycle, remove the gel from the staining chamber. The gel is now ready for examination.

INTERPRETATION OF RESULTS

The majority of monoclonal proteins migrate in the cathodic, gamma region of the protein pattern, but due to their abnormal nature, they may migrate anywhere within the globulin region on protein electrophoresis. The monoclonal protein band on the immunofixation pattern will occupy the same position and shape as the abnormal band on the serum protein pattern. The abnormal protein is identified by the antiserum type it reacts with.

When low concentrations of abnormal protein are present, the abnormal band may appear as a band within the normal polyclonal immunoglobulin. A band can also be seen within a polyclonal background when there is a large polyclonal immunoglobulin presence also.

The publication 'Immunofixation for the Identification of Monoclonal Gammopathies' is available from Helena BioSciences on request.

Visual interpretation of bands present on the gel by immunofixation:

Type of Proteinuria	Bands Observed On Gel	Proteins Present
Normal urine	Small albumin band	albumin
Glomerular	Albumin, alpha-1, beta, gamma	albumin, alpha-1 antitrypsin, transferrin, gamma globulins
Tubular	alpha-1, alpha-2, beta	Retinol Binding Protein, beta2-microglobulin, alpha-2 microglobulin
Overflow	gamma or variable	immunoglobulins, free light chains

LIMITATIONS

1. Antigen Excess

Antigen excess will occur if there is not a slight antibody excess or antigen/antibody equivalence at the site of precipitation. Antigen excess in IFE is usually due to an excess of the immunoglobulin in the patient sample. Antigen excess is characterised by prozoning (unstained areas in the centre of the immunofixed protein band, with staining around the edges). A higher dilution of the sample should be used in this event to optimise the immunoglobulin concentration.

2. Band In Cathodic End of Gamma Region Showing No Reactivity With IFE Antisera.

C Reactive Protein (CRP) may be detected in patients with acute inflammatory response^{6,7}. CRP appears as a narrow band at the cathodic end of the serum protein pattern. Elevated Alpha1-Antitrypsin and Haptoglobin are supportive evidence for CRP. Patients with a CRP band will probably have an elevated level when assayed for CRP.

3. Non-Reactivity With Kappa and Lambda Antisera

Occasionally a sample will have a reaction with a heavy chain antiserum but no light chain reaction is obvious. In this situation, the following need to be ruled out - a) Heavy chain disease, b) Very high concentrations of light chains, leading to antigen excess, c) Low concentrations of light chains, d) Atypical light chain molecule that does not react with the antiserum, e) Light Chains with 'hidden' light chain determinants (as sometimes seen with IgA and IgD). To obtain definitive results, testing may include a) A higher or lower dilution of the sample to optimise the antibody/antigen equivalence, b) Antisera from more than one manufacturer to aid in the identification of atypical immunoglobulins, and c) Treat the sample with b-2-mercaptoethanol to 'reveal' the light chains.

PERFORMANCE CHARACTERISTICS

a) Reproducibility

Band	Within-Run (n=24)		Between-Run (n=4)	
	Mean (%)	CV (%)	Mean	CV (%)
Band 1	54.9	2.7	53.6	3.0
Band 2	5.0	8.3	5.2	9.5
Band 3	8.9	9.1	9.3	9.2
Band 4	9.7	7.3	9.8	7.5
Band 5	6.8	8.9	7.1	10.4
Band 6	14.6	6.8	12.6	6.9

b) **Sensitivity:** 50ng per band (equivalent to 10mg/L in the sample), determined as the lowest concentration of protein which was evident as a discrete band on the completed gel.

c) **Linearity** Linearity is a function of densitometer specification as well as gel performance. It is recommended that each customer determine the linearity of the method based upon the densitometer in use in the laboratory.

BIBLIOGRAPHY

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox. Sang., 1969; 17 : 445-452.
7. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
8. Microprotein - PR™ (Procedure No. 611) Sigma Diagnostics, September 1995.

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BIBLIOGRAFIA

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.

5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox. Sang., 1969; 17 : 445-452.
7. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
8. Microprotein - PR™ (Procedure No. 611) Sigma Diagnostics, September 1995.

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