(en) SAS-MX IFE Antisera to IgD & IgE

INTENDED PURPOSE
SAS-MX IFE Antisera to IgD and IgE are intended for use in identifying monoclonal gammopathies using immunofixation.

WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

- SAS-MX IFE Antisera to IgD (Cat. No. 9409)**
The antisera contains monospecific antisera to human immunoglobulin heavy chains IgD. The antisera has been prepared in horse, sheep or goat. Contains sodium azide as preservative. The antisera is ready for use as packaged.
- SAS-MX IFE Antisera to IgE (Cat. No. 9410)**
The antisera contains monospecific antisera to human immunoglobulin heavy chains IgE. The antisera has been prepared in horse, sheep or goat. Contains sodium azide as preservative. The antisera is ready for use as packaged.
- Other kit components**
Each kit contains Instructions For Use.

STORAGE AND SHELF-LIFE

- SAS-MX IFE Antisera to IgD (Cat. No. 9409)**
The antisera should be stored at 2...6°C and is stable until the expiry date indicated on the label. Deterioration may be indicated by extremely cloudy antisera indicative of bacterial contamination.
- SAS-MX IFE Antisera to IgE (Cat. No. 9410)**
The antisera should be stored at 2...6°C and is stable until the expiry date indicated on the label. Deterioration may be indicated by extremely cloudy antisera indicative of bacterial contamination.

STEP-BY-STEP PROCEDURE

The IgD and IgE antisera can be used with the appropriate Helena BioSciences Europe kits: Refer to the appropriate Instructions For Use for the STEP-BY-STEP PROCEDURE.

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1. Antiseri SAS-MX IFE anti-IgD (Cod. N. 9409)
2. Antiseri SAS-MX IFE anti-IgE (Cod. N. 9410)
Antisueros frente a IgD e IgE SAS-MX IFE
1. Antisueros frente a IgD SAS-MX IFE (Nº Cat. 9409)
2. Antisueros frente a IgE SAS-MX IFE (Nº Cat. 9410)