**Instructions For Use**

**SAS-MX IFE Antisera Free Kappa / Lambda Light Chain**  
Cat. No. 9412 / 9413

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20090 Assago (Milano)

Helena Biosciences Europe

Other Helena Biosciences Europe offices:

HL-2-139FP 2003/11 (2)

**(en) SAS-MX IFE Antisera Free Kappa / Lambda Light Chain**

Cat. No. 9412 Free Kappa Light Chain Antisera  
Cat. No. 9413 Free Lambda Light Chain Antisera

**INTENDED PURPOSE**

SAS-MX IFE Antisera to Free Kappa and Free Lambda Light Chains are intended for use in identifying monoclonal gammopathies using immunofixation electrophoresis.

The IFE Free Kappa and Free Lambda Light Chain Antisera provide an excellent method for identifying monoclonal immunoglobulins in urine samples. These antisera can be used with the SAS-MX IFE Trivalent Antisera for screening urine samples. The reaction provides information about the excretion of the whole molecule of the monoclonal immunoglobulin. Both serum and urine can be tested on the same gel when using the SAS-MX Immunofix-Plus Kit thus providing a complete work-up on the patient.

**WARNINGS AND PRECAUTIONS**

For in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Refer to the product safety data sheet for risk and safety phrases and disposal information.

**COMPOSITION**

**1. SAS-MX IFE Free Kappa Light Chain Antisera, Cat. No. 9412**

Each vial contains monospecific antiserum to human immunoglobulin light chains. The antiserum was prepared in goat and contains sodium azide as preservative.

**2. SAS-MX IFE Free Lambda Light Chain Antisera, Cat. No. 9413**

The antisera are ready for use as packaged. Each vial contains monospecific antiserum to human immunoglobulin light chains. The antiserum was prepared in goat and contains sodium azide as preservative.

**3. Other kit components**

Each kit contains Instructions For Use.

- 1. SAS-MX IFE Free Kappa Light Chain Antisera, Cat. No. 9412**  
The antisera should be stored at 2...6°C and are stable until the expiration date indicated on the label. Deterioration may be indicated by extremely cloudy antiserum indicative of bacterial contamination.
- 2. SAS-MX IFE Free Lambda Light Chain Antisera, Cat. No. 9413**  
The antisera should be stored at 2...6°C and are stable until the expiration date indicated on the label. Deterioration may be indicated by extremely cloudy antiserum indicative of bacterial contamination.

**ITEMS REQUIRED BUT NOT PROVIDED**

The IFE Free Kappa and Free Lambda Light Chain Antisera are optional antisera which can be used with the following kits:

Cat. No. 100300 SAS-MX Immunofix Kit  
Cat. No. 100400 SAS-MX Immunofix-Plus Kit  
Cat. No. 100600 SAS-MX Urine Immunofix Kit

**STEP-BY-STEP PROCEDURE**

Refer to the appropriate kit Instructions For Use for detailed instructions.

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## (es) Cadenas sueltas y ligeras Kappa/Lambda SAS-MX IFE
**1. Antisuero para cadenas sueltas y ligeras Kappa SAS-MX IFE, Nº Cat. 9412**
**2. Antisuero para cadenas sueltas y ligeras Lambda SAS-MX IFE, Nº Cat. 9413**
Nº Cat. 100300 Kit ImmunoFix SAS-MX
Nº Cat. 100400 Kit ImmunoFix-Plus SAS-MX
Nº Cat. 100600 Kit ImmunoFix de orina SAS-MX