Instructions For Use

Protein C Rocket EID

Cat. No. 5357

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INTENDED PURPOSE

The Protein C Rocket EID (electroimmunodiffusion) procedure is intended for the quantitative determination of plasma protein C antigen by Laurell rocket electrophoresis^{1,2}.

Protein C is a vitamin K dependent plasma protein that functions as a regulator of fibrin formation. In its activated form, it inhibits thrombin formation by the inactivation of activated factors V and VIII^{3,4}. A deficiency of protein C constitutes a thrombotic risk factor^{5,6} of which superficial thrombophlebitis is the most common clinical feature. The virtual absence of plasma protein C has led to fatal thrombosis in neonates⁸.

The Helena Protein C Rocket EID Procedure is performed in a 1% agarose gel medium containing an antiserum specific for protein C. After the plasma specimens are applied to the wells in the agarose, electrophoresis is used to migrate the proteins into the antibody field. A rocket-shaped precipitin pattern forms along the axis of migration. The length of this rocket pattern is proportional to the antigen concentration.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in-vitro* diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody, however they should be handled with the same precautions as a human plasma sample.

COMPOSITION**1. Protein C Antigen Rocket Plates (Cat. No. 5357)**

Contains sheep or goat antibody to human protein C incorporated into agarose in Tris-tricine buffer⁸ and sodium azide as preservative. To prevent the formation of toxic vapors, sodium azide should not be mixed with acidic solutions.

Preparation: Remove the plate from the protective packaging and allow 5-20 minutes for the agarose to reach 15...30°C.

2. Tris-tricine Buffer (Cat. No. 5358 - not included)

Diluted buffer contains 0.08 M Tris and 0.024 M tricine.

Preparation: Dilute one package of buffer to 1000ml with deionized or distilled water. The buffer is ready for use when all material is completely dissolved.

3. Rocket Stain (Cat. No. 5362 - not included)

Rocket Stain is Coomassie Brilliant Blue stain.

Preparation: Dissolve the contents of the vial in 450ml deionized water, 450ml methanol and 100ml acetic acid. Mix thoroughly and filter before use if necessary.

4. Other kit components

Each kit contains Instructions For Use, a rocket ruler and report form.

STORAGE AND SHELF LIFE

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

1. Protein C Antigen Rocket Plates (Cat. No. 5357)

Rocket Plates **MUST** be stored at 2...6°C and maintained in moist condition within the bag.

DO NOT FREEZE. The plates are stable until the expiry date indicated on the package.

Signs of Deterioration: Discard the plate if dry in appearance or if the wells are not round. A crystalline appearance indicates the agarose has been frozen.

2. Tris-tricine Buffer (Cat. No. 5358)

The packaged buffer should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted buffer is stable for 2 months at 15...30°C.

Signs of Deterioration: Discard packaged buffer if the material shows signs of dampness or discoloration. Discard diluted buffer if it becomes turbid.

3. Rocket Stain (Cat. No. 5362)

The stain should be stored at 15...30°C and is stable until the expiry date indicated on the label.

Signs of Deterioration: If methanol evaporation occurs, a metallic sheen will be visible on the stain surface. Discard the stain if it does not adequately stain protein rockets as described in this procedure.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 5358 Tris - Tricine Buffer

Cat. No. 5362 Rocket Stain

Cat. No. 5185 SARP reference plasma

Cat. No. 1283 / 4063 Zip Zone Chamber / Titan Gel Chamber

Cat. No. 5014 Development Weight

Cat. No. 1520 EWS Power Supply

Cat. No. 9015 Sponge Wicks

Cat. No. 5037 Blotter Pads

Cat. No. 1558 TITAN GEL Multi-staining Set

Cat. No. 9025 IEP VuBox

Cat. No. 5116 I.O.D. - Incubator, Oven, Drier

Cat. No. 6210, 6211 Microdispenser and Tubes (10µl)

Lint free tissue

Destain solution: Mix 200ml deionized water, 200ml methanol and 50ml glacial acetic acid.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000-3000xg for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP-BY-STEP PROCEDURE

NOTE: Remove the plate from the refrigerator, allow approximately 5-20 minutes for the plate to reach 15...30°C and for excess moisture to be absorbed before use.

1. Reconstitute one vial of S.A.R.P. with 1.0ml deionized or distilled water. Make dilutions for preparation of the Standard Curve as follows:

Percent Activity	Dilution	Parts S.A.R.P.	Parts 0.85% Saline
100%	Use reconstituted S.A.R.P. undiluted		
50%	1:2	1	1
25%	1:4	1	3
12.5%	1:8	1	7

2. Dilute each patient sample and control with 0.85% saline. Prepare a 1:2 dilution (1 part patient sample and 1 part saline) and a 1:4 dilution (1 part patient plasma and 3 parts saline). Additional dilutions may be necessary depending on the patient history. Suspected abnormal samples may need to be tested undiluted.
3. **Preparation:**
 - If using Zip Zone Chamber:** Pour 200ml of Tris-tricine Buffer into each of the outer sections of the chamber (requires a total of 400ml buffer) and place one sponge wick in the buffer along each inner wall of the chamber.
 - If using TITAN GEL Chamber:** Pour 65ml Tris-tricine buffer into each inner section of the chamber.
4. Remove any excess buffer from the wells of the plate. **NOTE:** Excess moisture on the plate can result in poor rockets.
5. Apply 10 μ l of each dilution of the patient samples and controls to the designated wells. **NOTE:** Standard curve samples must be run on each plate. Duplicate applications of patient samples are advisable. When applying the samples to the plate wells, do not allow the pipette tip to touch the sides of the wells as this may cause damage.
6. Allow 5 minutes for specimens to diffuse into the agarose.
7. **Electrophoresis:**
 - If using Zip Zone Chamber:** Place the plate, agarose side down, into the chamber on the sponge wicks. Place the application point (wells) toward the cathode (negative) side of the chamber.
 - If using TITAN GEL Chamber:** Place the plate into the inner section of chamber, agarose side down, by gently squeezing the gel into place. Position the gel(s) so that the edges of the agar are in the buffer and the wells are toward the cathode (-) side of the chamber.
8. Electrophorese the plates at a constant current of 16mA per plate for 3 hours.
9. Following electrophoresis, remove the plates from the chamber. Discard the chamber buffer after each run. **NOTE:** Use the TITAN GEL Multi-Staining Set (Cat. No. 1558) as a staining and rinsing chamber.
10. Rinse the plate with deionized or distilled water and wash it in 0.85% saline overnight with gentle stirring.
11. Following the overnight wash, rinse the plate with deionized or distilled water.
12. Place the plate on a flat surface, agarose side up. Cover the agarose with a single, lint-free tissue.
13. Place 2-3 Blotter Pads and a Development Weight on the plate for 15 minutes and remove.
14. Dry the plate in a drying oven at 60...70°C for 10-20 minutes. **DO NOT** over dry plates. The plate will be transparent when completely dry. **NOTE:** If a dryer/oven is not available, the plates may be covered with wet lint-free tissues and allowed to dry at 15...30°C overnight or under a fan for 3 hours at 15...30°C as climate requires.
15. Following drying, stain the plate by immersing it in Rocket Stain for 20 minutes.
16. Place the plate in destain solution for 1-3 minutes. **NOTE:** Destaining is complete as soon as the background sufficiently clears in order to easily distinguish rocket peaks. **NOTE:** Excessive destaining may fade the rockets making correct measurements difficult. If over destaining does occur, repeat Step F.7. and stain the rockets again.
17. Rinse the plate twice in purified water for 5-10 minutes each rinse.
18. Dry the plates at 37°C for 5 minutes or at 15...30°C until dry.

- Place the plate on the Helena I.E.P. VuBox using a piece of white paper in the bottom of the VuBox for easier viewing of the rockets. Mark the apex of each rocket peak with marker.
- Using the Helena Rocket Ruler, measure the length of each peak in millimeters. The peak is measured from the top of each well to the apex of the rocket.
- Plot the values of the standard curve versus each rocket height on the Rocket Antigen Report Form or on 3 cycle semi-logarithmic paper. Draw the line of best fit for the four points. See Figures 1 and 2 in INTERPRETATION OF RESULTS for an example of a completed Rocket Plate and a standard curve drawn on a Rocket Antigen Report Form.

INTERPRETATION OF RESULTS

Read the patient values from the standard curve and multiply each by the appropriate dilution factor. If Rocket Reference Plasma is used to prepare the standard curve, the patient value read from that curve must be multiplied by the assigned Protein C Antigen value of the appropriate lot of Rocket Reference Plasma as well as the dilution factor.

Patient value from curve = 30%

Dilution factor = 2


S.A.R.P. assigned value = 98%

Actual Patient Factor


Protein C Antigen = $30\% \times 2 \times 0.98 = 58.8\%$

Patient samples with Protein C Antigen levels greater than the range of the standard curve, must be reassayed using the appropriate dilutions.

Figure 1: Rocket patterns on a Protein C Antigen Rocket Plate.

The lengths of the rockets (in millimeters) of the standard dilutions are used to prepare the standard curve. Patient results are read from the curve.

Figure 2: A representative standard curve prepared with S.A.R.P. on 3 cycle semi-logarithmic paper.**QUALITY CONTROL**

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 S.A.C.-1

Cat. No. 5302 S.A.C.-2

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range but be aware that there is currently no accepted international standard for Protein C antigen..

Protein C antigen values are usually expressed in relative percentages compared to a pooled normal plasma standard. The expected normal range for protein C antigen run by rocket electrophoresis has been reported at 65-129%⁸ with newborns showing less than 50% levels. Bertina et al. , reported a range of 65-145% in healthy individuals with diminished levels found following anticoagulant therapy. There is apparently no difference in protein C antigen levels between healthy males and females⁹.

Helena tested 39 plasmas of presumed healthy men and women. The results were as follows:

$$\bar{X} = 101\%$$

$$SD = 24$$

$$\text{Range} = 53 - 149\%$$

LIMITATIONS

Most patients with congenital protein C deficiency show diminished levels of both immunologic and functional activity, the so-called type I deficiency⁷. However, the finding of several patients with normal levels of protein C antigen and diminished functional protein C activity, type II deficiency⁷, makes the diagnosis more difficult. The Helena procedure will detect only the type I deficiencies.

PERFORMANCE CHARACTERISTICS

Precision Studies

Within-Run: A control plasma was tested in replicate on one plate with the following results.

$$\begin{aligned}n &= 8 & \bar{X} &= 91.8\% \\ & & SD &= 2.83 \\ & & CV\% &= 3.0\end{aligned}$$

Run-to-Run: A control plasma was tested in replicate on four different plates giving the following data.

$$\begin{aligned}n &= 26 & \bar{X} &= 85.2\% \\ & & SD &= 6.47 \\ & & CV\% &= 8.0\end{aligned}$$

Comparison Studies

Correlation studies were done on 43 normal and abnormal patient samples with the Helena Protein C method and the reference Protein C method. The study yielded an excellent linear regression equation and correlation coefficient.

$$\begin{aligned}n &= 43 & Y &= 0.614X + 27.006 \\ & & r &= 0.937 \\ & & X &= \text{Helena's Protein C} \\ & & Y &= \text{Reference Protein C method}\end{aligned}$$

BIBLIOGRAPHY

1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest 29 Suppl 124, 21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15:45-52, 1966.
3. Nesheim, M.E., Canfield, W.M., Kisiel, W., Mann, K.G. Studies of the Capacity of Factor Xa to Protect Factor Va from Inactivation by Activated Protein C. J Biol Chem, 257:1443-1447, 1982.
4. Marljar, R.A., Kleiss, A.J., Griffin, J.H. Human Protein C: Inactivation of Factor V and VIII in Plasma by the Activated Molecule. Ann N.Y. Acad Sci, 370:303-310, 1981.
5. Griffin, J.H., Evatt, B., Zimmermann, T.S., Kleiss, A.J. Wideman, C. Deficiency of Protein C in Congenital Thrombotic Disease. J Clin Invest, 68:1370-1373, 1981.
6. Broekmans, A.W., Veltkamp, J.J., Bertina, R.M. Congenital Protein C Deficiency and Venous Thrombo-embolism. A Study in Three Dutch Families. N Engl J Med, 309:340-344, 1983.
7. Bertina, R.M., Broekmans, A.W., Krommenhoek-Van Es, C., Van Wijngaarden, A. The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency. Thromb Haemostas, 51:1-5, 1984.

8. Seligsohn, U., Berger, A., Abend, M., Rubin, L., Attias, D., Zivelin, A., Rapaport, S.I. Homozygous Protein C Deficiency Manifested by Massive Venous Thrombosis in the Newborn. *N Engl J Med*, 310:559-562, 1984.
9. Pabinger-Fasching, I., Bertina, R.M., Lechner, K., Niesser, H., Korninger, Ch. Protein C Deficiency in Two Austrian Families. *Thromb Haemostas*, 50:810-813, 1983.

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BIBLIOGRAPHIE

1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest 29 Suppl 124, 21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15:45-52, 1966.
3. Nesheim, M.E., Canfield, W.M., Kiesel, W., Mann, K.G. Studies of the Capacity of Factor Xa to Protect Factor Va from Inactivation by Activated Protein C. J Biol Chem, 257:1443-1447, 1982.
4. Marlar, R.A., Kleiss, A.J., Griffin, J.H. Human Protein C: Inactivation of Factor V and VIII in Plasma by the Activated Molecule. Ann N.Y. Acad Sci, 370:303-310, 1981.
5. Griffin, J.H., Evatt, B., Zimmermann, T.S., Kleiss, A.J. Wideman, C. Deficiency of Protein C in Congenital Thrombotic Disease. J Clin Invest, 68:1370-1373, 1981.
6. Broekmans, A.W., Veltkamp, J.J., Bertina, R.M. Congenital Protein C Deficiency and Venous Thrombo-embolism. A Study in Three Dutch Families. N Engl J Med, 309:340-344, 1983.
7. Bertina, R.M., Broekmans, A.W., Krommenhoek-Van Es, C., Van Wijngaarden, A. The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency. Thromb Haemostas, 51:1-5, 1984.
8. Seligsohn, U., Berger, A., Abend, M., Rubin, L., Attias, D., Zivelin, A., Rapaport, S.I. Homozygous Protein C Deficiency Manifested by Massive Venous Thrombosis in the Newborn. N Engl J Med, 310:559-562, 1984.
9. Pabinger-Fasching, I., Bertina, R.M., Lechner, K., Niesser, H., Korninger, Ch. Protein C Deficiency in Two Austrian Families. Thromb Haemostas, 50:810-813, 1983.

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1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest 29 Suppl 124, 21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15:45-52, 1966.
3. Nesheim, M.E., Canfield, W.M., Kisiel, W., Mann, K.G. Studies of the Capacity of Factor Xa to Protect Factor Va from Inactivation by Activated Protein C. J Biol Chem, 257:1443-1447, 1982.
4. Marlar, R.A., Kleiss, A.J., Griffin, J.H. Human Protein C: Inactivation of Factor V and VIII in Plasma by the Activated Molecule. Ann N.Y. Acad Sci, 370:303-310, 1981.
5. Griffin, J.H., Evatt, B., Zimmermann, T.S., Kleiss, A.J. Wideman, C. Deficiency of Protein C in Congenital Thrombotic Disease. J Clin Invest, 68:1370-1373, 1981.
6. Broekmans, A.W., Veltkamp, J.J., Bertina, R.M. Congenital Protein C Deficiency and Venous Thrombo-embolism. A Study in Three Dutch Families. N Engl J Med, 309:340-344, 1983.
7. Bertina, R.M., Broekmans, A.W., Krommenhoek-Van Es, C., Van Wijngaarden, A. The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency. Thromb Haemostas, 51:1-5, 1984.
8. Seligsohn, U., Berger, A., Abend, M., Rubin, L., Attias, D., Zivelin, A., Rapaport, S.I. Homozygous Protein C Deficiency Manifested by Massive Venous Thrombosis in the Newborn. N Engl J Med, 310:559-562, 1984.
9. Pabinger-Fasching, I., Bertina, R.M., Lechner, K., Niesser, H., Korninger, Ch. Protein C Deficiency in Two Austrian Families. Thromb Haemostas, 50:810-813, 1983.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1424P 2007/05 (3)