von Willebrand Factor Antigen Rocket EIA

Instructions For Use

REF 5361

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INTENDED PURPOSE

The von Willebrand Factor (vWF) Antigen Rocket System is intended for the quantitative determination of von Willebrand Factor Antigen in human plasma by Laurell rocket electroimmunoassay^{1,2}.

The accurate quantitative measurement of von Willebrand Factor Antigen (vWF:Ag) in human plasma aids in differentiating patients with von Willebrand's disease from those with Hemophilia A and in detecting carriers of Hemophilia A³. Quantitation of von Willebrand Factor Antigen by rocket electroimmunoassay (EIA) is a procedure which combines the specificity of immunochemistry with the speed of electrophoresis⁴.

The Helena von Willebrand Factor Antigen Rocket EIA Method is performed in an agarose gel medium containing a monospecific antiserum to the von Willebrand Factor Antigen. After plasma samples are applied to the wells in the agarose, electrophoresis is used to migrate the proteins into the antibody field. A rocket-shaped precipitin pattern forms along the axis of migration with the length of this rocket pattern being proportional to the antigen concentration.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in-vitro* diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody, however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

1. von Willebrand Factor Antigen Rocket Plates (Cat. No. 5361)

Each plate contains antiserum to human von Willebrand Factor Antigen in barbital buffered agarose⁵ and sodium azide as preservative.

Preparation: Remove the plate from the bag, remove the plastic lid and allow 5-10 minutes for the agarose to reach 15...30°C.

2. Electra™ BI Buffer (Cat. No. 5016 - not provided)

When dissolved, the buffer contains barbital-sodium barbital buffer with sodium azide as preservative.

Preparation: Dissolve one package of buffer in 1000ml deionized water. The buffer is ready for use when all material is completely dissolved.

3. Rocket Stain (Cat. No. 5362 not provided)

Rocket Stain is a preparation of Coomassie Brilliant Blue stain.

Preparation: Dissolve the contents of one vial in 450ml methanol, 450ml deionized water and 100ml acetic acid.

4. Other Kit Components

Each kit contains Instructions For Use, a rocket ruler and report form.

STORAGE AND SHELF-LIFE

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

1. von Willebrand Factor Antigen Rocket Plates (Cat. No. 5361)

Rocket Plates MUST be stored at 2...6°C and maintained within the bag. **DO NOT FREEZE.**

The plates are stable until the expiry date indicated on the package.

Signs of Deterioration: Discard the plate if dry in appearance or if the wells are not round. A crystalline appearance indicates the agarose has been frozen and should be discarded.

2. Electra™ B1 Buffer (Cat. No. 5016)

The buffer should be stored at 15...30°C. Diluted buffer is stable for 2 months at 15...30°C.

Signs of Deterioration: Discard packaged buffer if the material shows signs of dampness or discoloration. Discard diluted buffer if it becomes turbid.

3. Rocket Stain (Cat. No. 5362)

The stain should be stored at 15...30°C.

Signs of Deterioration: If methanol evaporation occurs, a metallic sheen will be visible on the stain surface. Discard the stain if it does not adequately stain protein rockets as described in this procedure.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 5016 Electra B₁ Buffer

Cat. No. 5362 Rocket Stain

Cat. No. 5185 SARP reference plasma

Cat. No. 5037 Titan Blotter Pads

Cat. No. 9015 Sponge Wicks

Cat. No. 1283 Zip Zone Chamber

Cat. No. 5014 Development Weight

Cat. No. 9036 Humidity Chamber

Cat. No. 1520 EWS Power Supply

Cat. No. 5116 Incubator, Oven, Drier

Cat. No. 6210, 6211 Microdispenser and Tubes (10µl)

Laboratory rotator

Lint free tissue

Destain solution: Mix 200ml deionized water, 200ml methanol and 50ml glacial acetic acid.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000-3000g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. **DO NOT** keep at 37°C for more than 5 minutes.

STEP-BY-STEP PROCEDURE

NOTE: Remove the plate from the refrigerator, allow approximately 5-20 minutes for the plate to reach 15...30°C and for excess moisture to be absorbed before use.

1. Make dilutions of reconstituted S.A.R.P. for preparation of the Standard Curve as follows:

Percent Activity	Dilution	Parts S.A.R.P.	Parts 0.85% Saline
100%	* Use reconstituted S.A.R.P. undiluted		
50%	1:2	1	1
25%	1:4	1	3
12.5%	1:8	1	7

NOTE: If desired, S.A.R.P. may be prediluted 1:2 with 0.85% saline to enable the use of a reference value which is one-half the published reference value. This choice increases sensitivity when undiluted reference value exceeds 100%. Under those circumstances, the prediluted material will be used as reconstituted S.A.R.P. and one-half the reference value would be used in the calculations.

- Dilute each patient sample and control with 0.85% saline. Prepare a 1:2 dilution (1 part patient plasma and 1 part saline) and a 1:4 dilution (1 part patient plasma and 3 parts saline). **NOTE:** Additional dilutions may be necessary depending on the patient history. Suspected von Willebrand's samples may need to be tested undiluted.
- Pour 200ml of dissolved buffer into each of the outer sections of the chamber (requires a total of 400ml buffer).
- Place one sponge wick in the buffer along each inner wall of the chamber.
- Remove any excess moisture from the wells of the plate. **NOTE:** Excess moisture on the plate can result in poor rockets.
- Apply 10µl of each patient sample or dilution to the designated wells. **NOTE:** Standard curve samples must be run on each plate. Duplicate applications of patient samples are advisable. When applying the samples to the plate wells, do not allow the pipette tip to touch the sides of the wells as this may cause damage.
- Allow 5 minutes for specimens to diffuse into the agarose.
- Place the plate, agarose side down, into the chamber on the sponge wicks. Place the application point (wells) toward the cathode (-) side
- Electrophorese the plates at a constant current of 16 mA per plate (2 plates = 32 mA) for 4 hours.
- Following electrophoresis, remove the plates from the chamber. **NOTE:** Discard the chamber buffer after each run. **NOTE:** Use the Humidity Chamber (Cat. No. 9036) as a staining and rinsing chamber. A laboratory rotator should not be used during the rinsing process.
- Rinse the plate with deionized water and wash in 0.85% saline, agarose side up, overnight.
- After the overnight wash, rinse the plate again and then wash in deionized water for 15 minutes. Drain or shake the excess water from the plate.
- Place the plate on a flat surface, agarose side up. Cover the agarose with a single, lint-free tissue.
- Place 4-5 Titan Blotter Pads and then a Development Weight on the plate for 5 minutes. Remove the weight and the top 2-3 Blotter Pads. Replace with 2 fresh Blotter Pads and then a weight for an additional 10 minutes and remove.
- Score the edges of the plate with a sharp instrument to prevent peeling.

16. Dry the plate in a drying oven at 60...70°C for 10-20 minutes. Periodically inspect the agarose for transparency. The plate will be transparent when completely dry. DO NOT over-dry the plates. If the agarose starts to peel at the edges before the entire plate is dry, reduce the temperature. **NOTE:** If a dryer/oven is not available, the plates may be covered with the wet lint-free tissue and allowed to dry at 15...30°C overnight or under a fan for 3 hours at 15...30°C.
17. Following drying, pour the stain over the entire surface of the plate, agarose side up, for 20 minutes.
18. Place the plate in destain solution for up to 3 minutes. **NOTE:** Discontinue destaining as soon as the background sufficiently clears in order to easily distinguish rocket peaks. Excessive destaining may fade the rockets making correct measurements difficult. If over destaining does occur, repeat Step F.8. and stain the rockets again.
19. Rinse the plate twice in deionized water for 5-10 minutes each rinse.
20. Dry the plates at 37°C for 5 minutes or at 15...30°C until dry. After drying, tape the edges of the agarose to the plastic backing with clear tape. This will prevent peeling of the agarose and prolong the life of the plate.
21. Place the plate on a light box or white paper and mark the apex of each rocket peak with a sharp pencil or marker.
22. Using the Helena Rocket Ruler, measure the length of each peak in millimeters. The peak is measured from the top of each well to the apex of the rocket. Each box of plates contains one Helena Rocket Ruler.
23. Plot the percent activity of the reference curve versus each rocket height on the Helena Factor VIII Related Antigen Report Form or on 2 cycle semi-logarithmic paper. Draw the 'line of best fit' for the four points. Refer to Figures 1 and 2 in INTERPRETATION OF RESULTS for an example of a completed Rocket Plate and a standard curve drawn on 2 cycle semi-logarithmic paper.

INTERPRETATION OF RESULTS

Read the patient or control values from the standard curve and multiply each by the appropriate dilution factor. Use of a pre-diluted standard would be treated the same with no change in the dilution factors used. The only difference would be the final reference value used. If Coagulation S.A.R.P. is used to prepare the standard curve, the patient value read from that curve must be multiplied by the assigned or pre-diluted von Willebrand Factor Antigen value of the appropriate lot of Coagulation S.A.R.P. as well as the dilution factor.

Patient value from curve = 30%
 S.A.R.P. assigned value (from assay sheet) = 95%
 Dilution factor = 2
 Actual Patient Factor VIII
 Related Antigen = $30\% \times 2 \times 0.95 = 57.0\%$

Unknown patient samples with von Willebrand Factor Antigen levels greater than the range of the standard curve must be reassayed using the appropriate dilutions.

VON WILLEBRAND FACTOR ANTIGEN ROCKET EIA

Figure 1: Rocket patterns on a von Willebrand Factor Antigen Rocket Plate. The lengths of the rockets (in millimeters) of the Coagulation S.A.R.P. (the standard) dilutions are used to prepare the standard curve. Patient results are read from the curve. In this illustration, Patient #3 has no rocket formation indicating a von Willebrand Factor Antigen deficiency.



Figure 2: Standard Curve prepared with Coagulation S.A.R.P. on 2 cycle semi-logarithmic paper.

The rocket lengths of the standards are as follows:

100% activity = 26 mm

50% activity = 21 mm

25% activity = 16 mm

12.5% activity = 11 mm

From the curve, a patient with a rocket length of 17 mm for 1:2 dilution would have 30% activity. See Step G.4. of the STEP-BY-STEP PROCEDURE for complete patient calculations.

Severe von Willebrand's patients often show no detectable rocket precipitation by this method. Other patients with von Willebrand's disease have been reported with normal levels of von Willebrand Factor Antigen⁴. The diagnosis of von Willebrand's syndrome should not be made on the basis of this test alone. Factor VIII coagulant activity, ristocetin aggregation, template bleeding times and a thorough history and physical are essential for an accurate diagnosis of this disease. On the basis of the results from various tests, von Willebrand's patients may be classified into several categories. Suspected von Willebrand's patients should be tested on several different occasions due to the variability in the test results observed. There are fluctuations in the levels of the immunologically detectable protein as well as in the functional activity in the Factor VIII molecule^{5,7,8}. Each laboratory should establish its own reference range for this procedure.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 SAC1

Cat. No. 5302 SAC2

REFERENCE VALUES

von Willebrand Factor Antigen values are generally expressed in relative percents as compared to a standard or pooled normal plasma. The expected normal range for the Helena von Willebrand Factor Antigen Rocket EIA Method is 50% to 160%. This range was established by testing a total of 50 individuals, 25 males and 25 females. Testing on male donors versus female donors showed no significant difference in mean values. This normal range compared favorably with normal ranges reported by Zimmerman⁴ and Johnson⁵.

PERFORMANCE CHARACTERISTICS

Helena BioSciences or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Precision Studies

Precision of the Helena von Willebrand Factor Antigen Rocket EIA Method was determined by reconstituting seven vials of Coagulation S.A.R.P. on eight different days and performing the test in duplicate each day. The S.A.R.P. was run undiluted, and in dilutions of 1:2, 1:4 and 1:8. The mean coefficient of variation was 7.8%.

Sensitivity

The sensitivity limits are defined by the von Willebrand Factor Antigen Reference Standards (S.A.R.P.). Patient values greater than the highest value on the standard curve must be diluted and reassayed. Patient values less than the lowest value on the curve must be reported as such but no absolute value may be extrapolated.

BIBLIOGRAPHY

1. Laurell, C.B., Electroimmuno Assay, *Scan J Clin Lab Invest*, 29 Suppl 124:21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies, *Anal Biochem*, 15:45-52, 1966.
3. Eyster, M.E. et al., Carrier Detection in Classic Hemophilia by Combined Measurement of Immunologic (VIII AGN) and Procoagulant (VIII AHF) Activities, *AJCP*, 65:975-981, 1976.
4. Zimmermann, T.S., et al., Determination of the von Willebrand's Disease Antigen (Factor VIII-Related Antigen) in Plasma by Quantitative Immunoelectrophoresis. *J Lab Clin Med* 86:152-159, 1975.
5. Johnson, S.S. et al., Newborn Factor VIII Complex: Elevated Activities in Term Infants and Alterations in Electrophoretic Mobility Related to Illness and Activated Coagulation, *Brit J Haem* 47:597-606, 1981.
6. Peake, I.R., et al., Inherited Variants of Factor VIII-Related Protein in von Willebrand's Disease, *N Eng J Med* 291:113, 1974.
7. Moehring, C., A Rapid Factor VIII-Related Antigen Electroimmunoassay, *AJCP*, 72:21-28, 1979.
8. Brown, J.E. et al., Effect of Exercise on the Factor VIII Complex: A Correlation of the von Willebrand Antigen and Factor VIII Coagulation Antigen Increase, *Thromb Res* 15:61-67, 1979.

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BIBLIOGRAPHIE

1. Laurell, C. B., Electroimmuno Assay, Scan J Clin Lab Invest 29 Suppl 124 : 21-37, 1972.
2. Laurell, C. B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies. Anal Biochem, 15 : 45-52, 1966.
3. Eyster, M. E. et al., Carrier Detection in Classic Hemophilia by Combined Measurement of Immunologic (VIII AGN) and Procoagulant (VIII AHF) Activities, AJCP, 65 : 975-981, 1976.
4. Zimmermann, T. S., et al., Determination of the von Willebrand's Disease Antigen (Factor VIII-Related Antigen) in Plasma by Quantitative Immunoelectrophoresis. J Lab Clin Med 86 : 152-159, 1975.
5. Johnson, S. S. et al., Newborn Factor VIII Complex: Elevated Activities in Term Infants and Alterations in Electrophoretic Mobility Related to Illness and Activated Coagulation, Brit J Haem 47 : 597-606, 1981.
6. Peake, I. R., et al., Inherited Variants of Factor VIII-Related Protein in von Willebrand's Disease, N Eng J Med 291 : 113, 1974.
7. Moehring, C., A Rapid Factor VIII-Related Antigen Electroimmunoassay, AJCP, 72 : 21-28, 1979.
8. Brown, J. E. et al., Effect of Exercise on the Factor VIII Complex: A Correlation of the von Willebrand Antigen and Factor VIII Coagulation Antigen Increase, Thromb Res 15 : 61-67, 1979.

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LITERATUR

1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest, 29 Suppl 124:21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies, Anal Biochem, 15:45-52, 1966.
3. Eyster, M.E. et al., Carrier Detection in Classic Hemophilia by Combined Measurement of Immunologic (VIII AGN) and Procoagulant (VIII AHF) Activities, AJCP, 65:975-981, 1976.
4. Zimmermann, T.S., et al., Determination of the von Willebrand's Disease Antigen (Factor VIII-Related Antigen) in Plasma by Quantitative Immunoelectrophoresis. J Lab Clin Med 86:152-159, 1975.
5. Johnson, S.S. et al., Newborn Factor VIII Complex: Elevated Activities in Term Infants and Alterations in Electrophoretic Mobility Related to Illness and Activated Coagulation, Brit J Haem 47:597-606, 1981.
6. Peake, I.R., et al., Inherited Variants of Factor VIII-Related Protein in von Willebrand's Disease, N Eng J Med 291:113, 1974.
7. Moehring, C., A Rapid Factor VIII-Related Antigen Electroimmunoassay, AJCP, 72:21-28, 1979.
8. Brown, J.E. et al., Effect of Exercise on the Factor VIII Complex: A Correlation of the von Willebrand Antigen and Factor VIII Coagulation Antigen Increase, Thromb Res 15:61-67, 1979.

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BIBLIOGRAFIA

1. Laurell, C.B., Electroimmuno Assay, *Scan J Clin Lab Invest*, 29 Suppl 124:21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies, *Anal Biochem*, 15:45-52, 1966.
3. Eyster, M.E. et al., Carrier Detection in Classic Hemophilia by Combined Measurement of Immunologic (VIII AGN) and Procoagulant (VIII AHF) Activities, *AJCP*, 65:975-981, 1976.
4. Zimmermann, T.S., et al., Determination of the von Willebrand Disease Antigen (Factor VIII-Related Antigen) in Plasma by Quantitative Immunoelectrophoresis. *J Lab Clin Med* 86:152-159, 1975.
5. Johnson, S.S. et al., Newborn Factor VIII Complex: Elevated Activities in Term Infants and Alterations in Electrophoretic Mobility Related to Illness and Activated Coagulation, *Brit J Haem* 47:597-606, 1981.
6. Peake, I.R., et al., Inherited Variants of Factor VIII-Related Protein in von Willebrand Disease, *N Eng J Med* 291:113, 1974.
7. Moehring, C., A Rapid Factor VIII-Related Antigen Electroimmunoassay, *AJCP*, 72:21-28, 1979.
8. Brown, J.E. et al., Effect of Exercise on the Factor VIII Complex: A Correlation of the von Willebrand Antigen and Factor VIII Coagulation Antigen Increase, *Thromb Res* 15:61-67, 1979.

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1. Laurell, C.B., Electroimmuno Assay, Scan J Clin Lab Invest, 29 Suppl 124:21-37, 1972.
2. Laurell, C.B., Quantitative Estimation of Proteins by Electrophoresis in Agarose Gel Containing Antibodies, Anal Biochem, 15:45-52, 1966.
3. Eyster, M.E. et al., Carrier Detection in Classic Hemophilia by Combined Measurement of Immunologic (VIII AGN) and Procoagulant (VIII AHF) Activities, AJCP, 65:975-981, 1976.
4. Zimmermann, T.S., et al., Determination of the von Willebrand Disease Antigen (Factor VIII-Related Antigen) in Plasma by Quantitative Immuno-electrophoresis. J Lab Clin Med 86:152-159, 1975.
5. Johnson, S.S. et al., Newborn Factor VIII Complex: Elevated Activities in Term Infants and Alterations in Electrophoretic Mobility Related to Illness and Activated Coagulation, Brit J Haem 47:597-606, 1981.
6. Peake, I.R., et al., Inherited Variants of Factor VIII-Related Protein in von Willebrand Disease, N Eng J Med 291:113, 1974.
7. Moehring, C., A Rapid Factor VIII-Related Antigen Electroimmunoassay, AJCP, 72:21-28, 1979.
8. Brown, J.E. et al., Effect of Exercise on the Factor VIII Complex: A Correlation of the von Willebrand Antigen and Factor VIII Coagulation Antigen Increase, Thromb Res 15:61-67, 1979.

Helena Biosciences Europe
Queensway South
Team Valley Trading Estate
Gateshead
Tyne and Wear
NE11 0SD

Tel: +44 (0) 191 482 8440
Fax: +44 (0) 191 482 8442
Email: info@helena-biosciences.com
www.helena-biosciences.com

HL-2-1415P 2007/05 (3)